

ThinPrep® CytoLyt® Solution

Instructions for Use

**INTENDED USE**

CytoLyt® Solution is a water-based, buffered cell wash solution designed to lyse red blood cells, prevent protein precipitation, dissolve mucus, and preserve morphology of general cytology cellular samples prior to slide preparation with the ThinPrep® processor.

Note: CytoLyt Solution is not for use with samples processed with the ThinPrep 3000 processor.

SUMMARY AND EXPLANATION

CytoLyt Solution is designed for use with the ThinPrep® processor, a cytologic preparation device that produces slides for microscopic examination. CytoLyt Solution enables the transport and preservation of cells for up to eight days at 15°C–30°C.

PRINCIPLES OF PROCEDURE

CytoLyt Solution is a media used for collection and lysing of red blood cells. When used on samples processed with a ThinPrep processor, the process allows transfer of cells onto a microscope slide, providing a thin, uniform layer of cells suitable for cytologic evaluation.

COMPOSITION

Methanol-based, buffered preservative solution. CAS 67-56-1

WARNINGS

Danger. Flammable. Contains Methanol.

H302 - Harmful if swallowed.

H312 - Harmful in contact with skin.

H332 - Harmful if inhaled.

H370 - Causes damage to organs.

H226 - Flammable liquid and vapor.

For In Vitro Diagnostic use. Use with adequate ventilation. Not for external or internal use in humans or animals. Cannot be made non-poisonous.

PRECAUTIONS

P210 - Keep away from heat/sparks/open flames/hot surfaces.

Store in a well-ventilated place. Keep container tightly closed. Wear gloves and eye protection when handling solution. Wash hands thoroughly after handling. Do not use if primary packaging is damaged.

When transporting a CytoLyt Solution vial containing cells, make sure the vial is tightly sealed. Align the mark on the cap with the mark on the cup or tube to prevent leakage.

As with all laboratory procedures, universal precautions should be followed.

ThinPrep® CytoLyt® Solution

PRETREATMENT

No reconstitution, mixing or dilution is required.

STORAGE

Store CytoLyt Solution without cytologic samples at 15°C to 30°C (59°F to 86°F). Do not use CytoLyt Solution beyond the expiration date marked on the container. Close the 946 mL bottle after each use.

APPEARANCE AND INTEGRITY

Clear, non-sterile solution.

SPECIMEN COLLECTION AND PREPARATION

Collect non-gynecologic samples in a routine manner and refer to the ThinPrep processor operator's manual for preparation instructions. Record required patient information in the space provided. (See figure.)

Storage and Handling

CytoLyt Solution preserves cells for up to eight days at temperatures between 15°C (59°F) and 30°C (86°F).

PROCESSING INSTRUCTIONS

Cytologic specimens collected in CytoLyt Solution are further prepared with PreservCyt® Solution and then processed on a ThinPrep processor according to instructions in the ThinPrep processor's operator's manual.

LIMITATIONS OF PROCEDURE

CytoLyt Solution cannot be substituted with any other solution for specimen collection, preparation, or processing on any ThinPrep processor.

PERFORMANCE CHARACTERISTICS


Refer to the ThinPrep processor operator's manual.

DISPOSAL

Dispose in accordance with all applicable regulations.

FIRST AID MEASURES

IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. See www.hologicsds.com for the entire Safety Data Sheet.



Figure

ThinPrep® CytoLyt® Solution



Manufacturer



Consult instructions for use



In vitro diagnostic medical device



Authorised representative in the European Community



Catalogue number



Use by



Batch code



Temperature limitation



Flammable



Irritant, Respiratory Tract Irritation, Dermal Sensitizer



Respiratory Sensitizer, Target Organ Toxicity



Hologic, Inc. • 250 Campus Drive • Marlborough, MA 01752 USA
Tel: +1 508 263-2900 • www.hologic.com

EC REP Hologic BVBA • Da Vinciaan 5 • 1930 Zaventem • Belgium

Australian Sponsor: Hologic (Australia) Pty Ltd • Suite 302, Level 3 • 2 Lyon Park Road
Macquarie Park NSW 2113 • Australia • Tel: 02 9888 8000