FOR EXPORT ONLY

LDRCheck-D

For Quality Control of GEN-PROBE AccuLDR, LEADER 50, LEADER 50*i*, and LEADER I Luminometers

INTENDED USE

The GEN-PROBE LDRCheck-D reagent is a chemiluminescent DNA control used to qualitatively assess performance of GEN-PROBE AccuLDR, LEADER 50, LEADER 50*i*, and LEADER I luminometers.

SUMMARY AND EXPLANATION OF THE TEST

Calibration of GEN-PROBE AccuLDR, LEADER 50, LEADER 50*i*, and LEADER I luminometers is set using the GEN-PROBE Tritium Standard. The LDRCheck-D control is used to monitor performance of the luminometer every two weeks by measuring the control in the GEN-PROBE LEADER luminometer. On the first day of use, the test result is calculated as a ratio of the Observed Value of LDRCheck-D reagent to the Expected label value of LDRCheck-D reagent. The test result is compared to the quality control limits for the given LEADER luminometer model. The Observed Day 0 value is recorded and used for trending purposes. The instrument variation from test to test should fall within a ratio of 0.90 to 1.10 relative to the Observed Day 0 value.

PRINCIPLES OF THE PROCEDURE

The concentration of LDRCheck-D reagent has been determined and appears on the label as RLU / 100 μ L. LDRCheck-D reagent is measured in the luminometer being tested and the test results are compared against the quality control limits.

REAGENTS

GEN-PROBE LDRCheck-D Quality Control Reagent for LEADER Luminometers.

< 100 ng/bottle non-infectious DNA with a chemiluminescent label in succinate buffered solution containing detergent.

WARNINGS AND PRECAUTIONS

- A. For *in vitro* diagnostic use.
- B. Avoid contact of LDRCheck-D reagent with skin, eyes and mucous membranes. Wash with water if these reagents come into contact with skin or eyes. If spills of these reagents occur, dilute with water and wipe dry.
- C. Use only supplied or specified disposable laboratory ware.
- D. Avoid contact of Detection reagents I and II (GEN-PROBE Detection Reagent Kit) with skin, eyes and mucous membranes. **WARNING: IRRITANTS, CORROSIVES.** Wash with water if contact with these reagents occurs. If spills of these reagents occur, dilute with water before wiping dry.
- E. Use only GEN-PROBE Detection Reagent Kit (Cat. No. 1791, bioMérieux ref. 39300) with LDRCheck-D reagent.

STORAGE AND HANDLING REQUIREMENTS

- A. LDRCheck-D reagent should be capped when not in use to avoid changes in concentration due to evaporation.

- B. LDRCheck-D reagent is stable when stored at 2° to 25°C until the expiration date indicated. Storage above room temperature may result in product deterioration.
- C. The LDRCheck-D test is to be performed at room temperature.
- D. Avoid contaminating the LDRCheck-D reagent by contact with used or contaminated pipette tips.

MATERIALS

A. MATERIALS PROVIDED

Catalog No. 301042F, bioMérieux Ref. 39318. 20 Tests
GEN-PROBE LDRCheck-D Quality Control Reagent for AccuLDR (formerly PAL), LEADER 50, LEADER 50*i*, and LEADER I Luminometers 1 x 20 mL

B. MATERIALS REQUIRED BUT NOT PROVIDED

GEN-PROBE LEADER or AccuLDR (formerly PAL) luminometer
GEN-PROBE Detection Reagent Kit, Cat. No. 1791, bioMérieux Ref. 39300
GEN-PROBE Reaction Tubes, Cat. No. 2065, bioMérieux Ref. 39307
Calibrated micropipets capable of delivering 100 μ L or Eppendorf repeat pipettor
Micropipets capable of delivering 100 μ L and/or Eppendorf Repeating Pipettor with Combitips (0.5 mL or 1.25 mL)

C. MATERIALS AVAILABLE FROM GEN-PROBE

GEN-PROBE LEADER 50*i* luminometer, Cat. No. 3100i, bioMérieux Ref. 39400
GEN-PROBE Detection Reagent Kit, Cat. No. 1791, bioMérieux Ref. 39300
GEN-PROBE Reaction Tubes, Cat. No. 2065, bioMérieux Ref. 39307

TEST PROCEDURE

A. REAGENT PREPARATION

1. If stored at $\leq 15^{\circ}\text{C}$, prior to use, bring the LDRCheck-D reagent to room temperature by equilibrating at room temperature for one hour. Alternatively, swirl the bottle for 30 seconds in a water bath at $60^{\circ}\pm 1^{\circ}\text{C}$ and cool at room temperature for ≥ 10 minutes. Once equilibrated, LDRCheck-D reagent may be stored at room temperature.
2. Before pipetting sample, mix the LDRCheck-D reagent by swirling and gently inverting three times to mix. Do not vortex.

B. EQUIPMENT PREPARATION

1. Prepare the GEN-PROBE luminometer for operation. Make sure there is sufficient volume of Detection Reagents I and II to complete the tests.
2. Print out the instrument parameters as described in the Operator's Manual. On AccuLDR read and record instrument parameters.
3. Prepare the luminometer to read in the raw data mode using a 1 second count time as described in the instrument operator's manual.

C. CHEMICAL BACKGROUND CHECK

- Perform a chemical background check on the Detection Reagents.
1. Run the instrument wash protocol as described in the Operator's Manual.
 2. Read ten empty tubes for a one second count time in raw data mode. The exterior of each tube is to be wiped with a water dampened, lint-free wipe before reading in the luminometer.
 3. If the average reading is: ≥ 150 RLU for LEADER I, LEADER 50, LEADER 50*i*; or ≥ 100 PLU for AccuLDR, then replace the Detection Reagents with a fresh set of detection reagents and

repeat the test. If the retest is out, notify your Gen-Probe distributor Customer Service.

D. MEASURING LDRCHECK-D REAGENT

1. LDRCheck-D reagent should be read on the luminometer a minimum of every two weeks.
2. Precision pipette ten replicates of LDRCheck-D reagent (100 µL per tube) into clean polystyrene tubes.
3. Read the LDRCheck-D samples plus 10 empty tubes using a one second count time on raw data mode.
4. When the analysis is complete, remove the tubes from the luminometer.

PROCEDURAL NOTES

- A. CALIBRATION OF LUMINOMETER: The instrument calibration is set using Tritium Standard as described in the instrument Operator's Manual.
- B. Use calibrated micropipets to pipette LDRCheck-D reagent.
- C. Pipette all samples to the bottom of the tube taking care not to splash on the side of the tube in order to obtain accurate results.
- D. The exterior of each tube is to be wiped with a water dampened, lint free wipe before being read in the luminometer.
- E. The aliquots of LDRCheck-D reagent should appear as a clear liquid with no clumps. If microbial contamination is apparent, discard bottle. Obtain a new bottle of LDRCheck-D reagent and proceed.
- F. As in any reagent system, excess powder on some gloves may cause contamination of opened reagents or reaction tubes. Gen-Probe recommends that customers experiencing difficulty with the test avoid using this type of laboratory glove. Using powderless gloves (no talcum powder) will avoid this difficulty.

RESULTS

A. COEFFICIENT OF VARIATION (%CV):

The Coefficient of Variation (%CV) for LDRCheck-D reagent should be $\leq 6.0\%$. If the %CV is greater than 6.0% then prepare and read a new set of LDRCheck-D replicates. If the %CV is still greater than 6.0%, then contact your Gen-Probe distributor Customer Service.

For AccuLDR calculate the %CV as follows:

$$\%CV = \frac{\text{Standard Deviation of LDRCheck-D in RLU} \times 100\%}{\text{Average RLU} \times \text{LDRCheck-D}}$$

* For AccuLDR to convert PLU to RLU: $RLU = PLU \div 0.031$

B. CALCULATION OF RESULTS

1. DAY 0 - INITIAL OBSERVED/EXPECTED RATIO
Calculate the Observed Value / Expected Value ratio using the equation below:

$$\frac{\text{Observed Value}}{\text{Expected Value}} = \frac{(\text{Avg. RLU} \times \text{LDRCheck-D} - \text{Avg. RLU} \times \text{empty tubes})}{(\text{LDRCheck-D Value on label})}$$

2. SUBSEQUENT TEST DATES – OBSERVED/DAY 0 RATIO
Calculate the Observed Value / Day 0 Value ratio using the equation below:

$$\frac{\text{Observed Value}}{\text{Day 0 Value}} = \frac{(\text{Avg. RLU} \times \text{LDRCheck-D} - \text{Avg. RLU} \times \text{empty tubes})}{(\text{LDRCheck-D net RLU} \times \text{observed on Day 0})}$$

* For AccuLDR to convert PLU to RLU: $RLU = PLU \div 0.031$

QUALITY CONTROL RESULTS AND ACCEPTABILITY

- A. DAY 0 - INITIAL OBSERVED/EXPECTED RATIO
The Observed/Expected (O/E) ratio must be within the range indicated below. If the O/E ratio is within this range, the instrument is

functioning properly. If the O/E range is not met, perform a second quality control test. If the O/E range is still not met, contact your Gen-Probe distributor Customer Service.

Luminometer	O/E
AccuLDR	0.80 – 1.10
LEADER 50	0.90 – 1.20
LEADER 50i	0.90 – 1.20
LEADER I	0.80 – 1.10

B. SUBSEQUENT TEST DATES – OBSERVED / DAY 0 RATIO

The Observed / Day 0 ratio must be between 0.90 and 1.10. If the O/E ratio is within this range, the instrument is functioning properly. If the O/E range is not met, review all instructions and calculations and perform a second quality control test. If the O/E range is still not met, contact your Gen-Probe distributor Customer Service.

- C. This procedure should be followed whenever a new bottle of LDRCheck-D reagent is to be opened and used for checking instrument calibration.

LIMITATIONS

- A. LDRCheck-D reagent contains EDTA. The reagent is bacteriostatic, however the presence of fungal agents may affect the RLU value of the reagent. If contamination is observed or suspected, the reagent should not be used and should be replaced by a previously unopened bottle of LDRCheck-D reagent.
- B. Failure to follow the instructions given in this insert may result in erroneous results.

Manufactured by
Gen-Probe Incorporated
San Diego, California 92121 (USA)

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© 1999 Gen-Probe Incorporated