# ***Tumark<sup>®</sup> Vision for Brevera<sup>®</sup> Petite***

**REF** TUMARK-BREV-P-VISION

# ***Tumark<sup>®</sup> Vision for Brevera<sup>®</sup> Standard***


**REF** TUMARK-BREV-S-VISION

# ***Tumark<sup>®</sup> Vision for Eviva<sup>®</sup> Petite***

**REF** TUMARK-E13-P-VISION

# ***Tumark<sup>®</sup> Vision for Eviva<sup>®</sup> Standard***

**REF** TUMARK-E13-S-VISION

**INSTRUCTIONS FOR USE**  
GEBRUIKSAANWIJZING  
CONSIGNES D'UTILISATION  
GEBRAUCHSANWEISUNG  
ISTRUZIONI PER L'USO  
INSTRUÇÕES DE USO  
INSTRUCCIONES DE USO

**HOLOGIC<sup>®</sup>**

**CONTENT**

INHOUD / CONTENU / INHALT / CONTENUTO / CONTENÚDO / CONTENIDO:

ENGLISH - INSTRUCTIONS FOR USE .....	2
NEDERLANDS - GEBRUIKSAANWIJZING.....	6
FRANÇAIS - CONSIGNES D'UTILISATION.....	10
DEUTSCH - GEBRAUCHSANWEISUNG.....	14
ITALIANO - ISTRUZIONI PER L'USO .....	18
PORTUGUÊS - INSTRUÇÕES DE USO.....	22
ESPAÑOL - INSTRUCCIONES DE USO.....	26
SYMBOLS .....	30
ORDERING .....	32

**Compatibility table:**

Hologic biopsy device compatibility chart:

Biopsy Device		Com- patible with	Marker Application System		
Name*	REF		Name	REF	Length
Brevera® Biopsy Device	BREVDISP09, 12 mm variable aperture setting	→	Tumark® Vision for <b>Brevera® Petite</b>	TUMARK-BREV- P-VISION	123.5 mm
Brevera® Biopsy Device	BREVDISP09, 20 mm variable aperture setting	→	Tumark® Vision for <b>Brevera® Standard</b>	TUMARK-BREV- S-VISION	127.5 mm
Eviva® Biopsy Device (Petite)	Eviva 0913-12; Eviva 0913-12T	→	Tumark® Vision for <b>Eviva® Petite**</b>	TUMARK-E13-P- VISION	133 mm
Eviva® Biopsy Device (Standard)	Eviva 0913-20; Eviva 1213-20	→	Tumark® Vision for <b>Eviva® Standard**</b>	TUMARK-E13-S- VISION	129.5 mm

\* Biopsy device is sold with compatible biopsy introducer

\*\* Tumark Vision for Eviva® Standard and Tumark Vision for Eviva® Petite are not compatible with Eviva® 10 cm biopsy devices

**ENGLISH - INSTRUCTIONS FOR USE***Read instructions before use**Keep for future reference***Compatibility table:**

Hologic biopsy device compatibility chart:

Biopsy Device		Com- patible with	Marker Application System		
Name*	REF		Name	REF	Length
Brevera® Biopsy Device	BREVDISP09, 12 mm variable aperture setting	→	Tumark® Vision for <b>Brevera® Petite</b>	TUMARK-BREV- P-VISION	123.5 mm
Brevera® Biopsy Device	BREVDISP09, 20 mm variable aperture setting	→	Tumark® Vision for <b>Brevera® Standard</b>	TUMARK-BREV- S-VISION	127.5 mm
Eviva® Biopsy Device (Petite)	Eviva 0913-12; Eviva 0913-12T	→	Tumark® Vision for <b>Eviva® Petite**</b>	TUMARK-E13-P- VISION	133 mm
Eviva® Biopsy Device (Standard)	Eviva 0913-20; Eviva 1213-20	→	Tumark® Vision for <b>Eviva® Standard**</b>	TUMARK-E13-S- VISION	129.5 mm

\* Biopsy device is sold with compatible biopsy introducer

\*\* Tumark Vision for Eviva® Standard and Tumark Vision for Eviva® Petite are not compatible with Eviva® 10 cm biopsy devices

**Important Information:**

Read the instruction manual thoroughly and be familiar with its contents prior to use. Failure to read the entire manual and familiarize yourself with all instructions before using the device is unsafe and can result in life threatening or severe injury to the patient or user and to damage or malfunction of the device.

**Indications:**

The Tumark® Vision for Eviva® and Tumark® Vision for Brevera® are intended to attach a marker to soft tissue at the surgical site during a percutaneous procedure. The devices are indicated for use to radiographically and radiologically mark the surgical location in breasts following a percutaneous procedure. They are not indicated to be used with magnetic resonance imaging (MRI) techniques.

**Contraindications:**

The products are not intended for use except as indicated above.

The use of the products is contraindicated in patients who suffer from a severe nickel allergy.

**A summary of safety and clinical performance in accordance with Article 32 of Regulation 2017/745 (implantable devices) is available on [www.somatex.de](http://www.somatex.de).**

**Warnings:**

- Consult instructions for use of the biopsy device.
- Only qualified physicians with knowledge, experience and training in percutaneous soft tissue marking should use the product.
- This manual does not include descriptions or instructions for surgical techniques. It is the responsibility of the physician performing any procedure to determine the appropriateness of the procedure to be performed and of the use of this device and to determine the specific technique for each patient.
- DO NOT use the system in patients with breast implants.
- The cannula has a blunt tip. DO NOT use the cannula without a compatible biopsy introducer after a percutaneous biopsy procedure. Consult the compatibility table.
- The product is only sterile if the indicator on the package is green, if used before the expiration date and if package is unopened and undamaged. DO NOT use if indicator is not green, after the expiration date or if package is open or damaged.
- Single patient use only. DO NOT reuse or resterilize. We do not assume any liability for the use of this product or its components in case of re-sterilisation or reuse. This product may not be reused after a single application. The quality of the materials, coats and adhesive joints could degrade. Safe use is not guaranteed any longer. The product that is already used once is not designed for the required cleaning and sterilisation processes. The sterility of the reprocessed disposable products is therefore not guaranteed. The risk of unwanted injuries and infections, especially cross-infections between patient and medical staff inappropriately increases.
- The marker **application system** is **NOT** suitable for deployment under **MRI**.

**Precautions:**

- Make sure that the slider does not change its backward position while placing the cannula.
- The marker must be placed by pushing the slider forward as far as possible to the stop position.
- Consult the compatibility table for use together with introducers after a percutaneous biopsy procedure.
- Do use the system only with the referenced products.
- **Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician.**
- In rare cases the expansion of the marker may be delayed. Visibility in radiological imaging might be compromised until full expansion.

**Information about materials used:**


The implantable clip marker is made from a nickel-titanium alloy (Nitinol).

**Any serious events that occur in relation to the product should be reported to SOMATEX Medical Technologies GmbH as well as the competent national authority.**


**Device description:**

The product is a sterile, single use, preloaded tissue site marking systems consisting of a non-absorbable nickel-titanium marker (1), an introducer cannula (2) and a plastic handheld applicator (4) with deployment mechanism.

The introducer cannula has a blunt tip and can only be used together with an introducer. The handle is equipped with a slide-button (3) which allows for a one handed placement of the marker by pressing it forward. A safety catch system prevents the slide-button to inadvertently move forward and therefore prevents a premature deployment of the marker.

*Schematic illustration*

*Dimensions of clip marker*

**MRI Safety Information application system:**

MR unsafe

The *Tumark® Vision for Eviva®* and *Tumark® Vision for Brevera®* as an application system for clip markers are **not** suitable for use in MRI.

**MRI Safety Information clip marker:**

MR conditional

Clip markers, which have already been placed inside a patient, are conditionally MR safe. A patient with a clip marker can be safely scanned in an MR system meeting the following conditions:

- static magnetic field up to 3.0 Tesla with
- theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg.

Non-clinical tests were performed on the following systems:

- 3 Tesla Siemens Magnetom Skyra (Siemens Medical, Erlangen, Germany) MRI with software Numaris 4, syngo MR (Version “D11, N4\_VB11D\_LATEST\_20110129\_P3”);
- 1.5 Tesla Siemens Magnetom Avanto (Siemens Medical, Erlangen, Germany) MRI with software Numaris 4, syngo MR (Version “B13, N4\_VB13A\_LATEST\_20060607\_P29”).

Under the scanning conditions defined above, it is expected that clip marker of *Tumark® Vision for Eviva®* and *Tumark® Vision for Brevera®* will produce the following maximum RF-related temperature rise:

- at 1.5 Tesla: < 1,0 °C (2 W/kg SAR) after 20 min of continuous scanning,
- at 3.0 Tesla: < 1,0 °C (2 W/kg SAR) after 20 min of continuous scanning.

Under the scanning conditions defined above, it is expected that clip marker of *Tumark® Vision for Eviva®* and *Tumark® Vision for Brevera®* will produce the following image artefacts:

- at 1.5 Tesla: 4.3 mm spin echo sequence; 5.2 mm gradient echo sequence;
- at 3.0 Tesla: 5.2 mm spin echo sequence; 5,8 mm gradient echo sequence.

**Do not expose the implanted clip marker to unconventional and non-standardized MRI techniques other than the ones listed above, because it has NOT BEEN TESTED for that purpose.**

### **Device Preparation and Use:**

**Warning:** The cannula has a blunt tip. DO NOT use the cannula without a compatible biopsy introducer after a percutaneous biopsy procedure. Consult the compatibility table.

1. Prior to use, ensure that the sterilization indicator is green and that the packaging has not been opened and / or damaged. Do not use the product if the packaging or the product appears to have been compromised. Check the sterilization expiration date. Do not use the product if it has expired.
2. Check the marker application system length and consult the compatibility table to make sure the cannula length is compatible with your vacuum biopsy system.
  - Note that the marker **application system** is **NOT** suitable for deployment under **MRI**.
  - Consider the size of the clip marker in relation to the area to be marked (see picture *dimensions of clip marker*).
3. Remove the product from the packaging.
4. Remove the biopsy device, but leave the introducer sheath in the breast.
5. Place the marker application system through the hub of the introducer sheath.
6. Advance the marker application system until it stops or meets the proximal end of the introducer sheath.
7. Deploy the marker by advancing the slider forward **as far as possible to the stop position**.
8. Verify the deployment and proper position of the marker prior to removal of the device with the appropriate imaging modality (ultrasound and X-ray). Document the location of the marker.
9. Slowly remove the introducer sheath and marker application system as one unit from the breast.
10. After use: dispose the application device properly, following internal guidelines if appropriate; however, at least one suitable container intended for contaminated cannulas should be provided to ensure safe disposal.

### **Storage Instructions:**

Keep dry.

Keep away from sunlight and heat (temperature between 41 °F – 86 °F / 5 °C - 30 °C).

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Hologic, Inc.