ATEC®

Breast Biopsy and Excision System Ultrasound Introducer

Instructions For Use

MAN-08968 Revision 001

- GR Σύστημα βιοψίας και εκτομής μαστού, Εισαγωγέας υπερήχων, Οδηγίες χρήσης MAN-08968, Αναθεώρηση 001
- NL Biopsie- en excisiesysteem voor de borst, Echografische inbrenger, Gebruiksaanwijzing MAN-08968, Revisie 001
- SV Bröstbiopsisystem och excisionsystem, Ultraljudsintroducerare, Bruksanvisning MAN-08968, Revidering 001
- FI Rintabiopsia- ja poistoleikkausjärjestelmä, Ultraäänisisäänviejä, Käyttöohjeet MAN-08968, Tarkistettu versio 001
- NO System til brystbiopsi og -eksisjon, Innføringshylse for ultralyd, Bruksanvisning MAN-08968, Oppdatering 001
- DA System til brystbiopsi og excision, Ultralydsintroducer, Brugervejledning MAN-08968, Revision 001
- FC Système de biopsie et d'excision mammaire, Introducteur à ultrasons, Mode d'emploi MAN-08968, Révision 001
- PB Sistema de Biópsia e Excisão Mamária, Introdutor do ultrassom, Instruções de uso MAN-08968, Revisão 001
- EI Sistema de biopsia y escisión de mama, Introducutor por ultrasonido, Instrucciones de uso MAN-08968, Revisión 001
- CS Systém pro biopsii a excizi prsu, Ultrazvukový zavaděč, Návod k použití MAN-08968, Revize 001
- HU Emlőbiopsziás és -exciziós rendszer, Ultrahangos bevezetőeszköz, Használati utasítás MAN-08968, Átdolgozás 001
- RO Sistem de biopsie și excizie mamară, Dispozitiv de introducere cu ultrasunete, Instrucțiuni de utilizare MAN-08968, Revizuirea 001
- SK Systém na biopsiu a excíziu prsníka, Ultrazvukový zavadzač, Návod na použitie MAN-08968, Revízia 001
- SL Sistem za biopsijo dojke in ekscizijo, Ultrazvočni uvajalnik, Navodila za uporabo MAN-08968, Revizija 001
- PL System do biopsji piersi i mamentomii, Introduktor ultradźwiękowy, Instrukcja obsługi MAN-08968, Wersja 001
- TR Meme Biyopsisi ve Eksizyon Sistemi, Ultrason İntrodüseri, Kullanım Talimatları MAN-08968, Revizyon 001
- SR Систем за биопсију и ексизију дојке, Ултразвучни уводник, Упутство за употребу MAN-08968, Ревизија 001
- PT Sistema de excisão e biopsia mamária, Introducutor ecográfico, Instruções de utilização MAN-08968, Revisão 001
- IT Sistema per biopsia e asportazione di tessuto mammario, Introduuttore a ultrasuoni, Istruzioni per l'uso MAN-08968, Revisione 001
- DE Brustbiopsie und Exzisionssystem, Ultraschall-Einführhilfe, Gebrauchsanweisung MAN-08968, Revision 001
- FR Système d'excision et de biopsie mammaire, Introducuteur à ultrasons, Mode d'emploi MAN-08968, Révision 001

HOLOGIC®

ATEC® Ultrasound Introducer

Instructions for Use

Please read all information carefully. Failure to properly follow the instructions may lead to unintended consequences.

Important: This package insert is designed to provide instructions for use (IFU) for the ATEC® ultrasound introducer. It is not a reference to surgical techniques.

Patient Target Group

The target patient population of ATEC Breast Biopsy and Excision System includes patients undergoing breast tissue removal for diagnostic sampling of breast abnormalities.

Indications

The ATEC ultrasound introducer is indicated to be used as an aspiration needle and/or as an introducer with another biopsy needle in percutaneous biopsies of the breast for diagnostic purposes only.

Contraindications

Sound professional judgement should be used when performing biopsies on patients who are using anticoagulant therapy and/or patients with known hemophilia problems.

Device Description

The ATEC ultrasound introducer is for use in conjunction with the ATEC Breast Biopsy and Excision System to obtain percutaneous core biopsy samples from soft tissue. The ATEC ultrasound introducer is designed specifically for use with the 9cm and 12cm length, 20mm aperture ATEC breast biopsy devices. The ATEC ultrasound introducer is a sterile single-use system of devices comprised of a stylet with an echogenic aperture and introducer. The ATEC ultrasound introducer provides single insertion access to the targeted area of interest and assists in deploying the appropriate biopsy site marker.

Intended user

The ATEC Breast Biopsy and Excision System should be used only by physicians trained in open or percutaneous biopsy procedures.

Device Preparation

(All steps to be performed per standard interventional technique)

1. Standard aseptic patient preparation should be employed prior to the use of the ATEC ultrasound introducer.
2. Prior to use of the ATEC ultrasound introducer, inspect the packaging and device to verify that neither has been damaged during shipment. If it appears that the packaging has been compromised, do not use the device.
3. Carefully remove the ATEC ultrasound introducer from its protective packaging and remove its protective sheath.
4. If a 9cm ATEC ultrasound introducer is being used with a CeleroMark, SecurMark for Celero, or SecurMark Ultrasound biopsy site marker, remove the spacer provided in the packaging and install it onto the biopsy site marker until it engages into place. (Figure 1)
5. Set-up the appropriate ATEC handpiece per standard operating procedures to perform an ultrasound-guided biopsy utilizing the ATEC Breast Biopsy and Excision System.


Figure 1

*Note: A luer taper or luer lock syringe may be attached to provide additional fluid and/or pain management.

Expected Clinical Benefit


The ATEC Breast Biopsy and Excision System enables physicians to use a minimally invasive system to extract small biopsy samples of potentially malignant breast tissue.

Instructions for Use

1. Identify the target area under ultrasound visualization.
2. Insert the tip of the stylet and advance to the desired location. The echogenic side aperture will be visible under ultrasound-guidance and should be placed in the same area of interest where the ATEC handpiece will be used.

Note: The Introducer may be attached directly to the ATEC breast biopsy device to gain access to the area of interest without use of the stylet.

3. Remove the stylet, leaving the introducer in place. The stylet may be detached from the introducer by depressing the thumb latch mechanism on the introducer hub.
4. Advance the ATEC handpiece to the targeted area by sliding it into the introducer hub until the thumb latch engages the biopsy device hub. (Figure 2)
5. Prior to performing a biopsy using the ATEC handpiece, verify aperture position using ultrasound-guidance.
6. Once the ATEC handpiece is inserted into the introducer hub and latched, follow the instructions for handpiece use in the ATEC Console Operator's Manual to perform the biopsy.

ATEC handpiece with introducer attached

Figure 2

7. Once core acquisition has been completed, return the console to biopsy mode and remove the ATEC handpiece. The handpiece may be detached from the introducer by depressing the latch mechanism on the introducer hub. The ATEC ultrasound introducer can be left in place for biopsy site marker placement, otherwise remove it with the ATEC handpiece.
8. For biopsy site marker deployment, insert the deployment device into the proximal end of the introducer and carefully advance the deployment device until it engages the introducer hub thumb latch. (Figure 3)


Figure 3

9. Deploy the biopsy site marker per the instructions for use provided with the biopsy site marker.

10. Remove the introducer and biopsy site marker as a unit and properly dispose of them.

Warnings and Precautions

- The ATEC ultrasound introducer is for diagnostic use only and is NOT indicated for therapeutic use.
- The ATEC ultrasound introducer is NOT intended to be used outside the standards, specifications, and limitations as outlined in this IFU and other related materials.
- The ATEC ultrasound introducer is not recommended for use with MRI.
- As with any medical procedure, please ensure that users wear appropriate personal protective equipment to guard against potential contact with body fluids.
- Observe caution while handling sharp components.
- Breast biopsies using the ATEC ultrasound introducer should be performed only by persons having adequate training and familiarity with this procedure. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- This device should be used only by physicians trained in percutaneous biopsy procedures.
- Sound professional judgment should be used when the ATEC ultrasound introducer is used on patients with breast implants.
- Avoid operator or instrument contact with the needle portion of the ATEC ultrasound introducer.
- Minimally invasive instruments and accessories manufactured or distributed by companies not authorized by Hologic, Inc. may not be compatible with the introducer set. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Dispose of all opened instruments whether used or unused.
- Do not resterilize the ATEC ultrasound introducer. Re-sterilization and/or reuse may compromise the integrity of the instruments. This may lead to potential failure of the device to perform as intended and/or cross-contamination associated with using inadequately cleaned or sterilized devices.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- **Rx ONLY** Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

How Supplied

The ATEC ultrasound introducer is sterilized by radiation and supplied sterile and packaged for single patient use. Discard into an appropriate container after use.

As Identified on Labels:

QTY Number of Devices Enclosed.

YYYY-MM-DD Expiration date is represented by the following:

YYYY represents the year

MM represents the month

DD represents the day

Product Complaints and Technical Support

Report any complaints or problems in the quality, reliability, safety, or performance of this product to Hologic. If the device has caused or added to patient injury, immediately report the incident to Hologic Authorized Representative and Competent authority of the respective member state or country. The Competent Authorities, for medical devices, are usually the individual Member States' Ministry of Health, or an agency within the Ministry of Health.

For More Information

For technical support or reorder information in the United States, please contact:

Hologic, Inc.
250 Campus Drive
Marlborough, MA 01752 USA
Phone: 877-371-4372
BreastHealth.Support@hologic.com

International customers, contact your distributor or local Hologic Sales Representative:

EC REP European Representative
Hologic BV
Da Vincilaan 5
1930 Zaventem
Belgium
Tel: +32 2 711 46 80

Symbols Used on Labeling

Symbol	Description	Standard
R_x ONLY	Prescription use only	FDA 21 CFR 801.109
EC REP	Authorized Representative in the European Community	ISO 15223-1, Reference 5.1.2
CE XXXX	CE Mark with notified body reference number	MDR Regulation (EU) 2017/745
Translations in Box	Translations in Box	Hologic
 www.hologic.com/package-inserts	Consult instructions for use	ISO 15223-1, Reference 5.4.3, Hologic

	Follow instructions for use	IEC 60601-1, Reference No. Table D.2, Safety sign 10 (ISO 7010-M002)
	Do not use if package is damaged	ISO 15223-1, Reference 5.2.8
REF	Catalog number	ISO 15223-1, Reference 5.1.6
LOT	Batch code	ISO 15223-1, Reference 5.1.5
QTY	Quantity	Hologic
	Manufacturer	ISO 15223-1, Reference 5.1.1
	Country of Manufacture	ISO 15223-1, Reference 5.1.11
Patents	Patents	Hologic
	Do not re sterilize	ISO 15223-1, Reference 5.2.6
	Do not re-use	ISO 15223-1, Reference 5.4.2
STERILE R	Sterilized using irradiation	ISO 15223-1, Reference 5.2.4
	Use-by Date	ISO 15223-1, Reference 5.1.4
	Caution	ISO 15223-1, Reference 5.4.4
	Not safe for magnetic resonance imaging	ASTM F2503 Reference no. Table 2, Symbol 7.3.3;7.4.9.1; Fig. 9
MD	Medical device	ISO 15223-1, Reference 5.7.7
	Single sterile barrier system with protective packaging outside	ISO 15223-1, Reference 5.2.14
	Single sterile barrier system	ISO 15223-1, Reference 5.2.11
YYYY-MM-DD	Date format: YYYY represents the year MM represents the month DD represents the day	Hologic
	Country code for translation	ISO 3166

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© Copyright Hologic 2024. Todos os direitos reservados. Hologic, ATEC e logotipos associados são marcas comerciais ou marcas registradas da Hologic, Inc. e/ou de suas subsidiárias nos Estados Unidos e/ou em outros países.