Stereotactic Guided Breast Biopsy System

Instructions For Use

MAN-08967 Revision 001

GR Σύστημα στερεοτακτικής καθοδηγούμενης βιοψίας μαστού, Οδηγίες χρήσης ΜΑΝ-08967, Αναθεώρηση 001 NL Stereotactisch geleid borstbiopsiesysteem, Gebruiksaanwijzing MAN-08967, Revisie 001 SV Stereotaktiskt styrt bröstbiopsisystem, Bruksanvisning MAN-08967, Revidering 001 FL Stereotaktinen ohjattu rintabiopsiajärjestelmä, Käyttöohjeet MAN-08967, Tarkistettu versio 001 NO Stereotaktisk styrt brystbiopsisystem, Bruksanvisning MAN-08967, Oppdatering 001 Stereotaktisk guidet system til brystbiopsi. Brugerveiledning MAN-08967. Revision 001 DA FC Système de biopsie mammaire quidée par stéréotaxie, Mode d'emploi MAN-08967, Révision 001 PB Sistema de Biópsia Mamária Guiada Estereotáxica, Instruções de uso MAN-08967, Revisão 001 ΕI Sistema de biopsia de mama con guiado estereotáctico, Instrucciones de uso MAN-08967, Revisión 001 CS Stereotakticky naváděný systém pro biopsii prsu, Návod k použití MAN-08967, Revize 001 ΗU Sztereotaxiás irányítású emlőbiopsziás rendszer, Használati utasítás MAN-08967, Átdolgozás 001 RO Sistem de biopsie mamară ghidată stereotactic, Instrucțiuni de utilizare MAN-08967, Revizuirea 001 SK Systém na stereotakticky navigovanú biopsiu prsníka, Návod na použitie MAN-08967, Revízia 001 Stereotaktično voden sistem za biopsijo dojke, Navodila za uporabo MAN-08967, Revizija 001 SL PL System do stereotaktycznej biopsji piersi, Instrukcja obsługi MAN-08967, Wersja 001 TR Stereotaktik Rehberli Meme Biyopsi Sistemi, Kullanım Talimatları MAN-08967, Revizyon 001 SR Стереотаксички навођен систем за биопсију дојке, Упутство за употребу МАМ-08967, Ревизија 001 PT Sistema de Biopsia Mamária Guiada por Estereotaxia, Instruções de utilização MAN-08967, Revisão 001 IT Sistema di biopsia mammaria a guida stereotassica, Istruzioni per l'uso MAN-08967, Revisione 001 DE Stereotaktisch geführtes Brustbiopsiesystem, Gebrauchsanweisung MAN-08967, Revision 001 FR Système de biopsie mammaire sous guidage stéréotaxique, Mode d'emploi MAN-08967, Révision 001

HOLOGIC

Eviva[®] Breast Biopsy Device Instructions for Use (IFU)

Please read all information carefully. Failure to properly follow the instructions may lead to unintended consequences.

Important: This package insert is designed to provide instructions for clinical use (IFU) for the Eviva breast biopsy device to be used with the ATEC[®] breast biopsy and excision system console. It is not a reference to surgical techniques.

Patient Target Group

The target patient population of EVIVA Breast Biopsy and Excision System includes patients undergoing breast tissue removal for diagnostic sampling of breast abnormalities.

Indications

The Eviva breast biopsy device is indicated to provide breast tissue samples for diagnostic sampling of breast abnormalities under stereotactic breast biopsy guidance. The Eviva device is intended to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality. The extent of histologic abnormality cannot be reliably determined from its mammographic appearance. Therefore, the extent of removal of the imaged evidence of an abnormality does not predict the extent of removal of histologic abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedure.

Contraindications

The Eviva device is for diagnostic use only and is NOT indicated for therapeutic use. The Eviva device is contraindicated for those patients who, based on physician's judgment, may be at increased risk or develop complications associated with core removal or biopsy. Patients receiving anticoagulant therapy or who may have bleeding disorders may be at increased risk of procedural complications.

Intended User

The EVIVA Breast Biopsy and Excision System should be used only by physicians trained in open or percutaneous biopsy procedures.

Device Description

The Eviva device is a vacuum-assisted biopsy device which is used to remove breast tissue in a minimally invasive manner using stereotactic imaging. The Eviva device is designed for biopsy and diagnosis only. The device is NOT intended to be used outside the standards, specifications and limitations as outlined in this IFU and other related materials. The device is for single patient use and the entire device is fully disposable. The device primarily consists of a hollow needle with a side aperture and a sharpened inner cannula that rotates at high speeds and extends across the aperture to acquire targeted tissue. Throughout this process, vacuum is created inside the device to assist in pulling tissue into the aperture. Saline is supplied through the device to lavage the cavity and deliver tissue to the tissue filter. An introducer system, specifically designed for use with the Eviva device, is included with the device. This introducer system maintains access to the targeted area of interest and allows for deployment of a biopsy site marker.

Refer to the "ATEC Breast Biopsy and Excision System Operator's Manual" for detailed console instructions, performance characteristic and accessory information.

Expected Clinical Benefit

The EVIVA Breast Biopsy and Excision System enables physicians to use a minimally invasive system to extract small biopsy samples of potentially malignant breast tissue.

Device Preparation – Firing or Non-firing Device (all steps to be performed per standard interventional technique)

1. Assemble canister, place in the canister holder on the top of the console and connect vacuum-line assembly from console to the top port of canister lid. Ensure canister lid and large port cap are secure and sealed to avoid vacuum leaks.

WARNING: Before use, inspect the protective packaging and needle to verify that neither has been damaged during shipment. If it appears that the packaging or needle has been compromised, do not use the needle.

2. Place device package on console workspace.

- 3. Open package. Leave protective sheath on the device tip.
- 4. Open the saline bag, remove cap, spike saline bag and place the bag on the hook on the side of the console. Place saline line through pinch valve.
- 5. Connect tubing set to console and canister. Attach red sleeve to red device receptacle on console, black sleeve to black device receptacle on console and blue vacuum into the "Vacuum" port on the canister lid. For more detailed console connection instructions, please refer to the "ATEC Breast Biopsy and Excision System Operator's Manual."
- 6. Power console on (console will default to "Biopsy" mode).
- 7. Ensure device is in post-fire position by confirming the green button on the remote valve is fully depressed (firing device only).
- 8. Select "Setup" to prime device (saline visible in aperture) and confirm adequate vacuum is achieved ("Vacuum Ready" light should be illuminated green). If ample vacuum is not achieved, the "Vacuum Ready" light will flash. If the light flashes, refer to the Troubleshooting section of this IFU and/or the Troubleshooting section of the "ATEC Breast Biopsy and Excision System Operator's Manual."
- 9. Select "Test" and ensure one diagnostic biopsy cycle is completed (console automatically switches to "Biopsy" after passing the "Test" cycle). If the "Return to Setup" or "Retest Handpiece" lights illuminate, then refer to the Troubleshooting section of this IFU and/or the "ATEC Breast Biopsy and Excision System Operator's Manual."
- 10. Depress the white button on the remote valve to place the device in the pre-fire position. Refer to Figure C. Then, depress the green button on the remote valve to ensure that the outer cannula travels to the post-fire position.
- 11. Attach adapter and needle guide to the stage of the imaging system. Some imaging systems require input dimensions for the Eviva device. In these cases, refer to the documentation included with the adapter for your system.

Performing a Biopsy (all steps to be performed per standard intervential technique)

- 1. Identify targeted area.
- 2. Prepare the site as required.

Firing Device (20mm or 12mm aperture)

- a. Carefully remove the device from tray.
- b. Remove protective sheath from the device.
- c. Remove introducer from tray and carefully assemble introducer onto device, (refer to Figure A).
- d. Confirm the console is in "Biopsy" mode, prior to manipulation of the remote valve. Note: If console is not in the "Biopsy" mode, then the device may not function appropriately when attempting to place it in the "pre-fired" position.
- e. Set the device to the pre-fire position by depressing the white button on the remote valve, (refer to Figure C). Note: When the Eviva device is in the pre-fire position, the ATEC console must be powered on and in "Biopsy" mode to prevent migration of device to the post-fire biopsy position.
- f. Slide device/introducer onto the adapter until it locks into place with an audible click when the device is securely engaged it will no longer slide on the adapter.
- g. Advance the device to the desired pre-fire coordinates provided by the imaging system.

h. Fire the outer cannula by fully depressing the green button of the remote valve. This will advance the outer cannula by the stroke noted on the product label.

Non-firing Device (12mm aperture)

- a. Remove disposable stylet from tray and remove protective sheath from the stylet.
- b. Place stylet into the stylet bracket, (refer to Figure B).
- c. Remove introducer from tray and carefully assemble introducer onto stylet.
- d. Slide stylet bracket assembly onto the adapter. Advance the stylet bracket forward until the bracket and introducer lock into place with an audible click.
- e. Advance to the desired coordinates provided by the imaging system.
- f. While depressing the button on the back of the stylet bracket, remove the stylet bracket assembly from the patient and slide off the adapter. Introducer will remain connected to adapter and in the patient.
- g. Confirm the console is in "Biopsy" mode.
- h. Slide the device onto the adapter while guiding the blunt needle into the introducer until it locks into place with an audible click. When the device is securely engaged it will no longer slide on the adapter.
- 3. To administer preferred anesthetic, attach 10cc needleless syringe of anesthetic to Y-valve. To initiate automatic delivery of anesthetic, inject 1-2cc of anesthetic manually. Note: The device will not allow Y-valve injections when in the pre-fire position. Attempting to inject in the pre-fire position may result in device malfunction.
- 4. To begin tissue acquisition, press down on the footswitch and hold in the depressed position throughout the biopsy procedure. Removing your foot from the footswitch will inactivate or stop the device after completing the current cycle.
- 5. To rotate the aperture on the outer cannula of the device, turn the rear thumbwheel (refer to Figure A). Turn the thumbwheel from one position to the next desired position as the console beeps until the desired target area has been sampled. The numbers in the window on the device and the line on the thumbwheel indicate position of the aperture.
- 6. Select "Lavage" mode to irrigate and aspirate the cavity and clear device of tissue. (It is not necessary to depress the footswitch.)
- 7. To vacuum cavity, disconnect luer lock of Y-valve to open saline line to atmosphere resulting in constant aspiration to the biopsy cavity.
- 8. Disconnect filter chamber to remove tissue cores from tissue filter. Note: If blood loss from patient is noted, then the device can be placed in the biopsy mode. This closes the aperture and prevents any fluid loss from the biopsy cavity.
- 9. Place console in "Biopsy" mode in order to close the aperture prior to removal of device from biopsy cavity.
- 10. Remove device from adapter by lifting rear latch and sliding the device back and free of the introducer. Introducer will remain connected to adapter and in the patient.
- 11. Power console off.
- 12. For deployment of the biopsy site marker, refer to the appropriate Instructions for Use (IFU).
- 13. Remove introducer from adapter by pinching tabs and sliding back and free of the adapter.
- 14. Appropriately dispose of all used and unused biopsy items.
- 15. Properly clean and disinfect adapter as necessary.

Warnings and Precautions

- Observe caution when attaching sharp devices to the adapter.
- When the Eviva device is in the pre-fire position, the ATEC console must be powered on and in "Biopsy" mode to prevent migration of device to the post-fire biopsy position. The Eviva device must be in the post-fire position, to perform "Setup", "Test", "Biopsy", "Lavage" or "Manual Aspiration".
- Do not cut tubing between device and remote valve in order to avoid inadvertent firing. This is indicated by the following symbol:

- As with any medical procedure, please ensure that users wear appropriate personal protective equipment to guard against potential contact with bodily fluids.
- The Eviva device is not recommended for use with MRI and ultrasound.

- The Eviva device biopsy procedure should be performed only by persons having adequate training and familiarity with this procedure. Consult medical literature relative to techniques, complications and hazards prior to performance of any minimally invasive procedure.
- The Eviva device should be used only by physicians trained in percutaneous biopsy procedures. The Eviva breast biopsy device is not recommended for use with MRI.
- RONLY Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
- Sound professional judgment should be used when the Eviva device is used on patients with breast implants.
- Avoid operator or instrument contact with the sheathed needle portion of the Eviva device.
- Minimally invasive instruments and accessories manufactured or distributed by companies not authorized by Hologic, Inc. may not be compatible with the Eviva device. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all opened instruments whether used or unused.
- Do not resterilize and/or reuse the Eviva device and/or introducer. Resterilization and/or reuse may compromise the
 integrity of the instrument. This may lead to potential risks of failure of the device to perform as intended and/or crosscontamination associated with using inadequately cleaned and sterilized devices.
- A complete and comprehensive preoperative medical history and physical examination are suggested. Radiographic evaluation and laboratory tests may be included.

How Supplied

The Eviva device is supplied sterile for single patient use. Discard into an appropriate container after use. The biopsy device is sterilized by radiation to sterility assurance level (SAL) 10⁻⁶.

Product Complaints and Technical Support

Report any complaints or problems in the quality, reliability, safety, or performance of this product to Hologic. If the device has caused or added to patient injury, immediately report the incident to Hologic Authorized Representative and Competent authority of the respective member state or country. The Competent Authorities, for medical devices, are usually the individual Member States' Ministry of Health, or an agency within the Ministry of Health.

For More Information

For technical support or reorder information in the United States, please contact:

Hologic, Inc. 250 Campus Drive Marlborough, MA 01752 USA Phone: 877-371-4372 BreastHealth.Support@hologic.com

International customers, contact your distributor or local Hologic Sales Representative:

European Representative Hologic BV Da Vincilaan 5 1930 Zaventem Belgium Tel: +32 2 711 46 80

Troubleshooting

No Gap	No Gap	Gap
Figure D Fully Fired Position	Figure E Fully Armed Position	Figure F Partially Fired Position

Remote Valve Position	Fully Fired	Fully Armed	Partially Fired
Needle Position	Needle in fired position. Aperture will open and close as intended.	Needle in armed position. Aperture will remain closed in all modes.	Needle in fired position. Aperture may not fully close during biopsy.
Symptom(s)	Device functions as intended.	No saline flow in all modes. Device sounds normal during biopsy but unable to acquire tissue.	Saline flow is normal in all modes. Device sounds normal during biopsy but acquires inadequate or no tissue.*
Remedy	Not Applicable	Fully fire the needle. See Figure D.	Fully fire the needle. See Figure D.

 * May also lead to failures during setup and test.

Failure Mode	Possible Cause	Possible Remedy
Device will not assume the pre-fire or post- fire position.	Console is not in "Biopsy" mode.	Press "Biopsy" on console. If device will still not pre-fire/post- fire, replace the device and repeat the setup procedure.
During "Setup," "Test," "Biopsy," "Lavage" or "Manual Aspiration" fluid and/or tissue is not being pulled to the tissue filter or aperture is not opening or closing as expected.	Remote valve is not in full post- fired position.	Place device in post-fire position by fully depressing green button on remote valve.
Unable to depress syringe to administer anesthetic.	Device is not in post-fire position.	Place device in post-fire position by depressing green button on remote valve.

Symbols Used on Labeling

Symbol	Description	Standard	
	Prescription use only	21 CFR 801.109	
EC REP	Authorized Representative in the European Community	ISO 15223-1, Reference 5.1.2	
\triangle	Caution	ISO 15223-1, Reference 5.4.4	
CE xxxx	CE Mark with notified body reference number	MDR Regulation (EU) 2017/745	

Symbol	Description	Standard
()	CE Mark	MDR Regulation (EU) 2017/745
	European Conformity	
[]]	Consult instructions for use	ISO 15223-1, Reference 5.4.3
www.hologic.com/package-inserts		
	Follow instrctions for use	IEC 60601-1, Reference No. Table D.2,
		Salety sign 10 (150 7010-10002)
Translations in Day	Translations in Box	Hologic
	Do not use if package is damaged	ISO 15223-1, Reference 5.2.8
	Catalog number	ISO 15223-1, Reference 5.1.6
REF	_	
	Batch code	ISO 15223-1, Reference 5.1.5
ΟΤΥ	Quantity	Hologic
	Manufacturer	ISO 15223-1, Reference 5.1.1
	Country of Manufacture	ISO 15223-1, Reference 5.1.11
	Detente	
Patents	Patents	HOIOGIC
	Do not re sterilize	ISO 15223-1, Reference 5.2.6
STERNIZE		
	Do not re-use	ISO 15223-1 Reference 5.4.2
$ (\dot{\chi})$		
(MR)	MR Unsafe	ASTM F2503 Reference no. Table 2, Symbol 7 3 3 7 4 9 1, Fig. 9
		Symbol 7.3.3;7.4.9.1; Fig. 9
STERILE R	Sterilized using irradiation	ISO 15223-1, Reference 5.2.4
	Lice by Date	ISO 15222 1 Deference 5 1 4
	Use-by Date	130 13223-1, Reference 5.1.4
MD	Medical device	ISO 15223-1, Reference 5.7.7
$\langle \bigcirc \rangle$	Single sterile barrier system with protective packaging	ISO 15223-1, Reference 5.2.14
	outside	
\bigcirc	Single sterile barrier system	ISO 15223-1:, Reference 5.2.11
	Expiration date is represented by the following:	Hologic
	YYYY represents the year	
	MM represents the month	
	DD represents the day	
	ישט וואס אישט אין אי	
	Do not cut this line	N/A
► >		
	Country code for translation	ISO 3166

© 2024 Hologic, Inc. Hologic, The Science of Sure, Eviva and associated logos are trademarks and/or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries.

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Hologic, Inc. 250 Campus Drive Marlborough, MA 01752 États-Unis Tél. : 877-371-4372 BreastHealth.Support@hologic.com
Hologic, Inc. 250 Campus Drive, Marlborough, MA 01752 USA 1-877-371-4372