TriMark[®] for eviva[®]

Biopsy Site Identification System
Titanium Biopsy Site Marker

Instructions For Use
Gebruiksaanwijzing
Mode d'emploi
Bedienungsanleitung
Istruzioni per l'uso
Instruções de Utilização
Instrucciones de uso

HOLOGIC[®]

TriMark® for Eviva® Biopsy Site Marker System

Instructions for Use

Please read all information carefully. Failure to properly follow the instructions may lead to unanticipated surgical consequences.

Important: This package insert is designed to provide Instructions for Use of the TriMark® biopsy site marker system for the Eviva® biopsy device system. It is not a reference to surgical techniques.

Indications


The TriMark for Eviva biopsy site marker system is indicated for use to mark an open or percutaneous biopsy site to radiographically mark the location of the biopsy site.

Contraindications

None known.

Device Description

The TriMark for Eviva biopsy site marker system is supplied as a sterile, single patient use system, comprised of a singular implant-grade titanium marker. The deployment device is a hand held device that delivers the marker from the distal tip. The deployment device consists of a cannula, handle, rigid push rod and plunger.

Device Preparation and Use

1. Prior to use of the TriMark for Eviva biopsy site marker system, inspect the protective packaging and device to verify that neither has been damaged during shipment. If it appears that the packaging has been compromised, do not use the device.
2. Remove the TriMark for Eviva biopsy site marker system from its protective packaging using sterile technique.
3. Remove the biopsy device from the introducer sheath.
4. Place the TriMark for Eviva biopsy site marker system through the hub of the introducer sheath.
5. Advance the TriMark for Eviva biopsy site marker system until the handle snaps to the hub.
6. Deploy the TriMark for Eviva biopsy site marker by advancing the deployment plunger all the way forward until it latches onto the handle.
7. Slowly remove the deployment device and introducer sheath as one unit from the breast and properly dispose.

Warnings and Precautions

- There are possible adverse reactions when an object is implanted in the body. It is the responsibility of the physician to evaluate any risk or benefit prior to the use of this device.
- Potential complications of marker clip placement consist of pain, seroma formation, inflammation, bruising, hematoma, hemorrhage, infections, hypersensitivity or allergic reaction, soft tissue damage, misdiagnosis (due to marker clip migration), perforation or scar tissue.
- The TriMark for Eviva biopsy site marker system is not recommended for use within the bore of an MRI magnet.
- The TriMark for Eviva biopsy site marker system is not intended for use in patients with breast implants.
- The TriMark for Eviva biopsy site marker system should be used only by persons having adequate training and familiarity with this procedure. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.

- This TriMark for Eviva biopsy site marker system should be used only by physicians trained in open or percutaneous biopsy procedures.
- **ONLY** Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
- The TriMark for Eviva biopsy site marker should be deployed into the cavity created during the biopsy procedure. Deployment into tissue outside of the biopsy cavity is not recommended.
- Marker position relative to established landmarks may change under mammography upon subsequent breast compressions.
- The TriMark for Eviva biopsy site marker is not intended to be repositioned or recaptured after it has been deployed.
- Excess hematoma within the biopsy device can lead to marker adhesion, increasing the risk of marker drag out.
- Care should be taken to avoid damaging the cannula. Avoid operator or instrument contact with the TriMark for Eviva biopsy site marker system or the distal end of the device. Contact with the distal end may result in loss of sterility.
- The implanted TriMark for Eviva biopsy site marker is magnetic resonance imaging (MRI) conditional. The implanted TriMark for Eviva biopsy site marker presents no additional risk to patient or operator from magnetic forces, torque, heating, induced voltage or movement but may affect MRI image quality.
- Minimally invasive instruments and accessories manufactured or distributed by companies may not be compatible with the TriMark for Eviva biopsy site marker system. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all opened instruments whether used or unused.
- Do not resterilize and/or reuse the TriMark for Eviva biopsy site marker system. Resterilization and/or reuse may compromise the integrity of the instrument. This may lead to potential risks of failure of the device to perform as intended and/or cross-contamination associated with inadequately cleaned and sterilized devices.

How Supplied

The TriMark for Eviva biopsy site marker system is sterilized by radiation and supplied preloaded for single patient use. Discard into an appropriate container after use.

As Identified on Labels:

QTY

Number of Devices Enclosed.

YYYY-MM-DD Expiration date is represented by the following:

YYYY represents the year

MM represents the month

DD represents the day

For More Information

For technical support or reorder information in the United States, please contact:

Hologic, Inc.
250 Campus Drive
Marlborough, MA 01752 USA
Phone: 877-371-4372
BreastHealth.Support@hologic.com

International customers, contact your distributor or local Hologic Sales Representative:

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Additional Symbols Used on Labeling

Symbol	Description	Standard
	Prescription use only	21 CFR 801.109
	Authorized Representative in the European Community	ISO 15223-1:2016, Reference 5.1.2
	CE Mark with notified body reference number	MDR Regulation (EU) 2017/745
 www.hologic.com/package-inserts	Consult electronic instructions for use	ISO 15223-1:2016, Reference 5.4.3
	Translations in Box	Hologic
	Do not use if package is damaged	ISO 15223-1:2016, Reference 5.2.8
	Catalog number	ISO 15223-1:2016, Reference 5.1.6
	Batch code	ISO 15223-1:2016, Reference 5.1.5

	Quantity	Hologic
	Manufacturer	ISO 15223-1:2016, Reference 5.1.1
	Patents	Hologic
	Do not re sterilize	ISO 15223-1:2016, Reference 5.2.6
	Do not re-use	ISO 15223-1:2016, Reference 5.4.2
	MR Conditional	ASTM F2503 Reference no. Table 2; 7.4.6.1; Fig 6,7
	Sterilized using irradiation	ISO 15223-1:2016, Reference 5.2.4
	Use-by Date	ISO 15223-1:2016, Reference 5.1.4
	Made in Costa Rica	Hologic
YYYY-MM-DD	Expiration date is represented by the following: YYYY represents the year MM represents the month DD represents the day	Hologic

www.hologic.com/package-inserts

© 2010-2022 Hologic, Inc. All rights reserved. Hologic, Eviva and TriMark are registered trademarks and/or trademarks of Hologic, Inc. and/or its subsidiaries in the United States and other countries.

MAN-03485-4370 Revision 007
July 2022

TriMark® for Eviva®

Biopsy Site Identification System

Product Code: TRIMARK-EVIVA-10, TRIMARK-EVIVA-2S-10,
TRIMARK-EVIVA- 13, TRIMARK-EVIVA-2S-13

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