## JustRight<sup>®</sup> 5mm Stapler

**REF** JR-ST25 -2.0

5 mm shaft diameter, 20 cm shaft length, 25 mm jaw length, 2.0 mm staple

## JustRight<sup>®</sup> 5mm Reload

**REF** JR-REL25-2.0

5 mm shaft diameter, 20 cm shaft length, 25 mm jaw length, 2.0 mm staple

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**sv – Bruksanvisning**

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






## Caution

Read all warnings, cautions, and instructions provided with this instrument before use.

## Symbols

|                                                                                   |                                                                                            |                                                                                   |                                                                             |
|-----------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------|
|  | Sterilized Using Ethylene Oxide                                                            |  | Do Not Use if Package is Opened or Damaged and Consult Instructions for Use |
|  | Catalog, Reorder or Reference Number                                                       |  | Caution                                                                     |
|  | Consult Instructions for Use                                                               |  | Store in Temperatures Between 18°C – 28°C                                   |
|  | Do Not Re-use                                                                              |  | Transport Conditions Must Not Exceed: -18°C - 60°C                          |
|  | Manufacturer                                                                               |  | Not Made with Natural Rubber Latex                                          |
|  | Do Not Resterilize                                                                         |  | Medical Device                                                              |
|  | Use-by Date                                                                                |  | Keep Dry                                                                    |
|  | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician |  | Lot Number                                                                  |
|  | Date of Manufacture                                                                        |  | European Community Authorized Representative                                |
|  | This package forms the sterile barrier                                                     |  | Plastics - Other                                                            |

## Technical Characteristics for the JustRight 5.0mm Stapler and Reload

|                                                    |                                                                                     |                                        |                                                                                     |
|----------------------------------------------------|-------------------------------------------------------------------------------------|----------------------------------------|-------------------------------------------------------------------------------------|
| Size of Staple: 2.0 mm                             |    | Length of Shaft: 20 cm                 |    |
|                                                    |                                                                                     | Distal Jaw Opening: ≤15°               |   |
| Compatible Trocar Sizes: ≥5 mm                     |   | Compressed Tissue Range: 0.75 - 1.0 mm |  |
|                                                    |                                                                                     |                                        |                                                                                     |
| Staple Pattern: Staggered Rows per Cartridge: 4    |  |                                        |                                                                                     |
| Indicated Firing Zone Length of Staple Line: 25 mm |  |                                        |                                                                                     |
| Maximum Number of Reloads per Stapler Handle: 15   |  |                                        |                                                                                     |
| Pre-fire Compression Time: ≥10 seconds             |  |                                        |                                                                                     |

## Description

The JustRight 5mm Stapler is a sterile, single patient use device that staples and simultaneously divides tissue. The JustRight 5mm Stapler fitted with a JustRight 5mm Reload places four 25 mm staggered rows of titanium staples while simultaneously dividing the tissue between the two inner staple rows. The unformed staple leg length is 2 mm and the staple back span dimension is 3 mm. The JustRight 5mm Stapler is provided with a stapler handle and one JustRight 5mm Reload staple cartridge. The JustRight 5mm Reload cartridges are also provided separately to allow for multiple staple line applications within one surgical procedure. One JustRight 5mm Stapler handle can be fired up to 15 times in a single procedure using 15 JustRight 5mm Reload cartridges. The JustRight 5mm Stapler is compatible for introduction and use through endoscopic cannulas as small as 5 mm. The stapler handle has a rotation knob to allow for 360° rotation of the JustRight 5mm Reload cartridge. The JustRight 5mm Stapler is used on tissue that compresses to a thickness of 0.75 to 1.0 mm.

Titanium, as found in JustRight Surgical surgical staples, is acceptable for use with MRI and NMR procedures as it is nonferromagnetic. The implanted staples are not subject to displacement or induced torque, and do not result in localized heating. The implantable staples are categorized as MR Safe.

## Indications for Use

The JustRight 5mm Stapler is intended for use in abdominal, gynecologic, pediatric, and thoracic surgery for resection, transection, and creation of anastomosis.

## Contraindications

- Do not use the JustRight 5mm Stapler 2.0 mm staples on any tissue that compresses to less than 0.75 mm in thickness, or on any tissue that cannot comfortably compress to 1.0 mm.
- Do not use the JustRight 5mm Stapler on tissue such as liver, hepatic vasculature, spleen, pancreas, or biliary structures where compressibility is such that closure of the instrument would be destructive.
- Do not use the JustRight 5mm Stapler where adequacy of hemostasis cannot be verified visually after application.
- Do not use the JustRight 5mm Stapler with a staple line reinforcement material. Adequate staple formation cannot be verified when using staple line buttressing materials such as PERI-STRIPS or SEAMGUARD products.
- Do not use the JustRight 5mm Stapler on the aorta or on tissues that are necrotic, friable, or have altered integrity.

## Warnings and Cautions

### Warning

This product is intended for single use only; it cannot be adequately cleaned or resterilized for safe reuse. Attempts to clean or sterilize may result in biocompatibility infection or product failure risks to the patient.

Place tissue in the center of the jaws. To avoid incomplete ligation, do not grasp structure beyond the staple line surface; do not place tissue in the jaw hinge. Any tissue extending beyond the cut mark will not be transected.

When using the stapler more than once during a SINGLE surgical procedure, be sure to remove the empty JustRight 5mm Reload single use unit and reload a new one. When the JustRight 5mm Reload is empty, the blue staple deployment and cut lever will not engage. This prevents inadvertent actuation of the integrated knife.

After staples have been deployed, do not reclamp an empty staple cartridge on tissue or use to manipulate structures. Inadvertent tissue damage may occur.

Preoperative radiotherapy may result in changes to tissue. These changes may, for example, cause the tissue thickness to exceed the indicated range for the staple size. Careful consideration with regard to staple size should be given to any pre-surgical treatment the patient may have undergone.

Device performance testing was conducted preclinically in animals. No human clinical testing was performed.

Always inspect the tissue thickness prior to application of the JustRight 5mm Stapler.

Endoscopic procedures should be performed only by physicians having adequate training and familiarity with endoscopic techniques. Prior to performing any endoscopic procedures, consult the medical literature relative to technique, complications, and hazards.

When using endoscopic instruments and accessories from different manufacturers in a procedure, be sure to verify compatibility and ensure that electrical isolation or grounding is not compromised.

Do not use the JustRight 5mm Stapler unless properly trained.

### Caution

After firing, always inspect the staple line for hemostasis. Minor bleeding or oozing may be controlled by electrocautery or manual sutures.


The anvil of the stapler must be completely visible (past the cannula sleeve) prior to opening the jaws within the body cavity.

|                                                                                                                              |
|------------------------------------------------------------------------------------------------------------------------------|
| The JustRight 5mm Stapler is intended to be used with a laparoscopic cannula with a minimum diameter of 5.6 mm.              |
| Always close the blue Top Lever on the stapler handle prior to introducing and removing the stapler from the cannula sleeve. |
| If packaging is exposed to environmental conditions outside of those specified, do not use.                                  |
| Inspect packaging for damage. If damaged, do not use.                                                                        |
| Avoid placing fingers between the lever and handle, or in the jaws when clamping. Injury to the user may result.             |
| This device is intended for use by medical professionals only.                                                               |


## Preparing the Stapler

### Handle - JustRight 5mm Stapler

1. To insert the staple reload, the top blue lever must be in the closed position.

### Reload - JustRight 5mm Reload

2. Remove the green ship lock from the JR-REL25-2.0 staple reload by depressing the ribbed grip area on the reload while rotating the green ship lock in the direction of the arrow and pulling on it.

**NOTE:** Do not remove the green ship lock on the Reload until you are ready to load into the handle.

3. Align the blue dot on the reload connector with the blue dot on the handle's rotation knob and insert until fully seated.


**NOTE:** The reload can only be inserted one way; the tabs on the reload connector are different widths and match the slots in the rotation knob.

4. Rotate the reload shaft clockwise until the blue dot on the outside of the reload shaft aligns with the blue dot on the rotation knob.


## Using the Stapler

1. Insert the instrument into the cannula sleeve.

- NOTE:** The top blue lever (1) must be in the closed position before attempting to insert the stapler into cannula sleeve.
2. Lift the top blue lever (1) to open the jaws

3. Turn the rotation knob until the jaws are in the required position.


### Warning

Establish proximal control of blood vessels prior to stapling where practical and implement methods of blood vessel control in the event of stapler failure

When positioning the stapler on the application site, ensure that no obstruction, such as clips, are incorporated into the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.

Take care when positioning the stapler on the application site to avoid the unintended stapling of other anatomic structures.

4. Position the tissue to be transected between the black lines on the outside of the instrument jaws.


5. Pull back on the Jaw Clamp Lever (2) until the jaws lock on tissue to be transected.

### Warning

If the jaws of the device do not lock closed, the target tissue is outside the compression specifications. Reposition the jaws on to tissue within compression specification range.

If the tissue does not compress to 0.75 to 1.0 mm in thickness it is outside the compression specifications. Do not staple or cut on this tissue. Reposition the jaws on to tissue within the compression specification range.

6. To release and reposition jaws, depress the blue Jaw Release button (5) on the back of the handle.

7. Repeat Step 5.

**NOTE:** After closing, pausing for 10 seconds prior to firing may result in better tissue compression and better staple formation.

8. Flip the green Safety Lock (3) down. Do not hold the safety lock down.

Safety Lock (green) - Up

Safety Lock (green) - Down

9. Pull back completely on the blue Staple Deployment and Cut Lever (4) and release to reset the tissue structure.

**NOTE:** Only pull back on the blue Staple Deployment and Cut Lever (4) once to completely fire a staple line.

10. Press the blue Jaw Release button (5) to open the jaws.

11. Push down on the blue Top Lever (1) (closed position) and remove the device from the cannula.

12. Remove and set aside a spent staple Reload after firing.

13. Reload the device if necessary:

- Remove the spent reload by depressing the ribbed grip area on the reload and rotating it until the alignment pins are disengaged, then pull the spent reload from the handle and discard it.
- Remove the green ship lock from a new JR-REL25-2.0 staple reload by depressing the ribbed grip area on the reload while rotating the green ship lock in the direction of the arrow and pulling on it.
- Insert a new reload by aligning the blue dot on the reload connector with the blue dot on the handle's rotation knob and inserting until fully seated.
- Rotate the reload shaft clockwise until the blue dot on the outside of the reload shaft aligns with the blue dot on the rotation knob.

14. The JustRight 5mm Stapler may be reloaded and fired up to 15 times in a single procedure.

## Warranty

Except as otherwise expressly stated in an agreement between Hologic and its original customer ("Customer"), Hologic equipment ("Equipment") is warranted to the original Customer to perform substantially in accordance with published product specifications for one (1) year starting from the date of shipment ("Warranty Period"). Replacement parts and remanufactured items are warranted for the remainder of the Warranty Period or 90 days from shipment, whichever is longer. Consumable supplies are warranted to conform to published specifications for a period ending on the expiration date shown on their respective packages. Licensed software is warranted to operate in accordance with published specifications. Services are warranted to be supplied in a workman-like manner. Hologic does not warrant that use of Products will be uninterrupted or error-free, or that Products will operate with non-Hologic authorized third-party products.

These warranties do not apply to any item that is: (a) repaired, moved, or altered other than by Hologic authorized service personnel; (b) subjected to physical (including thermal or electrical) abuse, stress, or misuse; (c) stored, maintained, or operated in any manner inconsistent with applicable Hologic specifications or instructions, including Customer's refusal to allow Hologic recommended Software upgrades; or (d) designated as supplied subject to a non-Hologic warranty or on a pre-release or "as-is" basis.

---

## Disposal

The used disposable device must be treated as biohazardous waste and disposed of according to standard practices of the hospital or clinic where the treatment is performed.

Consult your local regulations for disposal/electronics recycling. Do not place into a municipal waste system unless authorized to do so by local authorities.

## Technical Support and Product Return Information

Contact Hologic Technical Support for proper disposal if any part of the JustRight Stapler fails to operate as intended. If product is to be returned to Hologic for any reason, Technical Support will issue a Returned Materials Authorization (RMA) number and biohazard kit if applicable. Return the JustRight Stapler according to the instructions provided by Technical Support. Be sure to clean and sterilize the product before returning it and include all accessories in the box with the returned unit. Return used or opened product according to the instructions provided with the Hologic-supplied biohazard kit.

## Complaints and Serious Incident Reporting

Report any complaints or problems in the quality, reliability, safety, or performance of this product to Hologic. If the device has caused or added to patient injury, immediately report the incident to Hologic Authorized Representative and Competent authority of the respective member state or country. The Competent Authorities for medical devices are usually the individual Member States' Ministry of Health, or an agency within the Ministry of Health.

## Technical Support Contact information

Hologic, Inc.  
250 Campus Drive  
Marlborough, MA 01752 USA  
Phone: 1.800.442.9892

## After Surgery

### Warning

Do not reuse or resterilize instruments and accessories labeled "disposable" or "single use only" and dispose according to the facility's policy for biohazards and sharps.

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All rights reserved.