## CoolSeal® Trinity

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# CoolSeal® Trinity

A Maryland Laparoscopic Sealer, Divider, and Dissector

**REF** CSL-TR105-30, Trinity Sealer/Divider/Dissector, 5 mm shaft diameter, 30 cm shaft length

**REF** CSL-TR105-37, Trinity Sealer/Divider/Dissector, 5 mm shaft diameter, 37 cm shaft length

**REF** CSL-TR105-44, Trinity Sealer/Divider/Dissector, 5 mm shaft diameter, 44 cm shaft length

## Compatible Generators:

**REF** CSL-200-50, CoolSeal Generator SW v1.0.0 or later

**REF** 374897, Da Vinci E-200 Generator from Intuitive, with Trinity compatible software

## Caution

Read all warnings, cautions, and instructions provided with this instrument before use.

Read the warnings, cautions, and instructions provided with the compatible generator before system use. Specific warnings, cautions, and instructions for the generator use are not included in this manual.

Federal (USA) law restricts this device to sale by or on the order of a physician.

This device is intended for medical professional use only.

## Symbols

|                                                                                     |                                          |                                                                                     |                                                                                            |
|-------------------------------------------------------------------------------------|------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|
|    | Sterilized Using Ethylene Oxide          |    | Do Not Use if Package is Opened or Damaged                                                 |
| <b>REF</b>                                                                          | Catalog, Reorder or Reference Number     |    | Attention, Consult Accompanying Documents                                                  |
|    | Consult Instructions for Use             | <b>LOT</b>                                                                          | Lot Number                                                                                 |
|    | Single use only                          |    | Not Made with Natural Rubber Latex                                                         |
|    | Manufacturer of Record                   | <b>MD</b>                                                                           | Medical Device                                                                             |
|    | Do Not Resterilize                       |    | Keep Dry                                                                                   |
|   | Use-by Date                              | <b>EC REP</b>                                                                       | Authorized Representative in the European Community                                        |
|  | Date of Manufacture                      |  | This package forms the sterile barrier                                                     |
|  | Store in Temperatures Between 0°C – 28°C | <b>Rx ONLY</b>                                                                      | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician |
|  | Country of Manufacture                   |                                                                                     |                                                                                            |

The CoolSeal Trinity, a Maryland Laparoscopic Sealer, Divider, and Dissector, with a 5 mm diameter shaft is designed for use with the CoolSeal Generator or any generator with the CoolSeal technology. Please refer to the cover page for details on compatible generator models. The Trinity creates seals by application of radiofrequency (RF) electrosurgical energy to vascular structures (vessels and lymphatics) or tissue bundles interposed between its jaws. A blade within the instrument is surgeon-actuated to divide tissue. The double action jaws have been designed to dissect tissue, which includes separating tissue planes and widening openings as necessary for the surgical procedure. Multiple shaft lengths provide additional flexibility for surgical procedures. Maximum rated voltage: 190 V<sub>peak</sub>

## Indications for Use

The CoolSeal Trinity is a bipolar electrosurgical instrument intended for use in minimally invasive or open surgical procedures where ligation and division of vessels, tissue bundles, and lymphatics is desired. The CoolSeal Trinity can be used on vessels (arteries, veins, and vascular bundles) up to and including 7 mm in diameter. It is indicated for use in general surgery procedures including urologic, vascular, and gynecologic. It is indicated for use in adult and pediatric populations (infants, children, and adolescents). Procedures may include, but are not limited to, Nissen fundoplication, colectomy, cholecystectomy, adhesiolysis, hysterectomy, oophorectomy, etc. The CoolSeal Trinity has not been shown to be effective for tubal sterilization or tubal


coagulation for sterilization procedures. Do not use the CoolSeal Trinity for these procedures. The device is contraindicated for use in ENT procedures.

General Warnings

| Warning                                                                                                                                                                                                                                                                                  |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| This product is intended for single use only; it cannot be adequately cleaned or resterilized for safe reuse. Attempts to clean or sterilize may result in bio-incompatibility, infection, or product failure risks to the patient.                                                      |
| These instruments are intended for use only with generators with the CoolSeal technology. Use of these instruments with other generators may not result in the desired tissue effect, may result in injury to patient or surgical team, or may cause damage to the instrument.           |
| Do not use the CoolSeal system unless properly trained. Use of this equipment without proper training may result in serious unintended injury to the patient or surgical team.                                                                                                           |
| Do not use in patients who have electronic implants such as cardiac pacemakers without first consulting a qualified professional (e.g., cardiologist). A possible hazard exists because interference with the action of the electronic implant may occur, or the implant may be damaged. |

| Caution                                                                                                                                                                                                                                                                          |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Use caution during surgical cases in which patients exhibit certain types of vascular pathology (atherosclerosis, aneurismal vessels, etc.). For best results, apply the seal to unaffected vasculature.                                                                         |
| The performance of this single-use device has been tested according to the expected conditions of a single surgical procedure. Subjecting the device to process steps, tools, and/or chemicals commonly used by third-party re-processors may negatively affect its performance. |

Getting Started

1. Remove instrument from tray by firmly pulling on the handle (5). Do not pull on the instrument's jaws (1) or cable (8).
2. Insert the plug (9) into the receptacle on the generator. Follow the instructions in the generator user's guide to complete the setup procedure.

| Warning                                                                                                                                                                                                                                                                                 |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Electric Shock Hazard:</b> Do not connect wet accessories to the generator.                                                                                                                                                                                                          |
| <b>Electric Shock Hazard:</b> Do not wrap the instrument cords around metal objects. This may induce stray currents that could lead to shocks, fires, or injury to the patient or surgical team.                                                                                        |
| Examine all instruments and connections to the system before using. Improper connection may result in arcs, sparks, accessory malfunction, or unintended surgical effects.                                                                                                              |
| Inspect the instrument, instrument cords and generator cable for breaks, cracks, nicks, or other damage before use. Failure to observe this warning may result in injury or electrical shock to the patient or surgical team or cause damage to the instrument. If damaged, do not use. |
| Do not use in the presence of flammable anesthetics or oxidizing gases (such as nitrous oxide (N <sub>2</sub> O) and oxygen) or in close proximity to volatile solvents (such as methanol or alcohol), as explosion may occur.                                                          |

Because of concerns about the carcinogenic and infectious potential of electrosurgical by-products (such as tissue smoke plume and aerosols), protective eye wear, filtration masks, and effective smoke-evacuation equipment should be used in both open and minimally invasive procedures.

#### Caution

Inspect packaging for damage. If damaged, do not use.

### Using the CoolSeal Trinity

#### Warning

Avoid placing fingers between the lever, handle, trigger, or in the jaws. Injury to the user may result.

Place the vessel or vascular bundle in the center of the jaws. To avoid incomplete sealing, do not grasp structure beyond the electrode surface; do not place tissue in the jaw hinge.

Contact between the active instrument electrode and any metal object (hemostats, staples, clips retractors, etc.) may increase current flow and can result in unintended surgical effects such as an effect at an unintended site or insufficient energy deposition.

Do not activate the CoolSeal Trinity until the instrument has fully latched. Activating the generator before this is done may result in improper sealing and may increase thermal spread to tissue outside the intended surgical site.

If the instrument shaft is visibly bent, discard and replace the instrument. A bent shaft may prevent the instrument from functioning properly.

The CoolSeal Trinity is a rigid instrument and should not be inserted through a cannulated endoscope.

For laparoscopic procedures, be alert to these potential hazards:

- Do not use hybrid trocars that are comprised of both metal and plastic components. Capacitive coupling of RF current may cause unintended burns.
- Use appropriately sized trocar to allow for easy insertion and extraction of the instrument.
- Carefully insert and withdraw the instrument through the cannula to avoid damage to the device and/or injury to the patient.
- Close jaws using device lever before insertion/extraction in the trocar.

Activate the CoolSeal Trinity only when the instrument is in direct contact with the target tissue to reduce the possibility of unintended burns.

During a seal cycle, energy is applied to the area between the instrument jaws. This energy may convert water into steam, and this steam may cause unintended injury in close proximity to the jaws. Take care in surgical procedures occurring in confined spaces in anticipation of this possibility.

Keep the cord free from the jaw and latch area of the instrument.

Do not place instruments near or in contact with flammable materials (such as gauze, surgical drapes, or flammable gases.) Instruments that are activated or hot from use may cause a fire. When not using instruments, place them in a clean, dry highly visible area not in contact with the patient. Inadvertent contact with the patient may result in burns.

Avoid the accumulation of naturally occurring flammable gases that may accumulate in body cavities such as the bowel.

Aspirate fluid from the area before activating the instrument. Conductive fluids (e.g., blood or saline) in direct contact with or in close proximity to an active electrode may carry electrical current or heat away from target tissues, which may cause unintended burns to the patient.

#### Notice:

Do not overfill the jaws of the instrument with tissue, as this may reduce device performance.

### Tissue Manipulation and Dissection

#### Warning

Use caution when handling the instrument between uses to avoid accidental activation of the CoolSeal system. Do not place the instrument on the patient or drapes when not in use.

The instrument can be used to manipulate and dissect tissue with the jaws either opened or closed.

#### Knob Rotation

Turn the black rotation knob (3) on the handpiece until the jaws are in the required position.

---

**Notice:**

Do not turn the rotation knob (3) when the lever (6) is latched. Product damage may occur.

**Grasping**

To grasp tissue with the device, place the tissue in the jaws and pull back on the lever.

**Sealing****Warning**

Do not use this instrument on vessels larger than 7 mm in diameter.

Eliminate tension on tissue when sealing and cutting to ensure proper function.

Do not attempt to seal over clips or staples or contact metal objects (e.g., retractors). Contact between an active electrode and any metal objects may result in alternate site burns or incomplete seals.

**Notice:**

Verification testing of the CoolSeal Trinity has been conducted to confirm device performance for up to 140 sealing cycles.

The surgeon may inspect the seal before cutting the vessel or tissue. After inspecting the seal, the surgeon may create a second seal adjacent to the first seal before cutting, as described below.

1. Open the jaws by pushing forward on the movable lever (6).
2. Grasp the intended vessel or vascular bundle in the center of the jaws.
3. Close the movable lever until it latches in place.
4. To activate the instrument, press and hold the blue activation button (4) on the back of the instrument. An activation tone sounds to indicate that the vessel or vascular bundle is being sealed. When the activation cycle is complete, a seal complete tone sounds and RF output ceases.
5. Release the seal activation button on the instrument when the seal cycle is complete and the tone sounds.
6. Open the jaws to release tissue by squeezing the movable lever until it unlatches, then open the jaws by moving the lever forward.
7. To seal adjacent tissue, overlap the edge of the existing seal. The second seal should be distal to the first seal to increase seal margin.

**Cutting****Caution**

Energy based devices such as electrosurgical pencils or ultrasonic scalpels that are associated with thermal spread must not be used to transect seals.

**Notice:**

Do not engage the cutting mechanism over clips, staples, or other metal objects as damage to the cutter may occur.

To activate the cutting mechanism:

1. Grasp the intended tissue in the center of jaws.
2. Close the movable lever until it latches in place.
3. Pull back on the blue cutting trigger (7) to deploy the knife.
4. Release the blue cutting trigger to retract the cutting blade.
5. Open the jaws by squeezing the movable lever until it unlatches, then open the jaws by moving the lever forward.

**Cleaning the Instrument During Use****Warning**

Do not activate the instrument or cutting trigger while cleaning the jaws. Injury to the operating room personnel may result.

Inspect the instrument jaws prior to cleaning to ensure the blade is not deployed.

## Caution

Keep the instrument jaws clean. Buildup of eschar may reduce the seal effectiveness. Wipe the jaw surfaces and edges with a wet gauze pad as needed. Do not clean instrument jaws with a scratch pad or scalpel blade.

Do not use excessive force (torque or bend the instrument jaws) during use or cleaning. Damage to the device may occur. If damaged, do not use.

Wipe jaw surfaces and edges with wet gauze pad as needed.

## Troubleshooting

The following is a list of troubleshooting suggestions for situations encountered when using the instrument with the CoolSeal generator. For details on specific situations, refer to the corresponding generator user's guide. For troubleshooting the da Vinci E-200 generator, refer to the da Vinci E-200 Generator's User Manual.

### Alert Conditions for the CoolSeal Generator:

When an alert condition occurs, energy delivery stops. After the alert condition has been corrected, energy delivery will be immediately available.

#### Seal Complete

##### Indicated By:

- A **three** sequential tone alert
- RF energy delivery stops
- Activation display brightens **blue** for a half (0.5) second

##### Causes:

- Successful vessel seal

##### To Resolve:

- None, normal operation

#### Reactivate / Incomplete Seal

##### Indicated By:

- A **three-pulsed** alert tone
- RF energy delivery stops
- Activation display brightens **amber** for one (1) second

##### Causes:

- Seal time exceeds five (5) seconds OR
- User either opened the instrument jaws or released the activation button, which causes the seal cycle to be interrupted before seal was completed OR
- Current stays at the maximum current limit for longer than four (4) seconds, which indicates an electrical short between the jaws has occurred OR
- The instrument has been activated in open air

##### To Resolve:

1. Release the activation button
2. Press the activation button to reactivate the seal cycle without repositioning the instrument
3. Open the instrument jaws and inspect for a successful seal
4. If possible, reposition the instrument and regrasp tissue in another location, then reactivate the seal cycle
5. Visually inspect seal before cutting

##### Possible use conditions include:

|                                                              |                                                                                                                                                       |
|--------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| Grasping thin tissue or activating in open air               | Open the jaws and confirm that a sufficient amount of tissue is inside the jaws. If necessary, increase the amount of tissue and repeat the procedure |
| Grasping too much tissue between the jaws                    | Open the jaws and reduce the amount of tissue which is grasped, and reactivate the seal cycle                                                         |
| Grasping a metal object                                      | Avoid grasping objects, such as staples, clips or encapsulated sutures in the jaws of the instrument                                                  |
| Activating in excess pooled fluids around the instrument tip | Minimize or remove excess fluids<br>Reactivate the seal cycle without repositioning the instrument                                                    |
| Excessive tissue eschar on electrode tips                    | Use a wet gauze pad to clean surfaces and edges of instrument jaws                                                                                    |

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                                                                                                                                                                                                                                        |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Instrument Error</b><br><i>Indicated By:</i> <ul style="list-style-type: none"> <li>• A <b>three-pulsed</b> alert tone</li> <li>• Instrument receptacle display illuminates <b>amber</b> and flashes continuously</li> <li>• Generator will not allow RF energy delivery</li> </ul> <i>Causes:</i> <ul style="list-style-type: none"> <li>• Generator is receiving an activation request from the instrument</li> </ul> <i>To Resolve:</i> <ol style="list-style-type: none"> <li>1. Disconnect instrument from the generator</li> <li>2. Ensure instrument activation button is not being pressed</li> <li>3. Reconnect instrument to the generator</li> <li>4. Confirm the instrument receptacle display illuminates green</li> </ol> If the instrument error reoccurs: <ul style="list-style-type: none"> <li>• Do not use the instrument</li> <li>• Use a different CoolSeal instrument</li> </ul> |                                                                                                                                                                                                                                                        |
| <b>Possible use conditions include:</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                                        |
| Inadvertent depression of the instrument activation button during instrument connection                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   | Remove anything depressing the instrument activation button and reconnect instrument                                                                                                                                                                   |
| The instrument switch is malfunctioning                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   | Replace the instrument                                                                                                                                                                                                                                 |
| <b>Invalid Instrument</b><br><i>Indicated By:</i> <ul style="list-style-type: none"> <li>• A <b>single-pulse</b> alert tone</li> <li>• Instrument receptacle display illuminates and remains <b>red</b></li> <li>• Generator will not allow RF energy delivery</li> </ul> <i>Causes:</i> <ul style="list-style-type: none"> <li>• An unusable instrument has been connected</li> </ul> <i>To Resolve:</i> <ol style="list-style-type: none"> <li>1. Disconnect instrument from the generator</li> <li>2. Ensure instrument is CoolSeal compatible</li> <li>3. Reconnect instrument to the generator</li> <li>4. Confirm the instrument receptacle display illuminates green</li> </ol> If the instrument error reoccurs: <ul style="list-style-type: none"> <li>• Do not use the instrument</li> <li>• Use a different CoolSeal instrument</li> </ul>                                                     |                                                                                                                                                                                                                                                        |
| <b>Possible use conditions include:</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                                        |
| Instrument connected is not CoolSeal technology compatible                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | Confirm CoolSeal technology compatibility of the instrument from the instrument instructions for use                                                                                                                                                   |
| Instrument connected has been used previously                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | Discard instrument                                                                                                                                                                                                                                     |
| Instrument is not usable with software version                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | Confirm the software version required by the instrument is not greater than the software version labelled on the bottom of the CoolSeal Generator<br>For a software upgrade refer to the Software Upgrade section of the CoolSeal Generator User Guide |

## Warranty

Except as otherwise expressly stated in an agreement between Hologic and its original customer ("Customer"), Hologic equipment ("Equipment") is warranted to the original Customer to perform substantially in accordance with published product specifications for one (1) year starting from the date of shipment ("Warranty Period"). Replacement parts and remanufactured items are warranted for the remainder of the Warranty Period or 90 days from shipment, whichever is longer. Consumable supplies are warranted to conform to published specifications for a period ending on the expiration date shown on their respective packages. Licensed software is warranted to operate in accordance with published specifications. Services are warranted to be supplied in a workman-like manner. Hologic does not warrant that use of Products will be uninterrupted or error-free, or that Products will operate with non-Hologic authorized third-party products.

These warranties do not apply to any item that is: (a) repaired, moved, or altered other than by Hologic authorized service personnel; (b) subjected to physical (including thermal or electrical) abuse, stress, or misuse; (c) stored, maintained, or operated in any manner inconsistent with applicable Hologic specifications or instructions, including Customer's refusal to allow Hologic

recommended Software upgrades; or (d) designated as supplied subject to a non-Hologic warranty or on a pre-release or “as-is” basis.

### Disposal

The used disposable device must be treated as biohazardous waste and disposed of according to standard practices of the hospital or clinic where the treatment is performed.

Consult your local regulations for disposal and/or electronics recycling as applicable. Do not place into a municipal waste system unless authorized to do so by local authorities.

### Technical Support and Product Return Information

Contact Hologic Technical Support if the CoolSeal Trinity fails to operate as intended. If product is to be returned to Hologic for any reason, Technical Support will issue a Returned Materials Authorization (RMA) number and biohazard kit if applicable. Return the CoolSeal Trinity according to the instructions provided by Technical Support.

If applicable, return used or opened product according to the instructions provided with the Hologic-supplied biohazard kit.

### Complaints and Serious Incident Reporting

Report any complaints or problems in the quality, reliability, safety, or performance of this product to Hologic. If the device has caused or added to patient injury, immediately report the incident to Hologic Authorized Representative and Competent authority of the respective member state or country. The Competent Authorities for medical devices are usually the individual Member States' Ministry of Health, or an agency within the Ministry of Health.

### Technical Support Information

Hologic, Inc.  
250 Campus Drive  
Marlborough, MA 01752 USA  
Phone: 1.800.442.9892

## After Surgery

### Warning

Do not reuse or resterilize the CoolSeal Trinity.

## Pre-Clinical Study

### Caution

There is no animal data qualified to predict the effectiveness of this device in sealing vessels containing atherosclerotic plaque.

Product performance of the device was established in a chronic in-vivo porcine model. The results showed that no animals studied experienced any hemostatic complications related to the device during the minimum 21-day survival period. A variety of tissue types and vessels were evaluated to demonstrate effective sealing in arteries, veins, and vascular bundles up to and including 7 mm.

The United States FDA clearance of this device was not based on human clinical testing.

| Vessel Type | Tissue/Vessel Name | Vessel Size Range |
|-------------|--------------------|-------------------|
| A/V Bundle  | Mesentery          | < 2.5 mm          |
|             | Ovarian            | 2.0 mm – 6.5 mm   |
|             | Uterine            | 1.8 mm – 5.4 mm   |
|             | Splenic            | 1.8 mm – 7.0 mm   |
|             | Short Gastric      | 5.5 mm            |
| Artery      | Renal              | 3.5 mm – 7.0 mm   |
|             | Splenic            | 3.0 mm – 7.0 mm   |
|             | Gastrointestinal   | 4.0 mm – 5.0 mm   |
| Vein        | Renal              | 7.0 mm            |
|             | Splenic            | 2.0 mm – 7.0 mm   |
|             | Short Gastric      | 5.5 mm            |
|             | Gastrointestinal   | 5.0 mm – 6.0 mm   |

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All rights reserved.
© 2025 Hologic, Inc.