
Manufacturer: Gynesonics, Inc.

Symbols Glossary

SYMBOL	SYMBOL TITLE	EXPLANATORY TEXT	STANDARD REFERENCE	STANDARD TITLE
	DO NOT re-use	Indicates a medical device that is intended for one use or for use on a single patient during a single procedure.	ISO 15223-1 #5.4.2	Medical devices – Symbols to be used with medical device labels, labeling, and information to be supplied.
			EN 980 #5.2	Symbols for use in the labeling of medical devices.
			ISO 7000-1051	Graphical symbols for use on equipment.
	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1 #5.3.7	Medical devices – Symbols to be used with medical device labels, labeling, and information to be supplied.
			ISO 7000-0632	Graphical symbols for use on equipment.
			EN 980 #5.17.3	Symbols for use in the labeling of medical devices.
	Humidity limitation	Indicates the range of humidity to which the medical device can be safely exposed.	ISO 15223-1 #5.3.8	Medical devices – Symbols to be used with medical device labels, labeling, and information to be supplied.
			ISO 7000-2620	Graphical symbols for use on equipment.
	Importer	To indicate the entity importing the medical device into the locale	ISO 7000-3725	Graphical symbols for use on equipment

DESCRIPTION	
Name of Product	Sonata® IUUS Probe Return Shipper Kit, RTN-001
Product Description	The Sonata IUUS Probe Return Shipper Kit ("Return Shipper Kit") is a product accessory that shall enable shipment of a non-sterile (after use and cleaning) Sonata Intrauterine Ultrasound (IUUS) Probe (IUSP-001 or IUSP-002).
Contents of Kit, RTN-001	• IUUS Probe Backer Card (1) • IUUS Probe Backer Card Tube (1) • IUUS Probe Backer Card Strap (1) • Biohazard Bag, Red (1) • Leak Proof Bag (2) • Bubble Bag for IUUS Probe Packaging (1) • Box, Individual IUUS Probe Shipment (1) • Used Medical Device label (1) • Biohazard label (1) • IUUS Probe Return Shipper Kit Instructions (1) • Device Return Label (1)
Dimensions length x width x height	As packaged for shipment: 24" x 6" x 4"
Storage Conditions of IUUS Probe	Temperature: -20 - 40°C (-4 - 104°F) Relative Humidity: 10-90%
Transit Conditions	
CAUTIONS	
CAUTIONS	For use only with the Sonata® Intrauterine Ultrasound Probe. Adequacy of this packaging system has not been demonstrated for any other instrument. All materials are single-use only.  In the event of a serious incident related to use of the Sonata IUUS Probe Return Kit, notify Gynesonics. If located within the European Union, also notify the competent authority of the EU Member State in which the user and/or patient is established.
INSTRUCTIONS	
Probe Preparation for Packaging	Prior to packaging for return after use, ensure that the <i>Sonata Intrauterine Ultrasound Probe</i> has been cleaned and dried per the <i>Sonata System IUUS Probe Instructions for Use</i> .
Packaging Probe for Shipment	Appendix A
REFERENCES	
CE Mark	 MDR (EU) 2017/745

UKCA Mark	 2002 No. 618 with 2019 No. 791 and 2020 No. 1478
Related Gynesonics Product	Sonata Intrauterine Ultrasound Probe IUSP-002
Reference Document	Sonata System Operators' Manual/Instructions for Use; Sonata System IUUS Probe Instructions for Use
Manufacturer Contact	Gynesonics Inc. 200 Cardinal Way #250 Redwood City, CA, 94063 USA Telephone: +1.650.216.3860 www.gynesonics.com EU: customersupport@gynesonics.com US: USCustomerService@gynesonics.com
European Authorized Representative	 Obelis s.a Bd. Général Wahis 53 1030 Brussels, BELGIUM Tel.: +(32) 2. 732.59.54 Fax: +(32) 2. 732.60.03 Email: mail@obelis.net
UK Responsible Person	Obelis UK Ltd. Sandford Gate East Point Business Park Oxford, OX4 6LB United Kingdom +44.1491.378012 info@obelis.co.uk
Importer	 Hologic Ireland Ltd. Sir Rogerson's Quay 70, D02 R296 Dublin, Ireland
	 Hologic Ltd. Oaks Business Parks, Crewe Road, Wythenshawe, Manchester, M23 9HZ, United Kingdom
	 Hologic Suisse SA World Trade Center, Avenue de Gratta-Paille 2, 1018 Lausanne, Switzerland

Originally approved and supplied in the English language

Sonata® IUUS Probe Return Shipper Kit RTN-001 Instructions

English


Sonata, Gynesonics, and the logo forms of the foregoing are registered trademarks and/or service marks of Gynesonics Inc. All other trademarks are properties of their respective owners.

APPENDIX A


Installing the Sonata Intrauterine Ultrasound Probe on the backer card

Ensure that the tip of the probe is unarticulated / straight / at 0-degrees.

Insert the distal end of the IUUS probe into the tip protector tube on the backer card with the handle facing down.

Lay the probe handle on the backer card. Lift the tabs on both ends of the handle for proper positioning.

Fold the strap over the probe handle to secure device. Fold short end of strap first, lifting the tab. Then fold over the longer end of strap by pushing the tab out and under the short end tab.

Loop the probe cable around the backer card flap about three times. Cable must be taut without deforming the backer card.


Carefully insert the electrical connector into its constraint on the backer card

The connector will fit snugly into the constraint. One way to seat the connector into the constraint is to first insert one of its corners and then slide the rest of the connector in, bringing the other corner in last.


CONNECTOR
CONSTRAINT

Insert the strain relief into its holding tab to lock the connector in place.

Continue inserting the device cable into its tabs to secure it in place.

Fold the top half (flap) of the backer card over the device.

Place the device (loaded onto the backer card) inside the provided biohazard bag. Fold the bag around the device in such a way as to ensure that the biohazard symbol is visible.

Insert the bagged device into one leak-proof bag. Fold the bag around the device and remove as much air from it as possible. Ensure that the biohazard symbol is visible through the leak-proof bag. Seal the bag (ziplock). Insert the now double-bagged device into the second leak-proof bag, remove air and seal that bag. At this point, the device is inside one biohazard bag and two leak-proof bags, triple-bagged.

Affix the biohazard label on the outside of the second leak-proof bag.

Fill out the provided Device Return Label, LS 04796-001 and attach it to the second leak-proof bag next to the Biohazard label as shown below.

Place on Bag

RMA #: _____

Known infectious substance, if any: _____

Condition (clean/pre-clean/disinfect?): _____

Product (fill out as applicable)

☐ RFA Handpiece

☐ IUUS Probe

☐ Other: _____

Model#: _____

Lot/SN #: _____

LS 04796-001 Rev D

Place on Bag
RMA #:
Known infectious substance, if any:
Condition (clean/pre-clean/disinfect?):
Product (fill out as applicable)
RFA Handpiece
IUUS Probe
Other:
Model #:
Lot/SN #:
LS 04796-001 Rev D

IF it is known that the device was in contact with an infectious substance, indicate the microorganism as shown:

Place on Bag	
RMA #:	<u>1059</u>
Known infectious substance, if any:	<u>Hepatitis B</u>
Condition (clean/pre-clean/disinfect?)	<u>Pre-clean</u>
Product (fill out as applicable)	
☐ RFA Handpiece	
☒ IUUS Probe	
☐ Other: _____	
Model#:	<u>IUSP-002</u>
Lot #:	<u>10026</u>
LS 04796-001 Rev D	

Place on Bag
RMA #:
1059
Known infectious substance, if any: Hepatitis B
Condition (clean/pre-clean/disinfect?): Pre-clean
Product (fill out as applicable)
RFA Handpiece
IUUS Probe
Other:
Model #: IUSP-002
Lot #: 10026
LS 04796-001 Rev D

1. If shipping a Gynesonics probe for off-site reprocessing, complete and attach LS 04796-001 Device Return label to the second leak-proof bag next to the Biohazard label, as pictured above. (No RMA number would be included, in this case.) Ship the probe to a validated reprocessor for processing.
2. If shipping a Gynesonics probe for any other reason (e.g. complaint / investigation associated with a Return Materials Authorization (RMA)), complete and attach LS 04796-001 Device Return label to the second leak-proof bag next to the Biohazard label, as pictured above. Ship the device to the following address:

Gynesonics, Inc.
 Attention: Quality Assurance
 200 Cardinal Way #250
 Redwood City, CA 94063
 USA

(Telephone:650-216-3860)

Insert the now triple-bagged device into the provided bubble bag. Do not seal bubble bag. Fold bubble bag as shown.


Place the wrapped device inside the individual box. Securely seal the box with packaging tape. Affix the “Used Medical Device” label and ship device to appropriate destination.

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