

Faxitron® Pro+

Specimen Radiography System



User Guide

MAN-10285 Revision 001

HOLOGIC®

Faxitron[®] Pro+

Specimen Radiography System

User Guide

For Software Version 3.1

Part Number MAN-10285

Revision 001

July 2023

Faxitron[®] Pro+ was previously known as UltraFocusXL.

HOLOGIC[®]

Technical Support

USA: +1.877.371.4372

Europe: +32 2 711 4690

Asia: +852 37487700

Australia: +1 800 264 073

All Other: +1 781 999 7750

Email: BreastHealth.Support@hologic.com

© 2023 Hologic, Inc. Printed in the USA. This manual was originally written in English.

Hologic, Faxitron, and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. All other trademarks, registered trademarks, and product names are the property of their respective owners.

This product may be protected by one or more U.S. or foreign patents as identified at www.Hologic.com/patents.

Table of Contents

List of Figures	ix
-----------------------	----

List of Tables	xiii
----------------------	------

1: Introduction	1
------------------------------	----------

1.1 Intended Use.....	1
1.2 System Capabilities.....	1
1.3 User Profiles.....	1
1.4 Quality Control Requirements.....	1
1.5 Where to Find the Installation Instructions.....	1
1.6 Where to Find Technical Description Information.....	2
1.7 Warranty Statement.....	2
1.8 Technical Support	2
1.9 Product Complaints.....	2
1.10 Hologic Cybersecurity Statement	2
1.11 Symbols	3
1.12 Descriptions of Warnings, Cautions, and Notes	4

2: General Information	5
-------------------------------------	----------

2.1 System Overview	5
2.2 Safety Information	6
2.2.1 Radiation Safety and Emissions Limits	6
2.2.2 Interlocks	6
2.2.3 Laser Locator.....	7
2.3 Warnings and Precautions.....	9
2.4 Compliance.....	13
2.4.1 Compliance Requirements	13
2.4.2 Compliance Statements	14
2.5 DICOM Standard.....	15
2.6 Interoperable Connections.....	15
2.7 Network / User Settings.....	16
2.8 Label Locations	17

3: Components, Controls, and Indicators	19
--	-----------

3.1 System Components.....	19
3.1.1 Imaging Cabinet Components.....	20
3.2 System Connections.....	22
3.3 How to Move the Faxitron Pro+ System.....	23

4: Quick Start – Basic Operation of the Faxitron Pro+ System	25
---	-----------

4.1 System Start Up.....	25
4.1.1 Power Up the System.....	25

Faxitron Pro+ Specimen Radiography System User Guide

Table of Contents

4.1.2	Start the Software	26
4.1.3	Calibrate the System.....	26
4.1.4	Sample Data	27
4.1.5	Acquiring Images	29
4.2	How to Shut Down the System.....	31
4.2.1	How to Remove All Power from the System.....	31
5: VisionLS Software Operation and User Interface		33
5.1	Overview of VisionLS Software.....	33
5.2	Calibration	33
5.2.1	Calibration Introduction.....	33
5.2.2	Calibration Sequence.....	34
5.2.3	Advanced Calibration.....	35
5.3	Automatic Exposure Control (AEC)	38
5.3.1	Preshot AEC Setup	38
5.4	Sample Data.....	40
5.4.1	Manual Entry	41
5.4.2	Continue with Current Procedure	42
5.4.3	Select from Work List (Remote Modality Work List)	42
5.4.4	Select from Database	42
5.5	Image Acquisition.....	44
5.6	Tool Bars and Menus.....	44
5.6.1	Left Tool Bar	44
5.6.2	Right Tool Bar	54
5.6.3	Bottom Tool Bar	57
5.6.4	Top Menu Bar	59
5.7	Current Sample Images Browser	82
5.8	Database/Review Mode	84
5.8.1	THE REVIEW PANE.....	86
5.8.2	Database Top Menu Bar	87
5.8.3	Edit or Delete Sample or Procedure Information from the Database.....	93
5.9	DICOM Functionality.....	96
5.9.1	Modality.....	96
5.9.2	Default DICOM Servers.....	97
5.9.3	DICOM Save Local	98
5.10	User Annotations	99
5.10.1	Overview	99
5.10.2	Operations	99
5.10.3	Design Mode and Mouse Right-Click Menu	101
5.10.4	Show – Hide Annotations	105
5.11	Help Button	106
6: Quality Control		107
6.1	Required Quality Control Procedures	107
6.2	Service and Troubleshooting.....	107

6.2.1	Service	107
6.2.2	Troubleshooting.....	107
6.3	Cybersecurity	107
7: Maintenance, Cleaning, and Disinfecting		109
7.1	General Information	109
7.1.1	For General Cleaning	109
7.1.2	For Disinfecting.....	109
7.1.3	To Prevent Possible Injury or Equipment Damage.....	109
7.2	Component Specific Care	110
7.2.1	How to Clean the Image Display Monitor	110
7.2.2	Power Cord	110
7.3	Maintenance	111
7.3.1	Preventive Maintenance Schedule	111
Appendix A System Specifications		113
A.1	Product Measurements	113
A.2	Operation and Storage Environment	114
A.2.1	General Conditions for Operation.....	114
A.2.2	General Conditions for Transport and Storage	114
A.3	Technical Specifications	115
A.3.1	System Specs	115
A.3.2	Beam Coverage and Source to Object Distance Chart.....	115
A.3.3	X-Ray Tube	116
A.3.4	Computer Specifications.....	117
8: Glossary of Terms		119

List of Figures

Figure 1: Faxitron Pro+ System	5
Figure 2: Generator Timing Mismatch Error Message	6
Figure 3: Laser Caution Label	8
Figure 4: Label Locations	17
Figure 5: System Components	19
Figure 6: X-Ray Chamber with Magnification Levels	20
Figure 7: Mag Tray Markings	21
Figure 8: System Imaging/Display Monitor	22
Figure 9: Back Panel Connections	23
Figure 10: Key Switch	25
Figure 11: VisionLS Software Icon	26
Figure 12: Calibration Page	26
Figure 13: Calibration Dialog Box	26
Figure 14: Calibration Dialog Box	26
Figure 15: Home Menu Screen	27
Figure 16: Home Menu	27
Figure 17: Sample Data Entry Menu	28
Figure 18: Sample Information Editor Dialog Box – Required Fields Empty and Filled In	28
Figure 19: Image Acquisition Screen	29
Figure 20: Chamber Floor – Active Image Area	29
Figure 21: X-Ray Chamber with Mag Tray Installed	30
Figure 22: Start X-Ray Buttons – Control Panel Touch Screen and Software	30
Figure 23: Image Acquisition Screen	31
Figure 24: Calibrate Button	34
Figure 25: Calibration Screen	34
Figure 26: Calibration Progress	35
Figure 27: Calibration Progress	35
Figure 28: Tools Menu – Advanced Tools – X-Ray Calibration/Exposure Setup	35
Figure 29: X-Ray Generator and Exposure Settings Dialog Box	36
Figure 30: X-Ray Generator and Exposure Settings Dialog Box	36
Figure 31: Tools Menu – Advanced Tools – PRESHOT AEC Setup	38
Figure 32: Preshot AEC Settings Dialog	39
Figure 33: Home Menu	40
Figure 34: Sample Data Entry Options	41
Figure 35: Sample Information Editor Dialog Box – Blank on Left, Complete on Right	41
Figure 36: Select from Database Page	43
Figure 37: Select from Database – Using Search Parameters	43
Figure 38: Left Tool Bar	44
Figure 39: View/Change Sample Information Button	45
Figure 40: Sample Information Dialog Box - Imaging Mode: Before Image (Left), After Image (Right)	45
Figure 41: Sample Information Editing Dialog – Review Mode	46
Figure 42: Enhance Image Button	46

Faxitron Pro+ Specimen Radiography System User Guide

Table of Contents

Figure 43: Window Level and Contrast Adjustment Tool	47
Figure 44: Window Level and Contrast Adjustment Tool	48
Figure 45: Window Level and Contrast Adjustment Tool	49
Figure 46: Example of Inverting: Inverted on the Left.....	49
Figure 47: Digital Zoom Tool Button and Dialog Box	50
Figure 48: Image in Default Zoom Level.....	51
Figure 49: Same Image of Snake at 1.00x Zoom	51
Figure 50: Annotation Tools Button and Dialog Box	52
Figure 51: PACS Store Button	52
Figure 52: PACS Print Button.....	52
Figure 53: Exposure Control Button and Dialog Box.....	53
Figure 54: Back to Main Menu Button	54
Figure 55: Pin (Un-Docking) Icon.....	54
Figure 56: Right Side Menu Bar	54
Figure 57: Reveal/Hide Top Menu Button	55
Figure 58: Image Information Button and Box	55
Figure 59: Toggle Annotations Button.....	56
Figure 60: Sample Viewer On-Screen Window	57
Figure 61: Bottom Tool Bar – Imaging & Review Modes	57
Figure 62: Bottom Tool Bar – Middle Section (B)	58
Figure 63: Bottom Tool Bar – Review Mode, Right Third	58
Figure 64: Bottom Tool Bar – Buttons/Status Bar – Imaging Mode	59
Figure 65: Top Menu	59
Figure 66: File Menu.....	59
Figure 67: Edit Menu.....	59
Figure 68: View Menu.....	60
Figure 69: View Menu - Zoom Mode Options.....	60
Figure 70: View Menu - Image Averages Options	61
Figure 71: View Menu – Image Perspective.....	61
Figure 72: View Menu - Grid Options	62
Figure 73: Grid Sample, Pitch Set to 256.....	62
Figure 74: View Menu – User Window Levels.....	63
Figure 75: View Menu – Default Window Level.....	63
Figure 76: View Menu – Default Aspect for Raw Files.....	64
Figure 77: View Menu – Enhance Image.....	64
Figure 78: View Menu – Magnifier Size.....	65
Figure 79: View Menu – Review Layout	65
Figure 80: View Menu – Review Layout (1 + (2 X 1)) Triple	66
Figure 81: Tools Menu.....	66
Figure 82: Advanced Tools Menu	67
Figure 83: Advanced Tools – Image Enhancement Setup	68
Figure 84: Login Utility Dialog Box.....	69
Figure 85: Advanced Tools – Cyclic Automated Acquisitions.....	70
Figure 86: Advanced Tools – Other Advanced User Tools – Multi Scale Enhancement	71
Figure 87: Advanced Tools - Properties.....	71

Figure 88: Tools Menu.....	72
Figure 89: Tools – Save/Store Current Image.....	73
Figure 90: LIS Save File As Dialog.....	73
Figure 91: Save as DICOM File Dialog	74
Figure 92: Save As Dialog.....	75
Figure 93: Print Current Image	75
Figure 94: Image Tools Menu.....	76
Figure 95: Image Tools Menu - Magnifier	77
Figure 96: Image Tools Menu – ROI Contrast.....	78
Figure 97: Image Tools Menu – ROI Contrast.....	78
Figure 98: Image Tools Menu - Ruler.....	79
Figure 99: Image Tools Menu - ROI Statistics	79
Figure 100: Image Tools Menu - ROI Histogram	80
Figure 101: Image Tools Menu - ROI Profile.....	80
Figure 102: Image Tools Menu – Line Profile	81
Figure 103: DICOM Utilities Menu	81
Figure 104: Help Menu	82
Figure 105: About Vision.....	82
Figure 106: Current Sample Images Browser Button.....	82
Figure 107: Current Sample Images Thumbnail Strip	83
Figure 108: Current Sample Images Thumbnail Strip	83
Figure 109: Open Database Button.....	83
Figure 110: Review Mode Button in the Home Menu	84
Figure 111: Image Database – Main Sections	84
Figure 112: Change Review Layout Button and Menu.....	86
Figure 113: Review Panes	86
Figure 114: Database Top Menu Bar	87
Figure 115: Save Images Dialog	88
Figure 116: Delete Image(s) Confirmation	89
Figure 117: Additional Options Menu	89
Figure 118: Additional Options Menu - Open.....	89
Figure 119: Additional Options Menu - Print.....	90
Figure 120: Database View Menu.....	90
Figure 121: Database Search Menu.....	91
Figure 122: Database Start Search Dialog Box	91
Figure 123: Database Tools Menu.....	92
Figure 124: Import Records Confirmation Dialog.....	93
Figure 125: Edit or Delete Sample Information Menu.....	93
Figure 126: Sample Information Editor Dialog Box	94
Figure 127: Edit Sample Information Confirmation Dialog Box	94
Figure 128: Delete Sample Information Dialog Box.....	94
Figure 129: Edit or Delete Procedure Information Menu	95
Figure 130: Sample/Procedure Information Editor	95
Figure 131: Edit Procedure Information Confirmation Dialog Box	96
Figure 132: Sample Information Editor Dialog box	96

Faxitron Pro+ Specimen Radiography System User Guide

Table of Contents

Figure 133: Top Menu – DICOM Utilities	96
Figure 134: DICOM Modality Settings Dialog Box – Default Settings Shown.....	97
Figure 135: DICOM Networking Dialog Box.....	98
Figure 136: Image Right Click Menu – Save/Store Current Image	98
Figure 137: Annotation Tool Button and Dialog Box.....	99
Figure 138: Annotation Tool Pallet.....	99
Figure 139: Annotation Tool Right Click Menu and Properties Menu.....	101
Figure 140: Annotation Tool Properties – Fixed Dialog Box.....	102
Figure 141: Adding Preset Annotation Tools.....	103
Figure 142: Preset Annotation Tools Drop-Down Menu.....	103
Figure 143: Annotation Tool Right Click Menu.....	104
Figure 144: Annotation Tools - Right Click Menu – Default Properties Menu	105
Figure 145: Toggle Annotations Button.....	105
Figure 146: Help Button – Home Menu.....	106
Figure 147: Help Menu/Page.....	106
Figure 148: Path+ System Measurements	113

List of Tables

Table 1: Required Procedures107

Table 2: User Preventive Maintenance.....111

Table 3: Service Preventive Maintenance111

Chapter 1 Introduction

Read all this information carefully before operating the system. Follow all warnings and precautions as stated in this manual. Keep this manual available during procedures.

**Note**

Hologic configures some systems to meet specific requirements. Your system configuration may not have all the options and accessories included in this manual.

1.1 Intended Use

Faxitron® Pro+ is a fully shielded x-ray cabinet offering a large field of view and high-resolution for imaging a broad range of biological and light NDT samples.

1.2 System Capabilities

Faxitron Pro+ is a fully shielded x-ray cabinet offering a large field of view and high-resolution with a 43cm x 43cm detector for imaging a broad range of biological and light NDT samples.

The Faxitron Pro+ system provides geometric magnification of up to 10X.

The Faxitron Pro+ system is fully shielded and meets FDA radiation safety standards. It is self-contained with no additional x-ray shielding required.

The Faxitron Pro+ is easy to use. No additional training or specialized x-ray requirements are needed to operate the system.

Hologic recommends that users receive training in basic ionizing radiation safety before using the system.

1.3 User Profiles

- The Faxitron Pro+ system does not utilize multiple User Profiles.
- System access is managed through the use of a physical key which serves as an ON/OFF switch.
- Service personnel have password protected access to all system functions.

1.4 Quality Control Requirements

Perform all Quality Control tests within the correct time frame.

1.5 Where to Find the Installation Instructions

Installation instructions are available in the *Service Manual*.

1.6 Where to Find Technical Description Information

Technical description information is available in the *Service Manual*.

1.7 Warranty Statement

Except as otherwise expressly stated in the Agreement: i) Equipment manufactured by Hologic is warranted to the original Customer to perform substantially in accordance with published product specifications for one (1) year starting from the date of shipment, or if Installation is required, from the date of Installation (“Warranty Period”); ii) digital imaging mammography x-ray tubes are warranted for twenty-four (24) months, during which the x-ray tubes are fully warranted for the first twelve (12) months and are warranted on a straight-line prorated basis during months 13-24; iii) replacement parts and remanufactured items are warranted for the remainder of the Warranty Period or ninety (90) days from shipment, whichever is longer; iv) consumable Supplies are warranted to conform to published specifications for a period ending on the expiration date shown on their respective packages; v) licensed Software is warranted to operate in accordance with published specifications; vi) Services are warranted to be supplied in a workman-like manner; vii) non-Hologic Manufactured Equipment is warranted through its manufacturer and such manufacturer’s warranties shall extend to Hologic’s customers, to the extent permitted by the manufacturer of such non-Hologic Manufactured Equipment. Hologic does not warrant that use of Products will be uninterrupted or error-free, or that Products will operate with non-Hologic authorized third-party products. These warranties do not apply to any item that is: (a) repaired, moved, or altered other than by Hologic authorized service personnel; (b) subjected to physical (including thermal or electrical) abuse, stress, or misuse; (c) stored, maintained, or operated in any manner inconsistent with applicable Hologic specifications or instructions, including Customer’s refusal to allow Hologic recommended Software upgrades; or (d) designated as supplied subject to a non-Hologic warranty or on a pre-release or “as-is” basis.

1.8 Technical Support

Refer to the copyright page of this manual for contact information for product support.

1.9 Product Complaints

Report any complaints or problems in the quality, reliability, safety, or performance of this product to Hologic. If the device has caused or added to patient injury, immediately report the incident to Hologic.

1.10 Hologic Cybersecurity Statement

Hologic continuously tests the current state of computer and network security to examine possible security problems. When necessary, Hologic provides the updates to the product.

For Cybersecurity Best Practices documents for Hologic products, refer to the Hologic Internet site at www.Hologic.com.

1.11 Symbols

This section describes the Symbols on this system.

Symbol	Description
	Discard electrical and electronic equipment separately from standard waste. Send decommissioned material to Hologic or contact your service representative.
	Manufacturer
	Date of Manufacture
 or 	Caution—Radiation
	Caution – Potentially infectious materials
	Warning Electricity
	Caution
	Follow instructions for use
	Catalog number
	Serial number

1.12 Descriptions of Warnings, Cautions, and Notes

Descriptions of Warnings, Cautions, and Notes used in this manual:



WARNING!

The procedures that you must follow accurately to prevent possible dangerous or fatal injury.



Warning:

The procedures that you must follow accurately to prevent injury.



Warning:

Warns the reader that ionized radiation is emitted into the labeled area when the x-ray beam is energized. The Faxitron Pro+ system has safety interlocks to prevent the labeled area from being accessed while the x-ray beam is energized.



Warning:

Warns the reader to use appropriate protection, such as gloves or an approved specimen container, to avoid exposure to potentially infectious materials such as patient tissue samples.



Caution:

The procedures that you must follow accurately to prevent damage to equipment, loss of data, or damage to files in software applications.



Caution:

Cautions the reader to use care in lifting heavy equipment to avoid personal injury and/or damage to the equipment



Note

Notes show additional information.

Chapter 2 General Information

2.1 System Overview




Figure 1: Faxitron Pro+ System

Figure Legend

1. Control Panel
2. X-Ray Imaging Chamber
3. Image Display Monitor
4. Keyboard & Mouse

2.2 Safety Information

Read and understand this manual before you use the system. Keep the manual available during the operation of the unit.

Always follow all the instructions in this manual. Hologic does not accept responsibility for injury or damage from incorrect system operation. Hologic can arrange for training at your site.

The system has safety interlocks, but the user must understand how to safely operate the system and be aware of the health hazards of x-ray radiation.

2.2.1 Radiation Safety and Emissions Limits

- Radiation shielded cabinet.
- X-Ray chamber door equipped with dual-safety interlocks. *See Sections 2.2.2 of this Manual.*
- Radiation: Less than 0.1 mR/hr at 5 cm (2 in.) from exterior surface at maximum kV.
- Certified to comply with standards set forth by the US Food and Drug Administration, Center for Devices and Radiological Health, (21 CFR 1020.40).
- Designed to comply with standards set forth by Underwriters Laboratory the US and Canada.

2.2.2 Interlocks

The Faxitron Pro+ has two independent safety interlocks on the chamber door. Opening the door will immediately terminate x-ray generation.

X-ray generation stops automatically when there is an interruption of an x-ray exposure. An x-ray exposure interruption may be caused by opening the cabinet door, pressing the stop button (either physical or software), or loss of communication between x-ray subsystem and the computer. The interruption may, or may not, generate an error message. *See figure below.*



Figure 2: Generator Timing Mismatch Error Message

1. Click **OK**.
2. Close and/or verify cabinet door is fully closed.
3. Repeat acquisition.

**Note**

Contact Hologic if the error message is persistent.

2.2.3 Laser Locator

The Faxitron Pro+ utilizes a pair of Laser generated lines to aid in positioning samples in the middle of the imaging field.

Compliance

The lasers comply with:

- IEC 60825-1:2007
- 21 CFR 1040.10 and 1040.11 except for deviations pursuant to Laser Notice No. 50, dated June 24, 2007.

The lasers require no assembly.

**Note**

There is no authorized operator maintenance on the lasers. The lasers are to be serviced only by a Hologic Authorized Representative.

**Caution:**

CLASS 2 LASER RADIATION WHEN THE CHAMBER DOOR IS OPEN. DO NOT STARE INTO THE BEAM.

**Caution:**

A laser beam can be damaging to the eyes. DO NOT look into the laser aperture when the laser is in operation.

**Caution:**

The laser apertures are located in the ceiling of the X-Ray chamber.

**Caution:**

Be aware that laser light can also be dangerous when reflected off a mirror-like surface.

**Note**

For further safety information regarding lasers, refer to ANSI-Z136.1 STANDARD FOR THE SAFE USE OF LASERS, available from the Laser Institute of America, phone 1-800-345-2737.

Laser CAUTION Label

Located on the upper right corner of the door. *Refer to Figure 4: Label Locations.*




Figure 3: Laser Caution Label

2.3 Warnings and Precautions



WARNING!

Only trained Service Engineers authorized through Hologic can open any of the panels. This system contains lethal voltages.



WARNING!

Electrical equipment used near flammable anesthetics can cause an explosion.



WARNING!

The user must correct problems before the system is used. Contact an approved service representative for preventive maintenance.



WARNING!

The protection provided by the equipment is reduced if the equipment is used in a method not specified by Hologic.



WARNING!

No modification of this equipment is allowed.



WARNING!

Do not attempt to use any x-ray cabinet that could present a shock hazard. Immediately contact Hologic or your distributor.



WARNING!

To prevent shock hazard, do not use the system unless the blades of the polarized plug can be fully inserted in the electrical outlet.



WARNING!

Before each use, check that the power cord is in good condition and is fully inserted in the electrical outlet.



WARNING!

Make sure that all power cords for the system are in good condition. A damaged power cord can be an electrical shock hazard. When disconnecting the cabinet, or computer components, from power, always grasp the plug at the insertion point and pull lightly. NEVER pull on the cord to unplug the unit.



WARNING!

Before moving the system, make sure that the power cords and any cables are unplugged and secured.



WARNING!

Before adjusting the position of the system, make sure the power cords and any cables are safely positioned out of the way.



WARNING!

To prevent fire or shock hazard, do not expose the system to rain or moisture.



Warning:

This device contains dangerous material. Send decommissioned material to Hologic or contact your service representative.



Warning:

Control the access to the equipment according to local regulations for radiation protection.



Warning:

Do not use the system if the glass is cracked or damaged.



Warning:

This system can be dangerous to the user. Always follow the safety precautions for x-ray exposures.



Warning:

Cables and cords can pose a tripping hazard. Position cords safely out of the way.



Warning:

Do not move the cabinet with excessive force or speed. Do not make quick stops.



Warning:

When moving the cabinet, be aware of bumps, ramps, inclines, or declines. Use extra caution when moving the cabinet on an uneven or sloped surface.



Warning:

Make sure there is adequate space to completely open the Imaging Cabinet door and to safely insert or remove the specimen tray.



Warning:

Keep the Imaging Cabinet door closed when not in use. Leaving the door open can pose a tripping hazard.



Warning:

Use the Imaging Cabinet door handle to open and close the door. Closing the Imaging Cabinet door incorrectly can cause a pinch hazard.



Caution

Risk of data loss. Do not put any magnetic media near or on devices that create any magnetic fields.



Caution

To prevent possible damage from thermal shock to the Digital Image Receptor, follow the recommended procedure to turn off the equipment.



Caution

Use the least possible amount of cleaning fluids. The fluids must not flow or run.



Caution

To prevent damage to the electronic components, do not use disinfectant sprays on the system.



Caution

To prevent damage to the cabinet and its components, use only the recommended cleaning methods.



Caution:

The system is a laboratory device and not a normal computer. Do not make changes to the hardware or software that are not authorized. Install this device behind a firewall for network security. The computer virus protection or network security for this laboratory device is not provided (for example, a computer firewall). The network security and anti-virus provisions are the responsibility of the user.



Caution

To prevent possible damage to the system, follow the recommended procedure to turn off the equipment.



Caution

The system must be installed and put into service according to the guidance provided in these instructions to ensure its electromagnetic compatibility.



Caution

Portable and mobile RF communications can affect medical electrical equipment.



Caution

To prevent low image quality, perform system calibrations when prompted.



Caution

To prevent low image quality, make sure to remove all sample material from the Imaging Cabinet and specimen tray. Use the recommended cleaning methods to clean and disinfect the Imaging Cabinet and specimen tray.



Caution

To prevent damage, use care when handling and cleaning the specimen tray.

2.4 Compliance

This section describes the system compliance requirements and the responsibilities of the manufacturer.

2.4.1 Compliance Requirements

The manufacturer has the responsibility for the safety, reliability, and performance of this equipment with the following provisions:

- The equipment is used according to the instructions for use.
- The assembly operations, extensions, adjustments, changes, or repairs are performed only by authorized persons.
- The network and communication equipment must be installed to meet IEC Standards.



Caution:

This system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the equipment or shielding the location.



Caution:

The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required), this equipment might not offer sufficient protection to radio frequency communication services. The user may need to take mitigation measures, such as relocating or re-orienting the equipment.



Caution:

The Medical Electrical (ME) Equipment or ME System should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, make sure that the ME Equipment or ME System operates correctly in this configuration.



Caution:

Changes or modifications not expressly approved by Hologic could void your authority to operate the equipment.

2.4.2 Compliance Statements

The manufacturer states this device is made to meet the following requirements:

IEC, EN:

- IEC 61010-1: 2010, 3rd Edition - Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use; Part 1: General Requirements
- IEC 61010-2-091:2012, 1st Edition – Safety requirements for electrical equipment for measurement, control and laboratory use – Part 2-091: Particular requirements for cabinet X-ray systems
- IEC 61010-2-101:2015, 2nd Edition - Safety requirements for electrical equipment for measurement, control and laboratory use – Part 2-101: Particular requirements for IVD medical equipment
- IEC 61326-1 Ed. 2.0 b:2012 - Electrical Equipment for Measurement, Control, and Laboratory Use – EMC Requirements – Part 1: General Requirements
- IEC 61326-2-6 Ed. 2.0 b:2012 - Electrical equipment for measurement, control, and laboratory use – EMC requirements – Part 2-6: Particular requirements – IVD medical equipment
- EN 55011 and CISPR 11 (Class A) - Industrial, Scientific and Medical (ISM) Radio-Frequency Equipment – Electromagnetic Disturbance Characteristics - Limits and Methods of Measurement
- IEC 60825-1:2007 – Safety of Laser Products; Part 1: Equipment classification and requirements
- EN ISO 14971 – Medical Devices – Application of risk management to medical devices
- EN ISO 13485 Medical Devices – Quality management systems – requirements for regulatory purposes
- EN ISO 18113-1 In vitro diagnostic medical devices – Information supplied by the manufacturer (labeling) Part 1: Terms, definitions, and general requirements
- EN ISO 18113-3 In vitro diagnostic medical devices – Information supplied by the manufacturer (labeling) Part 3: In vitro instruments for professional use
- EN ISO 15223-1 Medical Devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements

CFR:

- 21 CFR 1020.40 - Cabinet X-ray Systems.
- 21 CFR 1040.10 – Laser Products.
- 21 CFR 1040.11 – Specific Purpose Laser Products.
- 47 CFR Part 15.5 - General Conditions of Operation.

CAN/CSA:

- CAN/CSA-C22.2 No. 61010-1, 3rd Edition - Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements

2.5 DICOM Standard

The Faxitron Pro+ system uses the latest DICOM Standard:

- NEMA PS3 / ISO 12052, Digital Imaging and Communications in Medicine (DICOM) Standard, National Electrical Manufacturers Association, Rosslyn, VA, USA (available free at <http://www.dicomstandard.org/>)
- See DICOM Conformance Statement for more information.

2.6 Interoperable Connections

The Faxitron Pro+ is capable of operating as a standalone device without any input or output to other devices. Images can be captured, stored, and reviewed directly on the system.

However, the system does come with several interfaces. The primary one is the network connection which allows the system to be connected to a Picture Archiving and Communication System (PACS). Images can then be transmitted from the device directly to the facilities network via our software in a DICOM format. See Section 5.9 for details on DICOM version. Only Hologic employees or trained IT personnel at the user facility should set up this communication.

The system's connections and their intended use are listed below. No other type of devices should be connected to the Pro+ without approval from Hologic.

Connection	Purpose	Device Type	Intended User
Ethernet port at rear	To store radiographic images to networked data storage	PACS or Network storage	Set up by Hologic or facility IT personnel. Used by clinicians.
Ethernet port at rear	To allow remote connection by the Hologic service team, so Hologic can diagnose and repair the system.	Network internet connection	Hologic or facility IT personnel at the direction of Hologic only.
USB port on the monitor	To transfer images from device. This should only be used if PACS or network is not available.	USB storage device (Flash Drive)	Clinicians or IT personnel with approved access to the system.
USB port on the monitor	To service and maintain the system. For example, files copied to and from the computer for software updates or to export log files.	USB storage device (Flash Drive)	Hologic or facility IT personnel
USB port at the rear	To allow connection of a printer. Vision software allows for connection of a DICOM printer or a standard windows compatible printer.	DICOM or windows printer	Clinicians or facility IT personnel

2.7 Network / User Settings

The detector on the Faxitron Pro+ uses an ethernet (network) connection to communicate. Adding users or setting additional network security can prevent the detector from communicating with the software.



Caution:

Set up additional users with local administration rights to ensure the software runs correctly. Work with the Hologic team during installation to make sure the system is configured properly.

2.8 Label Locations

Figure 4: Label Locations




Figure Legend/Label Descriptions

NOTE: Unless specified all labels are from P/N 7120-6300. Label Set, PathVision.

1. "CAUTION – THIS DEVICE PRODUCES X-RAYS WHEN ENERGIZED". (Label D1 or D2) D2 for Canada. All other use D1. Located top/center inside door.
2. "Magnification Level". (P/N 7120-6315, PathVision LX Mag Label Set). Located on left and right side of chamber in front of mag tray rails.
3. "CAUTION – DO NOT OPERATE THIS X-RAY MACHINE WITH THE X-RAY SHIELDING GLASS REMOVED OR DAMAGED". (Label F) Located bottom/center inside door.
4. "CAUTION – X-Rays produced when energized". (Label A1 or A2). A2 for Canada. All others use A1.
5. "WARNING – SAFETY INTERLOCK – Do not tamper with this interlock". (Label E). Located on the right side of the cabinet frame/chamber opening.
6. "CAUTION – Specimens are potentially infectious". (Label C). Located front/center of chamber floor.
7. "ID/Nameplate Label". (P/N LBL-09947). 2x required.
8. "WARNING – Fuse replacement and rating". (Label K). Located bottom of power entry panel.

"CAUTION – No operator serviceable parts inside". (Label O).

Chapter 3 Components, Controls, and Indicators

3.1 System Components




Figure 5: System Components

Figure Legend

- | | |
|---------------------------|--|
| 1. Control Panel | 9. Keyboard and Mouse |
| 2. Power Light (Green) | 10. Ethernet Port |
| 3. Touch Screen | 11. USB Port |
| 4. X-Ray On Light (Red) | 12. System Power Switch and Power Cord Connections |
| 5. Imaging/X-Ray Chamber | 13. Mag Tray Storage Slot (<i>Behind Door</i>) |
| 6. Stop Button | 14. Casters/Feet |
| 7. Key Switch | 15. Computer (<i>Internal, Not Shown</i>) |
| 8. System Imaging Monitor | |

3.1.1 Imaging Cabinet Components

The cabinet houses a shielded imaging/X-ray chamber. The placement of the specimen tray determines the magnification level of the image. A label on each side of the chamber, directly in front of the mag tray slots (or rails), identifies the magnification level for each slot. See figure 21.




Figure 6: X-Ray Chamber with Magnification Levels

Magnification Tray

The Faxitron Pro+ system is supplied with a magnification tray (mag tray) which contains markings that define the active image area margins for each magnification level. The outer most markings define the active image area when the mag tray is positioned on the chamber floor.




Figure 7: Mag Tray Markings

Figure Legend

- | | |
|--|--------------------------------|
| 1. 1.0x Magnification Margins
(Note: Mag Tray on Chamber Floor) | 5. 4.0x Magnification Margins |
| 2. 1.5x Magnification Margins | 6. 5.0x Magnification Margins |
| 3. 2.0x Magnification Margins | 7. 6.0x Magnification Margins |
| 4. 3.0x Magnification Margins | 8. 8.0x Magnification Margins |
| | 9. 10.0x Magnification Margins |



Caution

To prevent damage, use care when handling and cleaning the magnification tray.



Note

Store the magnification tray in the slot provided to the right of the chamber.

3.2 System Connections

Also refer to the Service Manual, P/N MAN-07520, Chapter 3, Installation and Setup, System Set Up.

1. The following cords/cables are connected to the system imaging/display monitor:

- AC Power cord
- Display Port cable (*Video signal input from the system computer*)
- USB-B cable (*extends the computers USB ports to the monitor*)
- Mouse, and
- Keyboard




Figure 8: System Imaging/Display Monitor

Figure Legend

1. AC Power Receptacle
2. Display Port (*Video Cable*) Connection
3. USB-B, from Computer
4. USB, Mouse
5. USB, Keyboard

2. A panel on the back of the Faxitron Pro+ cabinet provides connections for:
 - AC Power cord (*for connection to the facilities AC Mains power*)
 - USB cable (*to allow connection of a printer*)
 - Ethernet cable (*to connect to the facilities network*)

Figure 9: Back Panel Connections

Figure Legend

1. AC Power Connection
2. USB Port
3. Ethernet Port

3.3 How to Move the Faxitron Pro+ System



Caution:

The Faxitron Pro+ cabinet weighs 635kg (1400lbs). If the system must be moved use caution to avoid injury to personnel or damage to the equipment.



Note

The Faxitron Pro+ system was designed as a stationary system and normally should not be moved once installed.



Note

Contact your Hologic representative if the system must be moved.

Chapter 4 Quick Start – Basic Operation of the Faxitron Pro+ System

4.1 System Start Up



Note

Refer to the Faxitron Pro+ Quick Start Guide (P/N 5081-4017)



Note

The start-up sequence must be followed to ensure proper system functionality and image quality.



Note

Allow the system to perform a complete calibration daily prior to system use. The system is configured to perform this calibration if a calibration has not been performed in the previous 12 hours.

4.1.1 Power Up the System

1. Verify there is nothing in the Faxitron Pro+ x-ray chamber.
2. Connect the system to power. If desired, connect the system to the network.
3. Set the power switch on the back of the Faxitron Pro+ cabinet to the “ON” (I) position.
4. Insert the key into the key lock switch on the front of the cabinet and turn it to “ON” (clockwise to the 3 o’clock position). Note that the key is captive and can’t be removed in this position.



Figure 10: Key Switch

5. Turn on the monitor if it is not already on.

4.1.2 Start the Software

The VisionLS software should automatically start when the Windows operating system has completed booting. If the software doesn't automatically start "Double Click" the VisionLS icon on the Windows desktop.



Figure 11: VisionLS Software Icon



Note

Allow the system to perform a complete calibration daily prior to system use.

4.1.3 Calibrate the System

1. Upon initial start-up of the day the software should open to the Calibration page.

Calibration



Figure 12: Calibration Page

2. Click on the **Start Calibration** button and allow the system to complete the calibration. A dialog box will open showing calibration progress.




Figure 13: Calibration Dialog Box



Figure 14: Calibration Dialog Box

3. Upon completion the "Calibration" screen will close, and the "Home Menu" screen will open.

HOLOGIC®

Home Menu




Figure 15: Home Menu Screen



Note

If a calibration has been completed within 12 hours VisionLS software will open directly to the “Home Menu” screen.



Note

The user may calibrate the system at any time by clicking on the **Calibrate** button in the Home Menu.



Note

See Section 5.2, Calibration, for additional calibration information.

4.1.4 Sample Data

1. The **Start Procedure** button on the “Home Menu” page opens a menu which allows the user to select the method for entering *Sample Data*.

Figure 16: Home Menu

2. Click on the **Manual Entry** button which will open the “Sample Information” editor dialog box.



Figure 17: Sample Data Entry Menu

Two side-by-side screenshots of the 'Sample Information' editor dialog box. The left screenshot shows the form with required fields (marked with an asterisk) empty. A red box highlights the 'Operator Name' field, which is split into 'First' and 'Last' sub-fields. A red arrow points to the 'Accept' button, which is greyed out. A text box below the 'Accept' button states: 'The Accept button is Not Active (greyed out)'. The right screenshot shows the same form with the required fields filled in: 'Operator Name' (A, B), 'Sample ID' (C), and 'Accession No.' (D). The 'Accept' button is now active (blue). A text box below the 'Accept' button states: 'All required fields are filled in. The Accept button is Active (Not greyed out)'. Both screenshots show the 'Date of Birth' field with a placeholder 'yyyy/mm/dd' and a dropdown menu for 'Sex'. The bottom of both screenshots shows a grey bar with a back arrow and the text 'Start a New Sample or Procedure'.

Figure 18: Sample Information Editor Dialog Box – Required Fields Empty and Filled In

3. Complete all required fields, which are indicated by an asterisk. Additional fields may be completed at the users' discretion. The **Accept** button will not be active (greyed out) until the required fields are completed. *See figure above.*



Note

If the date is entered it must be entered as yyyy/mm/dd (e.g., 1962/08/13).



Note

The “required fields” may be altered or deleted depending on the application. Contact Hologic support for further information.

- Click on the **Accept** button. The software will open the Image Acquisition Screen. See figure on next page.




Figure 19: Image Acquisition Screen

4.1.5 Acquiring Images

This section provides a brief description of how to acquire images. Image acquisition is covered in greater detail in Section 5.5, Image Acquisition.

- Place the sample in the active image area, or on the mag tray and insert the mag tray into the desired magnification level slot.

Figure 20: Chamber Floor – Active Image Area

Figure 21: X-Ray Chamber with Mag Tray Installed

2. Close the chamber door.
3. Press the green **Start X-Ray** button (either in the software or the control panel touch screen) to initiate image acquisition/x-rays.




Figure 22: Start X-Ray Buttons – Control Panel Touch Screen and Software

4. The system will initiate the exposure and process the data. When complete the image will appear on the screen.




Figure 23: Image Acquisition Screen



Note

Image analysis, manipulation, saving, and retrieval are covered in Sections:

- 5.6 - Tool Bars and Menus
- 5.7 - Current Sample Images Browser
- 5.8 - Database / Review Mode
- 5.9 - DICOM Functionality and
- 5.10 - User Annotations

4.2 How to Shut Down the System

If more than one procedure is scheduled for the day, it is advisable to leave the system on. If the system is no longer needed for the day, begin the shutdown procedure.

1. Turn the key on the front panel to the OFF position, counterclockwise 90°.
2. Wait 30 seconds to ensure the computer powers down completely.
3. Remove the key and store in a secure location to prevent unauthorized use of the system.

4.2.1 How to Remove All Power from the System

1. Make sure that the system is completely powered down.
2. Disconnect the system power cord from the facilities mains receptacle.

Chapter 5 VisionLS Software Operation and User Interface

5.1 Overview of VisionLS Software

Hologic's VisionLS Software is a complete image acquisition and processing package, designed exclusively to work with Hologic's Faxitron line of Life Science X-Ray Generators (such as Faxitron Pro+ and Faxitron Pro systems). In addition to the acquisition and processing functions, the software offers an extensive database module that allows user control and maintenance of image archiving, storage, filing and retrieval. The software is Digital Imaging and Communications in Medicine (DICOM) compliant, allowing saving and printing of sample images in DICOM format 3.0.



Note

Refer to Section 4.1, System Start Up, for system and software startup.

5.2 Calibration

5.2.1 Calibration Introduction

To optimize image quality the Faxitron Pro+ system requires periodic calibration, which requires the acquisition of reference images. These reference calibration images are used to calibrate a "raw" X-Ray image on a pixel-by-pixel basis to generate a "corrected" image for display. Image corrections performed include gain, offset, and linearity.

An x-ray tube warm-up sequence, which helps prolong the life of the tube, is incorporated into the system calibration. The warm-up sequence exercises the x-ray tube through the system's full energy range.

Calibration is a critical component of the system operation. Upon initial start-up of the system, or at any time desired by the user, a system calibration may be requested and performed. Additionally, the software will prompt for a calibration whenever the VisionLS software is started, and a calibration has not been performed with 12 hours.



Note

Perform a system calibration whenever the Faxitron Pro+ X-Ray System is started after prolonged periods of inactivity.



Note

Perform a system calibration if the Faxitron Pro+ X-Ray System is moved to another room with potentially different ambient temperature/humidity settings.

5.2.2 Calibration Sequence

For convenience, VisionLS software combines the numerous steps of calibration into a single unified step.

To calibrate the system at any time, click the **Calibrate** button on the Home Menu (Start-up Screen), and then press **Start Calibration**. See figures below. The Dialog will display a status message during calibration. See figures 26 & 27.




Figure 24: Calibrate Button

Calibration



Figure 25: Calibration Screen



Figure 26: Calibration Progress



Figure 27: Calibration Progress

5.2.3 Advanced Calibration

At times it may be beneficial to incorporate the X-Ray attenuation characteristics of a uniform attenuator in the flat-field reference calibration images. For example, if all samples are imaged on a thin plastic plate, the X-Ray attenuation characteristics of the plate will be observed in every image acquired. However, the image quality could be improved by placing the plate in the X-Ray beam during calibration, since the characteristics of the plate are captured in the reference images.

A summary of the advanced calibration protocol for “Both” (Flat Field and Image) follows:

1. Place the uniform attenuator on the detector. The entire detector surface should be covered by the attenuator. For best results, the attenuator should be larger than the detector’s imaging area.
2. From the Home Menu screen press the **Reveal/Hide Top Menu** button [^] in the Right Tool Bar (the top-right corner of the screen).
3. Select “Tools”, “Advance Tools”, then “X-Ray Calibration/Exposure Setup”. The “X-ray Generator and Exposure Settings” Dialog opens. See figure below.




Figure 28: Tools Menu – Advanced Tools – X-Ray Calibration/Exposure Setup




Figure 29: X-Ray Generator and Exposure Settings Dialog Box




Figure 30: X-Ray Generator and Exposure Settings Dialog Box



Note

The values displayed in the example screens above are not necessarily correct for your system. kV and Exposure Time settings may vary from system to system. These settings may also change as the x-ray tube ages.



Note

The alphabetical references in the figures above correspond to the steps in the following procedure.

Refer to the figures on the previous page when performing the following procedure.

- a. Under “X-Ray Generator and Exposure Settings”, enter the generator kV setting in the appropriate field. This will normally be the Calibration kV value. (25.0 in the example shown.)
- b. Enter the “Magnification” level or leave the default setting of 1.000.
- c. Under “EXPOSURE SETTINGS” > “Set For”, click “Both”.
- d. Under “User Time/ADU Target” verify that “ADU” is selected. Leave the default “user ADU target” setting. (Typically 10,000).
- e. Leave the default “Test Exposure” time. 2.00 is typical.
- f. Click the **Acquire** button. The system will automatically begin an x-ray exposure and obtain an image. Note that this image is not displayed on the screen.
- g. When the exposure is complete, click on the **Estimate/Calculate** button.
- h. Record the calculated exposure time “Time [s]” value under “User Time/ADU Target”. This time is calculated by the software and the number is greyed out.
- i. Enter the calculated exposure time reading “Time [s]” in the “Test Exposure” time box.
- j. Press the **Acquire** button to begin another x-ray exposure and obtain an image.
- k. When the exposure is complete, click on the **Estimate/Calculate** button.
- l. Enter the calculated exposure time reading in the “Test Exposure” time box once again.
- m. Click on the **Accept Settings** button.
- n. Perform a standard calibration procedure, as described in Section 5.2.2, “Calibration Sequence”. *However, make sure you leave the attenuator plate in the chamber.*

Additional Notes on Advanced Calibration

Experienced users may wish to determine optimum “EXPOSURE SETTINGS” for “X-Ray Cal (Flat-Field)” and “Manual Image” independently, or “Both” (simultaneously), as previously described. By selecting the “Flat Field” radio button, the system will determine the optimum exposure time for acquiring Flat Field (calibration) images at the desired ADU level. selecting “Manual Image” will determine the optimum exposure time for the system with an object which will attenuate the X-Ray signal in the chamber. For best results, the entire imaging area of the detector must be covered by the object.

After the selection of either “X-Ray Cal (Flat-Field)” and “Manual Image” and the ADU level has been set, select the **Acquire** button. For the “Manual Image” selection, ensure the object or X-Ray attenuator is in the X-Ray beam. Press the **Acquire** button to generate X-Rays and obtain an image. When the X-Ray exposure is complete, click on the **Estimate/Calculate** button. The system calculates optimum exposure time, and the number is greyed out. Click the **Accept Setting** button to save the “Exposure Settings” to the system’s configuration file. **Note: Only the settings for the selected radio button are**

saved. Next perform a standard calibration procedure, as described in Section 5.2.2, “Calibration Sequence”.

5.3 Automatic Exposure Control (AEC)

The Faxitron Pro+ system utilizes Automatic Exposure Control (AEC) to provide the best possible image quality with a single press of the **Start X-Ray** button. The software analyses a *Preshot* x-ray exposure, typically 2-3 seconds, to calculate the optimal exposure kV and time settings needed to produce the best image quality. The Preshot AEC Exposure is configured at the factory and is turned on by default.

5.3.1 Preshot AEC Setup

“Preshot AEC Exposure Setup” allows the user to adjust the system for particular imaging requirements, or to compensate for system aging.

To access the “Preshot AEC settings” dialog go to the Top Menu Bar, select “Tools”, “Advanced Tools”, then “PRESHOT AEC Exposure Setup”.




Figure 31: Tools Menu – Advanced Tools – PRESHOT AEC Setup

The “Preshot AEC Settings” dialog will open.

Preshot AEC Settings

PRESHOT - AEC ON/OFF

☒ TURN ON PRESHOT - AEC

PRESHOT ANALYSIS METHOD

☒ DYNAMIC METHOD
(Acquire Preshot Image Before Every Image, Analyze and Tune Settings)

☐ STATIC METHOD
(Analyze the Current Loaded Image, Tune Settings)

DYNAMIC PRESHOT PREFERENCES

PRESHOT kV DESIRED [kV] 35.0

PRESHOT TIME DESIRED [s] 3.00

PRESHOT COMBINATION [1] NONE, DISCARD PRESHOT, SHOW AEC SHOT ALONE

STATIC PRESHOT PREFERENCES

PRELOADED IMAGE kV [kV] 34.0

PRELOADED IMAGE TIME [s] 5.05

PRELOADED IMAGE MEAN INTENSITY UNDER SAMPLE [ADU] 9320.7

☒ AUTO DETECT SAMPLE AREA AND ITS MEAN VALUE

EXPOSURE TUNING METHOD

☒ AUTO kV / AUTO TIME

☐ AUTO kV / FIXED TIME
ENTER TIME [s] 5.03

☐ FIXED kV / AUTO TIME
ENTER kV [kV] 45.0

THICKNESS/BEAM HARDENING /OVERBURN

☒ AUTO DETECT SAMPLE THICKNESS FROM PRESHOT/PRELOADED IMAGE

☐ ENABLE OVERBURN PROTECTION

ENTER SOFT TISSUE EQUIVALENT SAMPLE+TRAY THICKNESS [0] 0.00" = 0.00 cm

* Preshot and AEC shot are assumed to have similar sample+tray thicknesses

TARGET INTENSITY

ENTER DESIRED MEAN INTENSITY UNDER SAMPLE [ADU] 10000.0

* GENERAL RECOMMENDATION: Use 60-90% of Mean Flat Field Value For Satisfactory Results

PRESHOT ANALYSIS FOR STATIC PRESHOT METHOD

OPTIMUM EXPOSURE PARAMETERS: CALCULATE kV= [kV] TIME= [s]

ACCEPT CANCEL

Figure 32: Preshot AEC Settings Dialog

- Verify that “TURN ON PRESHOT – AEC” is selected. It should be on by default.
- Verify that “DYNAMIC METHOD” is selected.
- Enter “PRESHOT kV DESIRED [kV]” and “PRESHOT TIME DESIRED [s]”. Exposure time and kV shown will vary depending on the particular system. These values may also need to be adjusted as the system ages.

- Verify “AUTO kV / AUTO TIME” is selected. The software will calculate the best exposure kV and Time to provide optimum image quality based on the “Preshot” exposure. Alternately you may choose to select and enter a “FIXED” exposure time or kV. The software will then calculate the best exposure kV or time.
- Verify that “AUTO DETECT SAMPLE THICKNESS FROM PRESHOT/PRELOADED IMAGE” is selected.
- And finally, “ENTER DESIRED MEAN INTENSITY UNDER SAMPLE (ADU)”. This is typically 10000 for the Faxitron Pro+ system. However, this value may vary depending on the type of sample being imaged.
- Click on the **ACCEPT** button.



Note

If the “STATIC METHOD” is used instead of the “DYNAMIC METHOD” the software will use Hologic predefined parameters.



Note

The user may also deselect “AUTO DETECT SAMPLE THICKNESS FROM PRESHOT/PRELOADED IMAGE” and enter the sample thickness. This may result in more accurate images.

Place the sample in the chamber and press the **Start** button in the VisionLS software (or on the Control Panel). The VisionLS software will perform a pre-shot exposure, create a Region of Interest (ROI) for obtaining the exposure parameters, then take the image exposure.

5.4 Sample Data

Sample data must be entered into the Sample Information Editor to initiate a procedure.

Clicking the **Start Procedure** button, *see figure below*, opens the Sample Data Entry page which provides methods to enter sample data. *See figure 34*. If a Work List has not been configured, the option will be greyed out and only the **Manual Entry** and **Select from Database** options will be accessible. These options are explained in the next sections.




Figure 33: Home Menu




Figure 34: Sample Data Entry Options

5.4.1 Manual Entry

Also refer to Section 4.1.4.

1. Click on the **Manual Entry** button to manually enter a new Procedure. The “Sample Information Editor” will open.

The figure consists of two side-by-side screenshots of the 'Sample Information Editor' dialog box. The left screenshot shows the 'Blank' state where the 'Operator Name' (with sub-fields 'First' and 'Last'), 'Sample ID', and 'Accession No.' fields are marked with an asterisk (*) to indicate they are required. These three fields are enclosed in a red rectangular box. The 'Accept' button at the bottom is greyed out. A red arrow points to the 'Accept' button with a text box stating: 'The Accept button is Not Active (greyed out)'. The right screenshot shows the 'Complete' state where the same fields are filled with 'A', 'C', and 'D'. A red box highlights a message: 'All required fields are filled in. The Accept button is Active (Not greyed out)'. The 'Accept' button is now active and blue. Both screenshots have a bottom bar with a left arrow and the text 'Start a New Sample or Procedure'.

Figure 35: Sample Information Editor Dialog Box – Blank on Left, Complete on Right

2. Complete all required fields, which are indicated by an asterisk. Additional fields may be completed at the users' discretion. The **Accept** button will be greyed out (inactive) until the required entries are completed.



Note

If the date is entered it must be entered as yyyy/mm/dd (e.g., 1962/08/13).



Note

The “required fields” may be altered or deleted depending on the application. Contact Hologic support for further information.

3. Click the **Accept** button. The Image Acquisition screen will open.

5.4.2 Continue with Current Procedure

The **Continue with Current Procedure** button is not active in the VisionLS software at this time.

5.4.3 Select from Work List (Remote Modality Work List)

Work Lists are normally only utilized in clinical applications. However, if a Work List has been configured the user can click the **Select from Work List** button which will display a page similar to the Select from Database page. *See figure 36.*

The user can scroll through the list to find the desired procedure, or can use the *Search Parameters* window, on the left side of the page, to narrow down the list.

Click on the desired procedure then click on the **Submit** button, or simply double click on the procedure. This will open the Sample Information Editor Dialog Box shown in *figure 35.*

The Operator Name and Sample ID fields will be populated. Enter the Accession No. The **Accept** button will now be active.

Click the **Accept** button to open the Image Acquisition screen.



Note

The user may enter additional sample information as desired.

5.4.4 Select from Database

Clicking the **Select from Database** button will display the following page. *See figure on next page.*

Database

Sample ID	Operator Last Name	Operator First Name	Time Modified	Date of Birth	Sex
1	1	1	2018/04/18 15:24...	2016/02/12	
1-A	JONES	A	2021/03/08 16:24...		
102	EYE	MOON'S	2018/02/28 16:24...	2018/02/28	
104	SHELL	TULIP	2018/02/28 16:30...	2018/02/28	
11	HOWES	SCULPIN	2018/02/28 15:26...	2018/02/28	
111	SHELLS	GROUP	2018/02/28 16:48...	2018/02/28	
2	HOWES	SNAKE	2018/02/28 14:44...	2018/02/28	
2-A	JONES	A	2021/03/08 16:29...		
24	KINGSLEY	PUFFER	2018/03/01 15:08...	2018/03/01	
3-A	GOOD	JOHNY	2021/03/08 16:38...		
47	SHELL	THIN	2018/03/02 12:31...	2018/03/02	
5	HOWES	FISH	2018/02/28 15:00...	2018/02/28	
8282_3 4	FOXFO	EDA CRISPR	2018/03/05 13:24...	2018/03/05	
SEAROBIN1	SEAROBIN1	SEAROBIN1	2018/02/28 11:46...	2018/02/28	
SH	SF	W	2013/04/25 15:25...	2013/04/25	

Figure 36: Select from Database Page

1. Select the sample from the table on the right side of the page.
2. Click on the **Accept** button (or double-click the sample entry in the table). The Sample Information Editor will open. *See figure 35.*
3. Complete the required fields, then click the **Accept** button. *Refer to Section 5.4.1, Manual Entry.*

The “Search Parameters” window allows the user to narrow the sample list by entering specific information in one or more fields. Only those samples matching all search parameters will be displayed. *See figure below.*

Sample ID	Operator Last Name	Operator First Name	Time Modified	Date of Birth	Sex
1	1	1	2018/04/18 15:24...	2016/02/12	
1-A	JONES	A	2021/03/08 16:24...		

Figure 37: Select from Database – Using Search Parameters

5.5 Image Acquisition

X-rays are started by pressing the green **X-Ray** button on the front panel touch screen, or clicking the **Start X-Ray** button in the lower right corner of the screen. *See Section 4.1.5, Acquiring Images.* The system will take a short “sample” x-ray, or pre-shot (which is explained in Section 5.3) followed by a longer exposure to acquire the image. The software will process the image before it is displayed on the monitor screen. Once the image is displayed, buttons on the left and right side of the screen, *covered in Section 5.6, Tool Bars and Menus*, will be enabled.

5.6 Tool Bars and Menus

VisionLS software provides various tools and menus to control various image attributes or parameters, x-ray acquisition, image views, image tools, and DICOM settings/access. The tools most frequently used are in Tool Bars at the sides and bottom of the screen. The Tool Bars and Menus are described in more detail in the following pages.

5.6.1 Left Tool Bar

VisionLS software provides single click access to commonly used tools or functions with a vertical tool bar located on the left side of the screen. More detailed descriptions follow in the next pages. *See figure below.*




Figure 38: Left Tool Bar

View/Change Sample Information

Click the **View/Change Sample Information** button to view or change the sample information after it has been input. *See figure below.* Note that sample information can only be changed before an image has been taken. The button will turn green when the mouse is hovered over it and while the Sample Information Editor is open.



Figure 39: View/Change Sample Information Button


The Sample Information Editing Dialog box will open. This editing window will be populated with all the information that was previously entered/selected. Note that the dialog box will now include a **Cancel** button. *See figure on left below.*

After acquiring an image, the user may view, but not edit, the sample information. The Sample Information Dialog will now only have options to **Close Dialog** or **Start a New Sample or Procedure**. *See figure on right below.*

The figure shows two versions of the 'Sample Information' dialog box. Both have a green header and a white body. The left version, labeled 'Before Image', includes fields for Operator Name (A, B), Sample ID (C), Accession No. (D), Date of Birth (yyyy/mm/dd), Sex (dropdown), Procedure Name (dropdown), Laterality (dropdown), Position View (dropdown), Body Part (dropdown), Refer. Physician (text), Technician (text), Technician ID (text), Study Name (text), and Study ID (text). At the bottom, there is a blue 'Accept' button with a right arrow, a grey 'Cancel' button with a left arrow, and a grey 'Start a New Sample or Procedure' button with a left arrow. The right version, labeled 'After Image', has the same fields but replaces the 'Accept' button with a blue 'Close Dialog' button with a right arrow. The 'Cancel' button is replaced by a grey 'Start a New Sample or Procedure' button with a left arrow.

Figure 40: Sample Information Dialog Box - Imaging Mode: Before Image (Left), After Image (Right)

Clicking the **View/Change Sample Information** button while in Review Mode allows the user to view the sample information. The Sample Information Dialog will now only have the **Close Dialog** button.



The dialog box is titled "Sample Information" and has a green header bar. It contains three sections: "Sample Information", "Procedure Information", and "Study Information". The "Sample Information" section includes fields for Operator Name (A, B), Sample ID (C), Accession No. (D), Date of Birth (2021/05/05), and Sex (dropdown). The "Procedure Information" section includes fields for Procedure Name (dropdown), Laterality (dropdown), Position View (dropdown), Body Part (dropdown), Refer. Physician (text), Technician (text), and Technician ID (text). The "Study Information" section includes fields for Study Name (text) and Study ID (text). At the bottom is a blue bar with the text "Close Dialog" and a right arrow.

Figure 41: Sample Information Editing Dialog – Review Mode

Enhance Image

The **Enhance Image** button will process the image on the display to emphasize the smaller details, such as micro-calcifications, and small masses in the image.



Figure 42: Enhance Image Button

The algorithm, commonly known as “Peripheral Equalization” or “Thickness Equalization” in advanced imaging, is applied to the whole image. In addition to this algorithm, Hologic employs a second, internally developed algorithm called the “Calcification Emphasis” algorithm. This secondary algorithm was created to bring out small details in the image and make them much easier to see and analyze.

The system will be idle (in Standby) during the processing. Once complete the button will stay depressed/active. The user can undo the process by pressing the button again.




Note

The system may have shipped to automatically process “Image Enhancement” after every image. In that case pressing the button will load the “Un-Enhanced” image.

Window Level and Contrast Adjustment and Inverting the Image

The Window Level and Contrast Adjustment tool allows the user to adjust the image's window level and contrast, restore the images original contrast, or invert the image. Each is described below:



Use a pair of sliders (arrows) to adjust the window level and contrast. Clicking between the sliders and moving up or down will maintain the window width (range) while moving the center of the window up or down. Moving the sliders individually will adjust just the upper or lower window levels individually.

You may also type upper and lower window levels into the white text fields.

The **Original Contrast** button restores the images original contrast.

The **Invert** button flips the grayscale values, i.e., converts the image from a positive to a negative, or vice versa.

Figure 43: Window Level and Contrast Adjustment Tool

Click the **OK** button to accept the changes. The **Cancel** button restores the previous settings.

Built-in Windowing Tool

The user may also adjust the Window Level and Width using mouse movements. This "Built-in Windowing Tool" can be used on any image, whether just acquired or opened in review mode.

Position the mouse over the displayed image then click and drag to adjust the brightness and contrast. Moving the mouse up and down adjusts brightness, left and right adjusts contrast. Release the mouse button when satisfied with the results.

The figures below show the same image with different window levels. Note the position of the sliders and the corresponding upper and lower limits displayed in the text fields above and below the sliders. Upper/lower limits of 13191/1912 versus 10649/5223.




Figure 44: Window Level and Contrast Adjustment Tool




Figure 45: Window Level and Contrast Adjustment Tool

An example of inverting is shown below.





Figure 46: Example of Inverting: Inverted on the Left

Digital Zoom Tool Button

The **Digital Zoom Tool** button opens the Zoom Adjustment tool, which provides options for adjusting an image's zoom level.



Use the slider to manually adjust the zoom level from 0.10x to 3.00x

The **Fit to Specimen** button will attempt to detect and zoom in onto the specimen to ensure the sample is in full view. If a specimen is not detected, it will default to the previous mode (either “Fit to Screen” or “Zoom 1.0x”)

Fit to Screen fits the image to the screen.

Zoom 1.0x or “Full Resolution”. The image is adjusted for a 1:1 ratio with the monitor resolution.

Figure 47: Digital Zoom Tool Button and Dialog Box

Zooming in on an image will require panning to navigate around the image. A thumbnail will appear in the top left corner of the screen. A blue box in the thumbnail indicates the portion of the image displayed on the screen. Click in the blue box and move it around the thumbnail to move the image in the main screen.

Traditionally, this feature is also called “Pan View” and/or “Pan and Scan”.

The default mode is “Fit to Screen”. The actual zoom level will vary depending on the size of the sample.



Note

Depending on the size of the specimen there may be minimal difference observed between **Fit to Specimen** and **Fit to Screen**.

The following images show the same sample at different zoom levels. The top figure is at the default zoom level. The bottom level is zoomed in to 1.00x.




Figure 48: Image in Default Zoom Level




Figure 49: Same Image of Snake at 1.00x Zoom

Note that when the image no longer fits on the screen the “Pan Window” opens which allows the user to move the sample around the screen.

Annotation Tools Button

The **Annotation Tools** button opens the annotation tools Dialog. These tools and their functions are discussed in detail in Section 5.10, User Annotations.




Figure 50: Annotation Tools Button and Dialog Box

PACS Store

The **PACS Store** button provides the ability for single click PACS storage. PACS Store only functions when the system is installed in a medical facility with a PACS network. To save files locally, please refer to Section 5.9, DICOM Functionality.



your

Figure 51: PACS Store Button

PACS Print

The **PACS Print** button provides the ability for single click print to a preconfigured PACS printer. PACS print will only function when the system is installed in a medical facility with a PACS network. Refer to Section 5.9, DICOM Functionality. For local printing please refer to Section 5.6.4, Top Menu Bar > Tools Menu > Print Current Image, and Section 5.8.2, Database Top Menu Bar > Additional Options > Print.



Figure 52: PACS Print Button

Exposure Control

The **Exposure Control** button opens the Exposure Mode Dialog box, which allows the user the convenience of controlling the X-Ray system through the software. The user can select the Exposure Mode, adjust X-Ray Tube kV and Exposure Time, adjust the number of images to acquire for Image Averaging (from 1 to 8), and initiate the X-Ray exposure.




Figure 53: Exposure Control Button and Dialog Box

The status, or availability, of the Exposure Settings is indicated by the up/down arrows adjacent to the setting. The arrows will be white when the setting is active and light blue when inactive. The Exposure Settings that are active is dependent on which Exposure Mode is selected.

- X-Ray tube current (mA) is not adjustable via software in the Faxitron Pro+ system and will always be inactive (light blue).
- X-Ray tube kV and Exposure Time settings are active/in-active dependent on which Exposure Mode is selected. KV will only be active in “Auto Time” and “Full Manual” modes. The Exposure Time will only be active in “Auto kV” and “Full Manual” modes.
- Image Averages is always active.

The figure above shows that “Full Auto” mode has been selected and that kV and Time are inactive. Only Image Averages is active in Full Auto Mode.

When finished click on the **Request X-Rays** button to acquire an image, or the **Set** button to save the settings and close the Exposure Control Dialog. All x-ray exposures initiated with the **Start X-Ray** button (either the physical or software button) will be taken with the saved settings.

Clicking the **Cancel** button will close the Exposure Control Dialog and discard any changes.

Back to Main Menu

The button to the right exits the image acquisition mode and returns to the Home Menu screen



Figure 54: Back to Main Menu Button

Undocking a Dialog

When a dialog is opened, it is “docked” to the left side of the screen. A user may un-dock the dialog box by clicking on the “Pin” icon in the upper-right corner of the box. They can now drag the dialog box anywhere on the screen. When the dialog box is undocked the “pin” then turns into an “X” icon allowing for quick closing of the dialog.



Figure 55: Pin (Un-Docking) Icon

5.6.2 Right Tool Bar

The Right Tool Bar provides the user with one click ability to reveal or hide the Top Menu, reveal or hide current image information, toggle annotations on or off and turn the sample viewer (x-ray chamber webcam) on or off. Each button is described below.




Figure 56: Right Side Menu Bar

Reveal/Hide Top Menu Button

The top menu is hidden by default when the VisionLS software is open. Clicking on the **Reveal/Hide** button brings the top menu into view. The top menu is hidden when the

button is pointing down and revealed when the button is pointing up. This button is available in all screens and will be white or black, depending on the background.




Figure 57: Reveal/Hide Top Menu Button

The Top Menu is discussed in detail in Section 5.6.4.

Image Information Button

Clicking on the **Image Information** button will either display or hide the current sample image information, depending on the current state. When the button is orange, *as shown below*, the image information box is revealed. The box is hidden when the button is white.

The Image Information box contains information such as the magnification level the image was taken at, date, time, and day the image was taken, dosage information, file name of the image and the location of the image. It will also contain the name of the institution if it was entered.




Figure 58: Image Information Button and Box

Toggle Annotations Button

The **Toggle Annotations** button will turn any annotations that have been made on an image On or Off. Turning the annotations Off temporarily hides them, it does not delete them from the image. The **Toggle Annotations** button will display “ON” or “OFF” as appropriate.



Figure 59: Toggle Annotations Button

Sample Viewer

The sample viewer allows the user to monitor samples in the x-ray chamber. Clicking the **Sample Viewer** button turns a webcam, located at the top of the chamber, On or Off. When the Sample Viewer is turned on a live view of the chamber is displayed in a small window in the upper left corner of the image acquisition screen. The **Sample Viewer** button will turn orange when the mouse is hovered over it.



Figure 60: Sample Viewer On-Screen Window

5.6.3 Bottom Tool Bar

The Bottom Tool Bar will vary depending on whether the system is in Imaging or Review Mode.

The left third of the Bottom Tool Bar (see A in figure below) will always show current sample information.




Figure 61: Bottom Tool Bar – Imaging & Review Modes

The middle third of the Bottom Tool Bar (see B in above figure and next figure) will vary between Imaging and Review modes.




Figure 62: Bottom Tool Bar – Middle Section (B)

Refer to the above figure for the following descriptions:

- **1 & 3 - Left and Right Arrow** buttons (1 & 3 above): Open either the “Previous” or “Next” image in a procedure that has multiple images.
- **2 - Browser** button: Opens the “Current Sample Images” view, which is a thumbnail strip showing all the images in the current sample.
- **4 - Add Images to Review** button flags images in the thumbnail strip for review. The number in the orange circle indicates how many images have been flagged. Only available in “Image Mode”.
- **5 - Eye** button: Will be displayed while in the “Image Acquisition Mode”. Click this button to enter “Review Mode”.
- **6 - Change Review Layout** button: Pulls up a menu to select the review pane configuration. Only available in “Review Mode”.
- **7 - Radiation Symbol** button: Will be displayed while in “Review Mode”. Click this button to return the system to the “Image Acquisition Mode”.

The right third of the Bottom Tool Bar contains the system status bar and magnification info.



Figure 63: Bottom Tool Bar – Review Mode, Right Third

- The “Magnification” icon displays the mag tray position while in “Image Acquisition” mode. *See figure above.*
- System status will always display “System in Standby” while in the “Review” mode. *See figure above.*
- The system status bar turns into the **Start X-Ray** button when the system is in “Image Acquisition” mode. When an x-ray is initiated, the button will turn into a **Stop X-Ray** button. When the exposure is complete the button turns into a **Processing** status bar. Once the image is processed and displayed the button/status bar will return to the **Start X-Ray** button. *See figure on next page.*



Figure 64: Bottom Tool Bar – Buttons/Status Bar – Imaging Mode

5.6.4 Top Menu Bar

The software provides access to “View”, “Tools”, “Image Tools”, “DICOM Utilities” and “Help” menus through the Top Menu Bar. See figure below. The Top Menu Bar is hidden by default and is displayed by clicking the **Reveal/Hide Top Menu** button in the Right Menu Bar. See figures 56 & 57. Each menu item is covered in detail in the following pages.



Figure 65: Top Menu

FILE MENU

This menu contains options commonly found in a Windows File Menu.




Figure 66: File Menu

EDIT MENU

This menu contains options commonly found in a Windows Edit Menu. However, all of the options are greyed out and unavailable to the user.




Figure 67: Edit Menu

VIEW MENU

The “View” menu provides the user the ability to modify the acquired images’ display format or characteristics. Each of the “View” menu items has sub-menus. Each menu item is described below.




Figure 68: View Menu

Zoom Mode




Figure 69: View Menu - Zoom Mode Options

- **NORMAL (FULL RES.)**

Displays the acquired image in its native format. In full resolution the image in the Faxitron Pro+ system exceeds the monitors display area. A thumbnail window will be displayed in the top left corner of the image which allows the user to “pan” the image. Refer to Digital Zoom Tool in Section 5.6.1. The same function is provided by the **Zoom 1.0x** button in the Digital Zoom Tool.

- **ADJUST ZOOM**

This option is greyed out and not available at this time.

- **FIT TO SAMPLE**

This option adjusts the acquired image to best fit the specimen sample within the monitors display area. This option provides the same function as the Digital Zoom **Fit to Specimen** button.

- **FIT ALL**

This option displays the entire acquired image, regardless of its native format. The Full Resolution options does not work in this view mode. Provides the same function as the Digital Zoom **Fit to Screen** button.

- **FIT WIDTH**

This option displays the acquired image across the available width of the entire display area.

- **FIT HIEGHT**

This option displays the acquired image across the available height of the entire display area.

Image Averages

This option allows the user to adjust the number of image averages to be taken per an exposure. The available options are 1x, 2x, 4x, or 8x shots. Default setting is 1x SHOT. If the user selects 2, 4, or 8 shots the software will acquire the number of images selected, one after the other, with a single push of the **Start** button, then average those images into a single image.




Figure 70: View Menu - Image Averages Options



Note

Image Averages may also be set in the Exposure Control Dialog Box accessed through the Side Tool Bar.

Image Perspective

This option allows the user to rotate the displayed image to 0, 90, 180, and 270 degrees.




Figure 71: View Menu – Image Perspective



Note

The default Image Perspective for the Faxitron Pro+ system is 270 degrees.

Grid Options

Grid Options allows the user to superimpose a “Grid” over an image and select the size of that grid. Available grid sizes are 64, 128, and 256 pixels.




Figure 72: View Menu - Grid Options




Figure 73: Grid Sample, Pitch Set to 256

User Window Levels

This allows the user to define up to three (3) pre-defined, user set Window Levels (WL).




Figure 74: View Menu – User Window Levels

- Take an image, or open an image from the database,
- Manually Adjust the image to the desired level window level.
- Select the “View” menu, then “User Window Levels”, then “RECORD CURRENT WL AS”. Then select User WL #1, #2, or #3”.
- The adjusted WL will be saved under the choice selected.

The user can now apply that saved window level to any other image by click on the “View” menu, then “User Window Levels”, then “APPLY USER WL #X”.

Default Window Level

There are 3 menu options under Default Window Level:

- “OPTIMIZE FOR WHOLE AREA” adjusts the window level for subsequent images for the entire image array.
- “OPTIMIZE FOR SAMPLE AREA” adjusts the window level for subsequent images to be optimized for the sample only. It ignores the image area outside of the sample. “Optimize for Sample Area” is the default setting.
- “Set Default Window Level to Current” applies the window level for the current image to all subsequent images. Adjust the window level for the enhanced image then select “Set Default Window Level to Current”. Repeat for the un-enhanced image. All subsequent images will now have those approximate window levels applied.




Figure 75: View Menu – Default Window Level

Default Aspect of Raw Files

“Default Aspect for Raw Files” helps when importing/opening previously saved image files with unknown file attributes, or binary image data without any header information. Since the exact imager array size is not available the software will import/open the image with “best fit” imager dimensions based on width being larger or smaller than height (default aspect information).




Figure 76: View Menu – Default Aspect for Raw Files



Note

These menu options are greyed out. Default setting is “Width >= Height and can’t be changed by the user.

Enhance Image (same as the Enhance Image Button)

Processes current Sample image and enhances visualization of tissue located near the periphery of the sample.

This menu option acts the same as the **Enhance Image** button in the Left Side Tool Bar.




Figure 77: View Menu – Enhance Image

Sample Viewer

The sample viewer allows the user to monitor samples in the x-ray chamber. *Refer to Sample Viewer in 5.6.2, Right Tool Bar.*

Magnifier Size

“Magnifier Size” is used to adjust the size of the *Magnifier* under the “Image Tools” menu. There are 3 available settings: Small, Medium, and Large.




Figure 78: View Menu – Magnifier Size

Review Layout

“Review Layout” is used to select the configuration of the image panes in Review (Database) Mode. See figure below for available options.

This same function is available using the **Change Review Layout** button in the Bottom Tool Bar.




Figure 79: View Menu – Review Layout

An example of the (1 + (2 X 1)) Triple layout is shown below.




Figure 80: View Menu – Review Layout (1 + (2 X 1)) Triple



Note

See Section 5.8, Database/Review Mode for additional information on using multiple review panes.

TOOLS MENU

The “Tools” menu provides the user with a number of tools is located in the Top Menu Bar, found in the upper left-hand corner of the screen. These tools apply to the image currently displayed.




Figure 81: Tools Menu

Preferences

The “Preferences” menu opens a configuration management utility which can be used to change various options or parameters.



Note

The “*Preference*” tool option is password protected and should only be accessed by an authorized Hologic service representative.



Caution:

Changing configuration parameters without proper knowledge or authorization can cause loss of image quality, system damage, and system malfunctions.

Service Mode

“*Service Mode*” is password protected and turned on by Hologic authorized technicians for troubleshooting purposes.

Service Tools

“*Service Tools*” is normally greyed out. It is only turned on when the system is in “*Service Mode*”.

Save Reference Files

“*Save Reference Files*” is normally greyed out. It is turned on when the system is in “*Service Mode*”.

Advanced Tools Options

The Advanced Tools menu offers various options to configure aspects of the software such as post-processing and configuration options.




Figure 82: Advanced Tools Menu

- **Auto Save On/Off**

“*AUTO SAVE ON/OFF*” allows the user to turn the auto save feature on or off, as desired. Menu options are “On” or “Off”.

When “Auto Save” is “On”, all subsequently acquired images are saved to the database. The files are written to a Database Directory defined in the configuration file.

“Auto Save” is set to “On” by default.

- **Edge Suppression**

Edge suppression is used to minimize the visual effects of the edges of CMOS sensors in the image. Menu options are “On” or “Off”.

Edge suppression is not used on the 43x43cm detector.

- **Tile Suppression**

Tile suppression is used to minimize the visual effect of the individual CMOS sensors in the image. Menu options are “On” or “Off”.

Tile suppression is not used on the 43x43cm detector.

- **X-Ray Calibration/Exposure Setup**

“X-Ray Calibration/Exposure Setup” aids the user in determining the calibration settings, proper exposure time for a selected kV setting for best image quality, or proper kV setting for a selected exposure time for best image quality. Please refer to Section 5.2.3, Advanced Calibration, for a more detailed description.

- **PRESHOT AEC Exposure Setup**

“Preshot AEC Exposure Setup” allows the user to adjust the system for particular imaging requirements or to compensate for system aging. Please refer to Section 5.3.1, Preshot AEC Setup, for a more detailed description.

- **IMAGE ENHANCEMENT Setup**

“Image Enhancement Setup” allows the user to adjust various parameters of the image enhancement algorithm.




Figure 83: Advanced Tools – Image Enhancement Setup

Record the settings before making any adjustments so you can return the system to the default settings if needed.

- **Login Utility**

The Login Utility allows the user to configure and activate *start-up password authentication* for VisionLS software. Activation of this feature is optional on the

Faxitron Pro+ system. If the user activates the *password authentication*, they will be prompted for the password each time the program is opened. The program will close if the correct password is not entered.

To activate *password authentication*:

- Click on “Login Utility” under “Tools” > “Advanced Tools”. The Login Utility Dialog will open. *See figure on next page.*




Figure 84: Login Utility Dialog Box

- Enter the Current Password. *Your authorized Hologic representative can provide you with your systems password.*
- Create a new password and enter it in the “NEW PASSWORD” and “RE-ENTER NEW PASSWORD” spaces.
- Click the checkbox for “Enforce User Authentication on Software Startup”.
- Click the **Apply Changes** button, then click the **Exit** button.

The system will now require the password each time VisionLS software is started.



Caution:

You can activate *password authentication* using your systems default password. However, Hologic highly recommends creating a new password for enhanced system security.

- **Cyclic Automated Acquisitions**

Cyclic Automated Acquisitions allows the user to automatically acquire sequential images with a user defined delay time between each image.




Figure 85: Advanced Tools – Cyclic Automated Acquisitions

- Default settings are 10 cycles with a 2.0-minute delay between cycles.
- Maximum delay between cycles is 1440 minutes.
- Maximum number of cycles is 1000.

Enter the desired delay time and the number of cycles. Click the **UPDATE** button, then click the **START CYCLIC EXPOSURE** button. The system will begin the exposures. Each exposure will include all processing as if taking a single exposure. The sequence is as follows:

1. System will take a “*preshot*” exposure,
 2. System will take the “*image*” exposure,
 3. Unenhanced image will be displayed,
 4. Software will process then display the enhanced image,
 5. Software will start timing the “*delay between cycles*”,
 6. When the defined delay time has elapsed, the system will start another “*cycle*” at step one.
 7. The process will be repeated until the defined number of cycles has been completed.
- All the images taken are saved and available for review from the Current Image Browser or the Database.

- **OTHER ADVANCED User Tools**

All menu items except “Multi Scale Enhancement” are greyed out and not available to the user.




Figure 86: Advanced Tools – Other Advanced User Tools – Multi Scale Enhancement

- **Multi Scale Enhancement** equalizes image details at multiple scales in the image, resulting in effective dynamic range compression and sharpening of faint image details.

PROPERTIES

“Properties” hides or shows the Image Properties Box.




Figure 87: Advanced Tools - Properties

EDIT Current Image Record

Currently not used.

DELETE Current Image Record

Deletes the current image from the file directory as specified in the configuration file.

ADD Current Image Record

Adds the current image to the file directory as specified in the configuration file.

CAPTURE Screen Shot JPG/BMP

“CAPTURE Screen Shot JPG/BMP” allows the user to capture and save a screen shot of the entire display. Clicking this menu option will open a “Save As” dialog. The dialog is populated with default/automated entries. The user can specify a different “Save in” directory, can create their own “File Name”, choose between JPEG and BMP file types, and select the “Sub Type”. *See figure on next page.*




Figure 88: Tools Menu

Save/Store current image

“Save/Store current image provides the user options for saving the current image:

- LIS Archive Image,
- As local DICOM image,
- As local image, and
- More options.




Figure 89: Tools – Save/Store Current Image

1. Selecting “LIS Archive Image” opens the dialog shown below.




Figure 90: LIS Save File As Dialog

- a. Click the **Browse** button to navigate to the directory you want the image saved in.

- b. Select the image format from the drop-down menu. Choices are: JPG, PNG, and TIF-8. JPG is selected by default.
 - c. Leave the “Selection Filter” at “All Images” or choose “Enhanced Images Only”, or “Un-Enhanced Images Only” from the drop-down menu.
 - d. Enter a “Custom Prefix” if desired.
 - e. Click the **Save** button. The displayed image will be saved to the specified directory.
2. Selecting “As Local DICOM Image” opens a “Save as DICOM File” dialog as shown below.




Figure 91: Save as DICOM File Dialog

- a. Navigate to the desired directory,
 - b. Enter a “File Name”,
 - c. Click the **Save** button. The image will be saved as a DCM file in the selected directory.

3. Selecting “As Local Image” opens the “Save As” dialog shown below.




Figure 92: Save As Dialog

- a. Navigate to the directory you want the image saved in,
 - b. Leave the system generated “File Name” or enter a “File Name” or your choice.
 - c. Select the “Save as type” from the drop-down menu. Options are: Win BMP (*.bmp), JPEG (*.jpg), TIF (*.tif), or RAW (*.raw). Default file type is BMP.
 - d. Leave “Sub Type” as “Uncompressed” or choose “RLE Compressed” from the drop-down menu.
 - e. Click the **Save** button.
4. Selecting “More Options” opens a “DICOM Networking Dialog”. See Section 5.9, *DICOM Functionality*, for more information.

Print current image

Prints current image to the WINDOWS printer.




Figure 93: Print Current Image

Selecting “To local printer” opens the Windows Print Dialog.

Selecting “More options...” opens a DICOM Networking Dialog box. See Section 5.9, *DICOM Functionality*, for more information.

IMAGE TOOLS MENU

The “Image Tools” menu provides the user tools to aid the analysis of the image. Click on any menu item to activate the tool. Click again to turn the tool off.




Figure 94: Image Tools Menu

Original Contrast

“Original Contrast” restores the images contrast levels. Behaves the same as the **Original Contrast** button in the Window Level and Contrast tool. See Section 5.6.1, *Left Tool Bar – Window Level and Contrast Adjustment and Inverting the Image*, and figure 43.

Invert

Clicking on “Invert” flips the grayscale values of the image, i.e., converts the image from a positive to a negative image, or vice versa. Click “Invert” to restore the image to its original state. Behaves the same as the **Invert** button in the Window Level and Contrast tool. See 5.6.1, *Left Tool Bar – Window Level and Contrast Adjustment and Inverting the Image*, and figures 43 & 46.

Magnifier

The “*Magnifier*” allows the user to magnify an area of interest, acting like a *Magnifying Glass*. Click on “*Magnifier*” to turn the feature on. Position the mouse cursor over the desired area to be magnified and click the left mouse button. The area of interest will be magnified. Click “*Magnify*” again to turn off the feature. *See figure below*. The size of the magnified area of interest can be adjusted by the user. *Refer to the “Magnifier Size” sub-menu in Section 5.6.4, Top Menu Bar - View Menu and figure 77.*



Figure 95: Image Tools Menu - Magnifier

Enhance

Clicking “*Enhance*” will create a sharper definition around the edges of regions with a contrast difference. This “*Edge Enhancement*” may amplify the noise in the image. Selecting the menu item again will turn off Edge Enhancement. Clicking on “*Original Contrast*” will also restore the image back to its original state. Note this is a different feature than the Enhance Image feature in Section 5.6.1, Left Tool Bar – Enhance Image and 5.6.4, Top Menu Bar – View Menu – Enhanced Image.

ROI Contrast

“*ROI Contrast*” will optimize the contrast levels within a user-defined region of interest (ROI). After selecting “*ROI Contrast*” click and hold the left mouse button at the starting point of the desired ROI, drag and release the mouse button at the desired ending point of the ROI. The contrast will automatically adjust to the best levels for that ROI. *See figures next page.*




Figure 96: Image Tools Menu – ROI Contrast




Figure 97: Image Tools Menu – ROI Contrast

ROI Window Level

The “ROI Window Level” function is not available in the current version of VisionLS software.

Ruler

“Ruler” allows the user to measure the distance between two points in the displayed image. Length measurements are in millimeters (mm). After selecting “Ruler” click and hold the left mouse button at the beginning of the point of measurement. Drag and

release the mouse button at the end point of the measurement. The Ruler is displayed as a red line on the image with the measurement displayed at the end of the ruler.




Figure 98: Image Tools Menu - Ruler

The Ruler remains active until the user selects another image tool or until another action is performed. Click “Ruler” again to turn off the feature.



Note

Hold the keyboard SHIFT key while drawing the ruler to get a perfectly horizontal, vertical, or 45-degree ruler.

ROI Statistics

Selecting “ROI Statistics” will enable the user to obtain the image statistics within a user defined ROI. Click and hold the left mouse button at the starting point of the desired ROI, drag and release the mouse button at the desired ending point of the ROI. A box will open displaying the images Statistical information. Statistics for the selected ROI are displayed on the left, and on the right for the whole image. The stats displayed are: beginning and end points of the ROI, Minimum, Maximum and Mean intensity values, and Standard Deviation (STD). To close the Statistics information box, simply click on the **Close** button. See figure on next page. Click “ROI Statistics” again to turn off the feature.




Figure 99: Image Tools Menu - ROI Statistics

ROI Histogram

Selecting “ROI Histogram” allows the user to create a histogram for a user defined ROI. Click and hold the left mouse button at the starting point of the desired ROI, drag and release the mouse button at the desired ending point. A “Graph Dialog” will open displaying Histograms for the selected ROI.




Figure 100: Image Tools Menu - ROI Histogram

The Black curve is the Histogram of the selected region. The Blue curve is the Histogram of the whole image. Both curves are normalized against the peak pixel value of the selected region of interest. Click “ROI Histogram” again to turn off the feature.

ROI Profile

Clicking “ROI Profile” allows the user to create a profile for a user defined ROI. Click and hold the left mouse button at the starting point of the desired ROI, drag and release the mouse button at the desired ending point of the ROI. A “Graph Dialog” will open displaying the profile of the selected ROI. By default, the averaging direction of the selected ROI is the short axis. Click the **SWITCH, HOR<->VERT** button to toggle the axis back and forth. Click “ROI Profile” again to turn the feature off.




Figure 101: Image Tools Menu - ROI Profile

Line Profile

Clicking “Line Profile” allows the user to create a profile for a user defined line. Click and hold the left mouse button at the starting point of the line, drag and release the mouse button at the desired end point. A “Graph Dialog” will open displaying a profile of the drawn line. *See figure on next page.*

The start and end point of the drawn line are displayed below the graph.

The starting point of the line will always be displayed on the left side of the graph.

Click “Line Profile” again to turn off the feature.




Figure 102: Image Tools Menu – Line Profile

DICOM Utilities Menu

There are two (2) drop down options in this menu:

- DICOM MODALITY Setup, and
- DUMP DICOM Header




Figure 103: DICOM Utilities Menu

DICOM Modality Setup is discussed in Section 5.9, DICOM Functionality.

Selecting “DUMP DICOM Header” creates a text file of the images DICOM header.

Help Menu

The Help Menu provides options to *open the Faxitron Pro+ User's Manual* and to view information about the VisionLS software and related components.




Figure 104: Help Menu




Figure 105: About Vision

5.7 Current Sample Images Browser

The **Current Sample Images Browser** button is located in the middle of the bottom tool bar. Clicking the **BROWSER** button opens a thumbnail strip of current sample images. *See figures below and on next page.*




Figure 106: Current Sample Images Browser Button




Figure 107: Current Sample Images Thumbnail Strip

The thumbnail strip displays the Procedure date, time, and day for each individual image in orange text along the top.



Figure 108: Current Sample Images Thumbnail Strip

The current image displayed in the image acquisition window will be highlighted in green. *See above figure.*

Images in the thumbnail strip can be “dragged and dropped” into the image display or review panes, as necessary.

The user can select (click on) multiple images and then Save the Images to a USB drive (or other location), send the images to PACS, or delete the images. The selected images are highlighted in orange.

There are arrows located to the left and right of the thumbnail strip which allows the user to scroll through images when they exceed the space provided in the thumbnail strip.

There is an **Open Database** button, which is used to enter the Review mode. The **Open Database** button is the equivalent of the **Review Mode** button in the Home Menu. *See figure below and figure 110.*



Figure 109: Open Database Button

5.8 Database/Review Mode

Review Mode allows the user to open images from the database. To enter Review Mode, click on the **Review** button in the Home Menu, or the **Open Database** button in the Current Image Browser. *See figures above and below.* Either button will open the Database, *figure 111.* The system will be in Standby while in review mode.




Figure 110: Review Mode Button in the Home Menu

The Database is comprised of a Top Menu Bar and three main sections or sub-windows:

- The Sample List/Table, on the left,
- The Procedure List/Table, on the top right, and
- The Image Thumbnail List/Viewer located on the bottom right.




Figure 111: Image Database – Main Sections

The sections, or lists/tables, function in a top-down hierarchy:

Sample List→Procedure List→Image Thumbnail List/Viewer

**Note**

The Procedure List and Thumbnail List/Viewer will be empty until a sample is selected from the Sample List.

When a sample is selected, the Procedure window becomes active and displays a listing of all the procedures performed under the selected Sample. Additionally, all images for the selected Sample will be displayed in the Thumbnail List/Viewer. Selecting a procedure filters the images to only those taken under that procedure.

- Multiple procedure can be selected using the keyboard CTRL or SHIFT keys with the left mouse button. SHIFT + Click is used to select a group of adjacent procedures. CTRL + Click is used to select non-adjacent procedures. All the images belonging to all the selected procedures will be displayed.
- Multiple images can be selected by simply clicking on the specific images. The CTRL or SHIFT keys have no affect when selecting images.

The database may be sorted by clicking on the top fields of the Sample List/Table or Procedure List/Table. Clicking on a field such as Patient ID, Last Name, or Accession Number will sort the List/Table in ascending or descending order. The Field title used to sort the List/Table will contain an arrow to indicate whether the sorting is ascending or descending. The arrow will point up for ascending order and down for descending.

Scroll bars are located in the database to allow for ease of navigation when the List/Table content exceeds to available space.

**Note**

“Thumbnail” images are separate files from the actual image data (DCM, DAT, DMM, and RAW) files. They have the same file name as the image file but have a PNG or BMP file extension. If the image does not have a thumbnail VisionLS will attempt to create a thumbnail to be loaded into the database. If the software fails to create the thumbnail the image will not be displayed in the Thumbnail List/Viewer.

5.8.1 THE REVIEW PANE

The VisionLS allows up to 4 review panes to be opened in the review screen. These panes are only open in Review Mode.

Configuring the Review Panes

The review panes can be configured using the **Change Review Layout** button in the middle section of the VisionLS Bottom Tool Bar. The review layout can be configured with up to 4 panes. This button is only available during Review Mode.




Figure 112: Change Review Layout Button and Menu



Note

The Review Pane Layout can also be configured using the “Review Layout” menu in the “Top Menu Bar – View Menu”. See Section 5.4.4.

Opening Images into Review Panes

- Open the Image Review screen.
- Click on the desired pane. The border of that pane will turn green indicating it is active.




Figure 113: Review Panes

- Click on the **View Current Sample Images** (Browser) button in the middle of the Bottom Tool Bar.
- Select an image from the Thumbnail Strip, or click on the **Open Database** button, select a sample, select a procedure, then double click on an image. That image will be loaded into the active review pane and you will be returned to the Image Review screen.
- Repeat the above process for each Review Pane and image to be loaded.
- The **Change Review Layout** button indicates how many images have been loaded.

5.8.2 Database Top Menu Bar

The Top Menu bar consists of 4 buttons and a TAB along the top of the database screen.



Figure 114: Database Top Menu Bar



Note

The left button will adjust from **CLOSE DATABASE** to **BACK TO MAIN MENU** depending on the workflow used to open the database. If the database was opened from the Acquisition Mode, the button will be **CLOSE DATABASE**. If the database was opened using the **Review** button in the Home Menu, the button will be **RETURN TO MAIN MENU**.

From left to right they are:

- **BACK TO MAIN MENU** or **CLOSE DATABASE**, depending on how Review Mode was opened.
- **SAVE IMAGE(S)**
- **SEND IMAGE(S) TO PACS**
- **DELETE IMAGE(S)**, and
- **ADDITIONAL OPTIONS** (*Opens a drop-down menu.*)

Prior to selecting an image, only the **BACK TO MAIN MENU/CLOSE DATABASE** button, and the **ADDITIONAL OPTIONS** tab will be active. The **SAVE IMAGE(S)**, **SEND IMAGE(S) TO PACS**, and **DELETE IMAGE(S)** buttons will only be active after an image is selected from the *Image Thumbnail List/Viewer*.

SAVE IMAGE(S)

The **SAVE IMAGE(S)** button is available when one or more entries from the Image List/Table or Thumbnail List is selected, otherwise it is disabled.

Click on the image(s) you want to save, then click the **SAVE IMAGE(S)** button. The "Save File As" dialog will open.

- Click the **Browse** button to navigate to the folder you want to save the images to.

- Enter a “Custom Prefix” if desired, or
- Select an “Automated Prefix” from a drop-down list. The options are: “Sample ID”, “Operator Last Name”, “Operator First Name”, or “Accession Number”.
- An example of the file name will be displayed towards the bottom of the dialog.
- Click the **Save** button to save the image(s). Note that the process may take a minute or more.
- The image is saved as a DCM file. A PNG thumbnail will also be created.




Figure 115: Save Images Dialog



Note

Entering or selecting a prefix is optional. Both Custom Prefix and Automated Prefixes entries can be left blank.



Note

The window level of the saved image is the default window level.

SEND IMAGE(S) TO PACS:

The **SEND IMAGE(S) TO PACS** button is used to send images to a pre-configured DICOM Server and will open a DICOM Networking Dialog Box. *Refer to Section 5.9, DICOM Functionality for additional information.*

DELETE IMAGE(S)

Allows the user to delete one or more images from a procedure. Select the image(s) to be deleted in the Image List/Viewer then click the **DELETE IMAGE(S)** button. The software will request confirmation that you want to delete the image(s). Click the **OK** button to continue deleting the image(s) from the database. *See figure on next page.*




Figure 116: Delete Image(s) Confirmation

ADDITIONAL OPTIONS

The **ADDITIONAL OPTIONS** tab opens a menu of additional database options or functions. These options are described in the following pages.




Figure 117: Additional Options Menu

Open Menu

The "Open" Menu has only one option "Open Image from File" which opens a standard Windows dialog box that allows the user to navigate to and open image files located on local or network drives.




Figure 118: Additional Options Menu - Open

Print Menu

The "Print" menu provides options to print an image to a local printer, or to a DICOM printer (if one has been configured). Neither print option is enabled until an image is selected from the Image List/Table or Thumbnail List. *See figure on next page.*



Figure 119: Additional Options Menu - Print

- “Print Image to Local Printer” sends the selected image to the local default Windows printer.

The image will be printed with the default window level which may not provide an ideal printed image. If you require the window level to be adjust for the printed image, open the image in Review Mode, adjust the window level, then print using the “Print Current Image” in the *Top Menu, Tools Menu*.

- Print Image(s) to DCM Printer

This function sends the image(s) to a pre-configured DICOM printer.

View Menu

The “View” Menu provides the user options affecting the Database view.




Figure 120: Database View Menu

- Show All DB Entries

Reloads the Sample List/Table from the Database file. Any new changes or additions to the Database will be added to all the lists/tables. Any active Search Filter which may have been applied to the tables will be cleared.

Search Filter is described under the Search menu item below.

- Refresh DB Entries

Refreshes the Sample List/Table but maintains any existing Search Filter that has been applied.

- Browse Data Directory (No DB)

Displays all thumbnail images (.png files) in the Data directory established by the configuration file. Default directory is C://Faxitron/Data. The files are displayed in the Thumbnail List/Viewer. They are not associated with the Database, or table lists. This tool is helpful to review un-recorded images, or to debug database problems with images in the current data directory.

- Browse Any Directory (No DB)

Opens a typical Windows “Browse for Folder” dialog, which allows the user to browse to any file directory. The files are displayed in the Thumbnail List/Viewer. They are not associated with the Database, or table lists. This tool is helpful to review an image directory which was never recorded in the database, or when database records are corrupted.

Search Menu

The database “Search” tool allows the user to perform a detailed and specific database search, or filtering.




Figure 121: Database Search Menu

- Start Search

“Start Search” opens the DB Search Dialog which allows the user to filter the database by the selected criteria. A Drop-down menu provides a list of *search categories*, such as Accession Number, First Name, Last Name, or Procedure Name. The user types in the search parameter for that search category, then clicks the **OK** button to start the search. The Sample List/Table will be filtered to display only those samples matching the search parameter. Subsequent searches can be conducted to further narrow down the samples, procedures, or images.




Figure 122: Database Start Search Dialog Box



Note

When the DB Search Dialog is opened, the Database page is temporarily closed, and the software opens the Image Review page. The Database page will re-open when the DB Search Dialog is closed.

- Stop Search

“Stop Search” clears all active Search Filter and re-loads the list of all Samples to Sample List Table.

Tools Menu

The “DB-Tools” menu provides per Image/Procedure/Sample functions as well as some useful functions to repair or reformat the database.




Figure 123: Database Tools Menu

- **Image Properties**

This sub-menu is only active if an image has been selected from the Image List. Clicking on Image Properties will open the “Image Properties” window, which is the same “Image Properties” window that is opened from the Top Menu Bar – Tools – Properties. *See Section 5.6.4 and Figure 94.*

- **Edit Image Record**

“Edit Image Record” is currently not available in VisionLS software.

- **Add Record to DB**

“Add Record to DB” is only active when an entry from the Image List/Table or Thumbnail List is selected. This tool is useful for adding individual records/images to the database and can be used in conjunction with the “Browse Any Directory” function from the View Menu.

Click on an image in the Image List/Viewer, then click on “Add Record to DB”. That image will now be added to the database.

Note



“Add Record to DB” only allows you to add 1 image/record at a time.

- **Import Records in DIR to DB**

“Import Records in DIR to DB” allows the user to import all the image files in a directory (including sub-directories) to the database. The image formats can be DCM, DAT, DMM and RAW. The software will create thumbnails (BMP files) if they are missing from any of the images.

Selecting “Import Records in DIR to DB” will open a “Browse for folder” which allows the user to click through (navigate) to the desired folder/directory. When the directory has been selected click the **OK** button. The software will prompt for confirmation before importing the images. *See figure on next page.*



Figure 124: Import Records Confirmation Dialog

- Reformat all DB
“Reformat All DB” tool helps to reformat Date [YYYY/MM/DD], Date/Time [YYYY/MM/DD HH:MM:SS] fields of tables and will also capitalizes letters of entries for some other fields. It is a useful tool when an older database file is imported and opened for the first time, or when the user alters database entries manually using MS Access.
- Close Database
 Closes the database and opens the Image Review Screen. The image area will be empty if an image had not been selected before closing the database.

5.8.3 Edit or Delete Sample or Procedure Information from the Database

The user may edit or delete SAMPLE or PROCEDURE information from within the database. Right click on the SAMPLE or PROCEDURE entry to open an on-screen menu.



Note

If the system has an “Authentication Password” activated (see Section 5.6.4, Top Menu Bar > Tools > Advanced Tools > Login Utility, and Figure 83), you will be prompted for the same password before sample or procedure records can be edited or deleted.

Edit Sample Information

1. Right Click on the sample in the Sample List to open the on-screen menu shown below.




Figure 125: Edit or Delete Sample Information Menu

2. Select “Edit Sample Information” which will open a “Sample Information Editor” dialog box. See figure on next page.


A dialog box titled "Sample Information" with a green header and a close button (X) in the top right corner. It contains four input fields: "Operator Name:" with a red asterisk, a text field containing "GEORGE", a dropdown menu with "C" selected, and a text field containing "SCOTT"; "Sample ID:" with a red asterisk, a text field containing "2021"; "Date of Birth:" with a text field containing "yyyy/mm/dd" in red; and "Sex:" with a dropdown menu. At the bottom, there is a blue "Save Changes" button with a right arrow and a grey "Cancel" button with a left arrow.

Figure 126: Sample Information Editor Dialog Box

3. Make the desired changes. You may edit the Operator Name, Date of Birth, and Sex. The Sample ID cannot be edited.
4. Click **Cancel** to close the dialog without making any changes.
5. Clicking the **Save Changes** button will open a "Confirmation" dialog box.


A confirmation dialog box with a light grey background and a large yellow question mark icon on the left. The title is "Confirmation". The text reads: "You asked to ADD/EDIT the record of the selected SAMPLE. This action may yield undesired results or harm your DB, if SAMPLE\PROCEDURE\IMAGE records already exist. Would you like to PROCEED?". At the bottom are three buttons: "Yes", "No", and "Cancel".

Figure 127: Edit Sample Information Confirmation Dialog Box

6. Click **Yes** to save the changes. Click **No** or **Cancel** to discard the changes.

Delete Sample Information

1. Right click on the sample in the Sample List to open the on-screen menu shown in figure 125.
2. Select "Delete Sample Information", which will open a "Confirmation" dialog box.

A confirmation dialog box with a light grey background and a large yellow question mark icon on the left. The title is "Confirmation". The text reads: "Do you want to DELETE ALL the records for the selected SAMPLE? Click Yes to CONFIRM.". At the bottom are three buttons: "Yes", "No", and "Cancel".

Figure 128: Delete Sample Information Dialog Box

3. Click **Yes** to delete the records. Click **No** or **Cancel** to retain the records and close the dialog.

Edit Procedure Information

1. Right click on the procedure in the Procedure List to open the on-screen menu shown below.




Figure 129: Edit or Delete Procedure Information Menu

2. Select "Edit Procedure Information" which will open a "Sample/Procedure Information Editor" dialog box.
3. Make the desired changes.
4. Click **Cancel** to close the dialog without making any changes.
5. Clicking the **Save Changes** button will open a "Confirmation" dialog box. See figure on next page.

Sample Information

Operator Name: * GEORGE S SCOTT

Sample ID: * 2021

Date of Birth: yyyy/mm/dd

Sex: Female

Procedure Information

Accession No.: * Y

Procedure Name:

Laterality:

Position View:

Body Part:

Refer. Physician:

Technician:

Technician ID:

Study Information

Study Name:

Study ID:

Save Changes

Cancel

Figure 130: Sample/Procedure Information Editor



Note

None of the Sample Information fields can be edited.



Note

All Procedure and Study Information, including Accession No., can be edited.



Figure 131: Edit Procedure Information Confirmation Dialog Box

6. Click **Yes** to save the changes. Click **No** or **Cancel** to discard the changes.

Delete Procedure Information

1. Right click on the Procedures List to open the on-screen menu shown in figure 128.
2. Select "Delete Procedure Information", which will open a "Confirmation" dialog box.



Figure 132: Sample Information Editor Dialog box

3. Click **Yes** to delete the records. Click **No** or **Cancel** to retain the records and close the dialog.

5.9 DICOM Functionality

5.9.1 Modality

The DICOM file settings for Modality are accessed from the Top Menu Bar by selecting "DICOM Utilities" then "DICOM Modality Setup".




Figure 133: Top Menu – DICOM Utilities

Selecting "DICOM MODALITY Setup" opens the DICOM File/Transfer Settings dialog. Here the user can select the DICOM Storage Type (Presentation or Processing), DICOM Image Modality (Digital X-ray (General) or Digital X-ray Mammography Image) and may select Default DICOM Servers. The Value Representation and Transfer Syntax are not user definable. They are always assigned to be VR_IMPLICIT and META HEADER PRESENT, LITTLE_ENDIAN.

The dialog box is titled "DICOM File/Transfer Settings" and contains the following sections:

- DICOM Storage Type:**
 - ☒ Storage For Presentation
 - ☐ Storage For Processing
- Value Representation:**
 - ☒ Implicit
 - ☐ Explicit*
- DICOM Image Modality:**
 - ☐ Digital X-ray Image (General)
 - ☒ Digital X-ray Mammography Image
- Transfer Syntax:**
 - ☒ Little Endian
 - ☐ Big Endian*

(*) IMPLICIT-VR and LITTLE-ENDIAN syntax is recommended. Other options may not be fully compatible with third party software or DICOM services. IMPLICIT-VR and BIG-ENDIAN cannot be simultaneously selected.

Default DICOM Servers:

- STORE 1: [Dropdown menu]
- STORE 2: [Dropdown menu]
- ☒ Enable Single-Click Auto Store
- PRINT: [Dropdown menu]
- ☒ Enable Single-Click Auto Print
- MWL: [Dropdown menu showing "Worklist 1"]
- ☒ Enable Auto Search When Connecting to MWL Server
- ☒ Enable Single-Click Auto MWL Connection
- ☐ Autoclose dialog after [5] seconds of successful transfer/printing.

[DICOM Servers...] Click DICOM Servers button to add, configure or delete DICOM servers.

[OK] [Cancel]

Figure 134: DICOM Modality Settings Dialog Box – Default Settings Shown



Note

Once the DICOM Modalities have been selected, they remain in effect even after the VisionLS software has been closed. The Modality settings are stored and will remain in the configured manner until they are change again using the DICOM File Settings. These settings (Storage Type & Modality, Transfer Syntax, and Value Representation) will be effective for all DICOM Local Save or Print functionalities.

5.9.2 Default DICOM Servers

Refer to Figure 133.

This section contains settings for default DICOM servers and to automate some of the processes when storing, printing, or searching PACS and MWL servers.

“Enable Single-Click Auto Store”, “Print” or “MWL Connection” checkboxes allow the user to store or print images or automatically connect to the MWL server. If you use a single PACS server for those actions or majority of your work involves a single PACS server, check a desired box to automate those actions.

“Enable Auto Search When Connecting to MWL Server” will tell the application to automatically perform a search with default search settings when looking up sample’s data.

“Auto close dialog after N seconds of successful transfer/printing” will automatically close the DICOM Networking dialog after the specified number of seconds passed. The dialog will not close automatically if there was an error with transfer or printing.

Clicking the **DICOM Servers** button will open DICOM Networking dialog, where a user can add, modify, or delete DICOM servers.




Figure 135: DICOM Networking Dialog Box

5.9.3 DICOM Save Local

To save the current image to a local drive/directory, right click anywhere in the image. This will open the menu shown below. Click “Save/Store Current Image” then “As Local DICOM Image”. This will bring up the standard Windows File Save Dialog. Browse or click through to the desired directory/folder, enter the file name and click **Save**.




Figure 136: Image Right Click Menu – Save/Store Current Image

Each “Save/Store Current Image” sub-menu item is discussed in detail in Section 5.6.4, Top Menu Bar, under Tools Menu > Save/Store Current Image.

5.10 User Annotations

5.10.1 Overview

User Annotations in the Faxitron VisionLS Software allows the user to apply a graphic “overlay” commonly called “annotations” to an image. These annotations take the form of text, lines, “free form” lines, arrows, ellipses, and boxes. This allows additional information and “markers” to be added to the image if desired. Once applied to an image the annotations can be easily turned on or off, i.e., revealed or hidden, at will.

5.10.2 Operations

To begin annotating, select the **Annotations Tool** button from the left tool bar. *See figure 38 and figure below.* This opens the Annotations Dialog box and puts the software into annotations design mode.




Figure 137: Annotation Tool Button and Dialog Box




Figure 138: Annotation Tool Pallet

The tool pallet, along the top of the dialog box, consist of 8 annotation tools, described on the next page.



Note

The method for drawing each of the annotation objects is basically the same: Position cursor at starting point, left click, drag, and release at the end point. All of the annotation objects, except the Freehand Tool, follow a linear path, starting point to end point. The Freehand Tool follows irregular movements of the mouser/cursor and allows the user to trace, or encompass, irregular features of the sample.

Selection Tool	The Selection Tool is used to select (highlight) one or more current annotation objects. • To select a single annotation object, click on the tool, move the cursor over the annotation object and click. The annotation object will be highlighted. • To select multiple objects, click on the selection tool, click and hold the left mouse button down near one of the objects, then drag the mouse to encompass the annotation objects you want to select. Release the mouse button. Any annotation object within the rectangle will be highlighted.
Line Tool	The Line Tool is used to draw a straight line on the image.
Rectangle Tool	The Rectangle Tool is used to draw a rectangle or square on the image.
Ellipse Tool	The Ellipse Tool is used to draw an ellipse or circle on the image.
Arrow Tool	The Arrow Tool is used to draw a straight line with arrowhead at one end.
Freehand Tool	The Freehand Tool is used to draw irregular shapes. The tool will follow the mouse movements.
Text Tool	The Text Tool is used to draw a text box in the image window, which allows the user to type text-based data, make notes on the image.
Ruler Tool	Creates a ruled line segment. This is similar to the Image Tools Ruler.

When all the desired annotations have been made click the **Save** button to *save* the annotations to the image. An additional image file with the extension “.ANN” will be created. The software will exit annotations mode and return to image acquisition mode.

Click the **Cancel** button to exit annotation mode without saving any annotations.

5.10.3 Design Mode and Mouse Right-Click Menu

Annotation Tool Properties

The user can modify the properties of each annotation tool, except for the Selection Tool.

- Select the desired annotation tool and draw the object on the image.
- Right click on the annotation object. A menu will open.




Figure 139: Annotation Tool Right Click Menu and Properties Menu

- Select Properties (i.e., “Rectangle Properties”, “Line Properties”, etc.). This will always be the bottom menu item.
- The properties menu will vary slightly depending on the annotation tool. But generally, the user will be able to change the line thickness, line type, fill color and pattern (if applicable), foreground and background colors, along with a few additional options. See figure above.
- The “Fixed” menu item allows the user to adjust the behavior of annotations when an image is zoomed or scrolled. Annotations that are not fixed move with the image. i.e., if the zoom level is changed, the annotations will be adjusted accordingly. When annotations are fixed, they are tied to the window. i.e., if the zoom level is changed, or the image is panned, the annotation does not adjust accordingly. It’s size and position are locked to the monitor screen. *See figure on next page.*




Figure 140: Annotation Tool Properties – Fixed Dialog Box

Preset Annotation Tools

The user may modify an annotation tools' properties and save those changes as a new preset. The new preset will retain all the tool properties applied at the time it was created. Additionally, the new preset will retain the exact size and location existing when it was created. The new preset may be applied to any image in the future.

- Select and draw the desired annotation tool.
- Modify the properties as described above.
- Assign a name to the modified annotation tool (i.e., "Red Rule", "Wide Line") and enter that name in the "New Preset Name Here" space provided in the Annotations Tools Dialog box. See figure on next page.




Figure 141: Adding Preset Annotation Tools

- Click on the **Add as New Preset** button.
- The new “line” will be added to the Preset Dropdown Menu.




Figure 142: Preset Annotation Tools Drop-Down Menu



Note

Creating a New Preset temporarily applies the new presets’ properties to all the annotation tools in the pallet. After creating a new preset, close the Annotation Tools Dialog box. Re-enter Annotation Mode by clicking the **Annotation Tools** button in the side tool bar. All Annotation Tool Pallet items will behave with the default properties.

Using Preset Annotation Tools

- To utilize preset annotation tools, click the down arrow in the “New Preset Name Here” box to display a list of the available presets.
- Select the desired preset, then click on **Apply**.

- The tool will be applied to the image. Note that the tool will be located in the exact position and orientation as when it was created.
- To move or re-orientate the tool click on the **Selection Tool** button, then click on the tool/object. Grab the tool/object handles and adjust as desired.
- Click the **Save** button to save the new annotation(s) to the image and return to the image acquisition mode.
- Click the **Cancel** button to delete any new image annotations and return to the image acquisition mode.

Universal Annotations Right-Click Menu

Right Clicking anywhere in the image screen will pull up a “Universal” Annotations menu.




Figure 143: Annotation Tool Right Click Menu

Most of the menu items are inactive (greyed out) and not available for use.

A description of the active menu items follows:

- **Select All:** Clicking on Select All will draw a box around all annotation tools on the image. The user can delete all the selected tools using the keyboard “Delete” key. The user may also grab one of the selection box handles and resize all the selected tools.
- **Lock All...:** Lock All is used to prevent any active annotation tools from being deleted or altered. Clicking on Lock All brings up a password prompt. Your Hologic representative will provide you with the password.
- **Unlock All...:** Unlock All is used to reverse the Lock All action. As with Lock All, Unlock All will bring up a password prompt.
- **Hide Rotate Handles for All Objects.** Hides the Rotate Handles on all objects/tools. If an object is selected the Rotate Handles will not be shown. When active this menu item changes to **Show Rotate Handles for All Objects**.
- **Default Properties:** This menu item allows the user to set the “Default Properties” for all Annotation Tool Pallet items.




Figure 144: Annotation Tools - Right Click Menu – Default Properties Menu

5.10.4 Show – Hide Annotations

The **Toggle Annotations** button is located in the Right Tool Bar (See Section 5.6.2 and Figure 56) and allows users to quickly toggle annotations on and off. By default, the annotations will be toggled on.



Figure 145: Toggle Annotations Button

5.11 Help Button

The **Help** button is located in the Home Menu and is used to open the Help Menu.




Figure 146: Help Button – Home Menu




Figure 147: Help Menu/Page

From the Help Menu you can:

- Click the **Back** button to return to the Home Menu,
- **Open** the **Quick Start Guide**,
- Get the contact information for help in using your system, and
- Get information on contacting Hologic's technical support team through the **Remote Connect** button.

Chapter 6 Quality Control

6.1 Required Quality Control Procedures

The following procedures are necessary for correct system operation.

Table 1: Required Procedures

Test	Frequency
Calibration	Daily

6.2 Service and Troubleshooting

6.2.1 Service

There are no serviceable parts in the Faxitron Pro+. Refer to the copyright page at the front of this manual.

6.2.2 Troubleshooting

The following are a few issues that you may be able to solve with some basic troubleshooting steps:

Computer or System will not turn on.

- No power to the computer or x-ray cabinet.
- Verify power cords are plugged in and power to the system is ON.

Image looks like there is another object in the x-ray.

- Verify there are no other objects in the chamber during calibration.
- Remove objects and repeat calibration.

Poor image quality.

- Verify that the image intensity values (I (x,y)) found in the upper right-hand corner of the display screen are near 10,000 ADUs when preshot AEC target is also set to 10,000 ADUs, and exposure mode is set to Auto kV/Auto Time.
- Repeat the imaging procedure if needed.

6.3 Cybersecurity

If a cybersecurity event occurs disconnect the system from the network and contact Hologic immediately. We will work with your IT department to resolve the issue.

A cybersecurity even can be detected in a number of ways:

- Slower than normal operation and/or network connection.
- Suspicious pop-ups or home page changed in internet browser.

- Password(s) no longer working.
- Unidentified programs in start menu or system tray.
- Missing, corrupt, or altered data.

Chapter 7 Maintenance, Cleaning, and Disinfecting

7.1 General Information

7.1.1 For General Cleaning

This procedure includes cleaning the unit and the surfaces that the specimens do not contact.

Use a lint-free cloth and apply a diluted dishwashing liquid to clean components and surfaces.



Caution:

Use the least possible amount of cleaning fluids. The fluids must not flow or run.

7.1.2 For Disinfecting

To disinfect the system's components and surfaces, Hologic recommends the following disinfecting solutions:

- 10% chlorine bleach solution and water with one part commercially available chlorine bleach solution (normally 5.25% chlorine and 94.75% water) and nine parts water. Mix this solution daily for best results.
- Commercially available isopropyl alcohol solution (70% isopropyl alcohol by volume, not diluted).
- 3% maximum concentration of hydrogen peroxide solution.

After you apply any of the above solutions, use a clean cloth to clean the system, observing appropriate contact (wet) times.



Caution:

To prevent damage to the electronic components, do not use disinfectant sprays on the system.

7.1.3 To Prevent Possible Injury or Equipment Damage

Do not use a corrosive solvent, abrasive detergent, or polish. Select a cleaning/disinfecting agent that does not damage the plastics, aluminum, or carbon fiber.

Do not use strong detergents, abrasive cleaners, high alcohol concentration, or methanol at any concentration.

Do not expose equipment parts to steam or high temperature sterilization.

Do not let liquids enter the internal parts of the equipment. Do not apply cleaning sprays or liquids to the equipment. Always use a clean cloth and apply the spray or liquid to the cloth. If liquid enters the system, disconnect the electrical supply and examine the system before returning it to use.



Caution:

Wrong cleaning methods can damage the equipment, decrease imaging performance, or increase the risk of electric shock.

Always follow instructions from the manufacturer of the product you use for cleaning. The instructions include the directions and precautions for the application and contact time, storage, wash requirements, protective clothing, shelf life, and disposal. Follow the instructions and use the product in the most safe and effective method.

7.2 Component Specific Care

7.2.1 How to Clean the Image Display Monitor

- Do not touch the display when working with specimens.
- Use care when you clean the outer surface of the display area.
- Always use a clean, soft, lint-free cloth to clean the display area. Microfiber cloths, available at most stores, are recommended.
- Strong chemicals and abrasives can damage the display and should never be used.
- Never use a spray or flow a liquid onto the display.
- Never apply pressure to the display area.
- Never use a detergent with fluorides, ammonia, alcohol, or abrasives.
- Never use bleach.
- Never use steel wool.
- Never use a sponge with abrasives.

There are many commercially available products to clean LCD displays. Any product that is free of ingredients or abrasive described above and used according to the directions of the manufacturer can be safely used.

7.2.2 Power Cord

How to Inspect the Power Cords

1. Once every quarter, visually inspect the power cords for cuts, damage to the cover, or strain relief.
2. If any of the power cords appears damaged, contact Hologic for replacement of the power cord.

7.3 Maintenance

7.3.1 Preventive Maintenance Schedule

Table 2: User Preventive Maintenance

Action	Each Use	Daily	As Required
Clean and disinfect the specimen tray	✓		
Clean and disinfect the contact surface	✓		
Perform System Calibration		✓	
Clean and disinfect the Imaging Cabinet chamber			✓
Clean and disinfect the cabinet door handle			✓
Clean and disinfect the control display surface			✓

Table 3: Service Preventive Maintenance

Action	Done Annually
X-Ray Tube Conditioning ^{Note 1}	✓
Operational/Safety Checks	✓
Image Quality Check	✓
Radiographic Survey (X-Ray Leakage Test) ^{Note 2}	✓

Note 1: Performed annually, or if system has been unused for 2 months or more.

Note 2: Performed annually or as required by local regulatory agencies.

Appendix A System Specifications

A.1 Product Measurements




Figure 148: Path+ System Measurements

Path+ Cabinet/System Measurements

A	Cabinet Height	153 cm (60.3 in)
A ₁	System Height, floor to top of monitor	168 cm (66 in)
B	Cabinet Depth	80 cm (31.4 in)

B ₁	System Depth, Monitor Bracket/Keyboard Shelf positioned against cabinet	82 cm (32 in)
B ₂	System Depth, Monitor Bracket/Keyboard Shelf fully extended.	103 cm (40 in)
C	Cabinet Width	86 cm (34 in)
C ₁	System Width, Monitor Bracket/Keyboard Shelf positioned against cabinet	137 cm (54 in)
C ₂	System Width, Monitor Bracket/Keyboard Shelf fully extended.	160 cm (63 in)
D	Height, floor to Keyboard Shelf	120 cm (47.4 in)
	Weight	635 kg (1400 lbs)

A.2 Operation and Storage Environment

A.2.1 General Conditions for Operation

This equipment is designed for safe and effective operation under the following conditions for Indoor Use:

<i>Altitude</i>	<i>To 3,048 m (10,000 feet)</i>
<i>Degree of Pollution</i>	<i>Pollution Degree 2</i>
<i>Temperature Range</i>	<i>15.5 °C (60 °F) to 30 °C (86 °F)</i>
<i>Relative Humidity Range</i>	<i><75%</i>
<i>Over-voltage Category</i>	<i>II</i>

A.2.2 General Conditions for Transport and Storage

<i>Temperature Range</i>	<i>0°C (32°F) to 45°C (113°F)</i>
<i>Relative Humidity Range</i>	<i><75%</i>

A.3 Technical Specifications

A.3.1 System Specs

<i>Spatial Resolution</i>	<i>100μm (10μm with geometric magnification, up to 50+ lp/mm))</i>
<i>Source to Imager Distance (SID)</i>	<i>36" (90mm)</i>
<i>Exposure Control</i>	<i>Fully Automatic, Manual Time-Auto kV, Auto-Time-Manual kV or Fully Manual control</i>
<i>Geometric Magnification</i>	<i>1.5:1, 2:1, 3:1, 4:1, 5:1, 6:1, 8:1 and 10:1</i>
<i>Imaging Area</i>	<i>43cm x 43cm</i>
<i>Energy Range</i>	<i>20kV to 100kV. Configurable.</i>
<i>Exposure Time Range</i>	<i>1 s to 60.0 s. Configurable.</i>
<i>Typical Image Acquisition time</i>	<i><1 min</i>
<i>Target Material</i>	<i>Tungsten (W)</i>
<i>Beryllium Window Thickness</i>	<i>254μm (0.010")</i>
<i>Tube Current</i>	<i>0.5mA max (Limited to 12W)</i>
<i>Focal Spot Size</i>	<i>< 15 microns</i>
<i>Input Power</i>	<i>100-240VAC $\pm$10%, 50/60Hz, 300VA Max</i>
<i>DICOM Interface</i>	<i>Store, Print, Annotate and Modality Worklist capabilities</i>

A.3.2 Beam Coverage and Source to Object Distance Chart

<i>Shelf Position/Mag. Factor</i>	<i>BC nominal diameter cm (inches)</i>	<i>SOD cm (inches)</i>
10	6.1 (2.4")	9.2 (3.6")
8	7.62 (3.0")	11.5 (4.5")
6	10.2 (4.0")	15.3 (6.0")
5	12.2 (4.8")	18.3 (7.3")
4	15.3 (6.0")	22.9 (9.0")
3	20.3 (8.0")	30.5 (12.0")
2	30.5 (12.0")	15.8 (18.0")
1.5	41.6 (16.0")	61.0 (24.0")
1 (Floor)	61.0 (24.0")	91.5 (36.0")

A.3.3 X-Ray Tube

Hologic P/N:	130-552027
Manufacturer P/N:	18032 (MXR)
<i>Max Anode Voltage</i>	100 kV
<i>Max Anode Current*</i>	0.3 mA
<i>Maximum Power</i>	12W
<i>Focal Spot Size</i>	<15 μ m
<i>Anode Material / (Target)</i>	Tungsten (W)
<i>Cooling</i>	Air Cooled, 100 CFM fan recommended
<i>Beryllium Window Thickness</i>	0.25mm (0.010")
<i>Window Material</i>	Be 98.5%
<i>Cone Angle</i>	45 Degrees
<i>Operating/Storage Temp Range</i>	10°C-30°C
<i>Operating/Storage Humidity Range</i>	10%-70% non-condensing
<i>Weight</i>	3 lbs



***Note**

Max Anode current is limited by the UltraFocus® controller. Tube can operate at 0.5mA.

A.3.4 Computer Specifications

This table provide minimum recommended specifications at time of publication.

Hologic P/N: 2009-0029

<i>Processor</i>	<i>Intel® Core™ i5-8500 (6 Cores)</i>	
<i>Operating System</i>	<i>Windows 10 (64 bit)</i>	
<i>System RAM</i>	<i>8 GB (2s4GB) 2666MHz DDR4</i>	
<i>Hard Drive</i>	<i>M.2 512GB NVMe Class 35 Solid State Drive</i>	
	<i>(NOTE: Early system may have a 2.5-inch 500GB 5400rpm SATA Solid State Hybrid Drive)</i>	
<i>Optical Drive</i>	<i>N/A</i>	
<i>Network Connectivity</i>	<i>Integrated Intel I219-LM Ethernet LAN 10/100/1000</i>	
<i>USB Ports</i>	<i>FRONT Panel</i>	<i>REAR Panel</i>
	<i>USB 2.0 – 2X</i>	<i>USB 2.0 – 2X</i>
	<i>USB 3.1 Gen 1 – 1X</i>	<i>USB 3.1 Gen 1 – 4X</i>
	<i>USB 3.1 Gen 1 Type-C – 1X</i>	
<i>Video RAM</i>	<i>2256 MB (on board)</i>	
<i>Graphics</i>	<i>Integrated Intel® HD Graphics 630</i>	
<i>Case/Chassis</i>	<i>Small Form Factor (SFF)</i>	
<i>Dimensions (H x W x D)</i>	<i>11.42" x 3.7" x 11.50" (29.0cm x 9.26cm x 29.2cm)</i>	
<i>Min. Weight</i>	<i>11.57 lbs (5.26 kg)</i>	
<i>Monitor</i>	<i>27"</i>	
<i>Hologic P/N: 2009-0026</i>	<i>Ultra HD (3840 x 2160) 4K Resolution</i>	
	<i>1300:1 Contrast Ratio (typical)</i>	
	<i>350 cd/m² Brightness</i>	
	<i>Display Port Input</i>	

Glossary of Terms

AEC

Automatic Exposure Control

Annotations

Markings on an image to indicate an area of interest.

DICOM

Digital Imaging and Communications in Medicine

PACS

Picture Archiving and Communications System. A computer and network system that transmits and archives digital medical images.

ROI

Region of Interest

SID

Source to Image Distance

HOLOGIC®



Hologic Inc.
600 Technology Drive
Newark, DE 19702 USA
1.800.447.1856
www.hologic.com



Refer to the corporate website for more facilities worldwide.
www.hologic.com