LDRCheck-D

(bioMérieux ref. 39318, Hologic Cat. No. 301042F/1042F)

FOR EXPORT ONLY

LDRCheck-D
(bioMérieux ref. 39318, Hologic Cat. No. 301042F/1042F)
For Quality Control of Hologic AccuDR, Leader 50, Leader 100 and Leader 1 Lumimeter

INTENDED USE
The Hologic LDRCheck-D reagent is a chemiluminescent DNA control used to qualitatively assess performance of Hologic AccuDR, Leader 50, Leader 100, and Leader 1 Lumimeters.

SUMMARY AND EXPLANATION OF THE TEST
The LDRCheck-D control is used to monitor performance of the lumimeter every two weeks by measuring the control in the Hologic Leader Lumimeter. On the first day of use, the test result is calculated as a ratio of the Observed Value to LDRCheck-D reagent to the Expected label value of LDRCheck-D reagent. The test result is compared to the quality control limits for the given Leader lumimeter model. The Observed Day 0 value is recorded and used for trending purposes. The instrument variation from test to test should fall within a ratio of 0.90 to 1.10 relative to the Observed Day 0 value.

PRINCIPLES OF THE PROCEDURE
The concentration of LDRCheck-D reagent has been determined and appears on the label as RLU / 100 µL. LDRCheck-D reagent is measured in the lumimeter by beads and the test results are compared against the quality control limits.

REAGENT INFORMATION
For information on any hazard and precautionary statements that may be associated with reagents, refer to the Safety Data Sheet library at www.hologic.com/dfs. Hologic LDRCheck-D Quality Control Reagent for Leader Lumimeters.
< 100 ng/dish non-infectious DNA with a chemiluminescent label in succinate buffered solution containing detergent.

WARNINGS AND PRECAUTIONS

A. For in vitro diagnostic use only.
B. Avoid contact of LDRCheck-D reagent with skin, eyes and mucous membranes. Wash with water if these reagents come into contact with skin or eyes. If spills of these reagents occur, dilute with water and wipe. Avoid contact of LDRCheck-D reagent with skin, eyes and mucous membranes. Wash with water if contact with these reagents occurs. If spills of these reagents occur, dilute with water before wiping dry.
C. Use only supplied or specified disposable laboratory ware.
D. Avoid contact of Detection Reagent Kit (Cat. No. 201791/1791, bioMérieux ref. 39300) with LDRCheck-D reagent.

STORAGE AND HANDLING REQUIREMENTS
A. LDRCheck-D reagent should be stored where not in use to avoid changes in concentration due to evaporation.
B. LDRCheck-D reagent is stable when stored at 2° to 25° until the expiration date indicated. Storage above room temperature may result in product deterioration.
C. The LDRCheck-D test is to be performed at room temperature.
D. Avoid contaminating the reagent by contact with used or contaminated pipette tips.

MATERIALS

A. MATERIALS PROVIDED
Hologic LDRCheck-D Quality Control Reagent for AccuDR (formerly PAL), Leader 50, Leader 100 and Leader 1 Lumimeters. (bioMérieux ref. 39318, Hologic Cat. No. 301042F/1042F)

B. MATERIALS REQUIRED BUT NOT PROVIDED
Hologic AccuDR (formerly PAL) lumimeter
Hologic Detection Reagent Kit (bioMérieux ref. 39300, Hologic Cat. No. 201791/1791)
Hologic Reaction Tubes (bioMérieux ref. 39307, Hologic Cat. No. 102065/2065)
Calibrated micropipets capable of delivering 100 µL or Eppendorf Repeating Pipettor with Corbittips (0.5 mL or 1.25 mL)

C. MATERIALS AVAILABLE FROM Hologic
Hologic Leader 50 Lumimeter (bioMérieux ref. 39400, Hologic Cat. No. 1031003/1000)
Hologic Leader 100 Lumimeter (bioMérieux ref. 39300, Hologic Cat. No. 201791/1791)
Hologic Reaction Tubes (bioMérieux ref. 39307, Hologic Cat. No. 102065/2065)

TEST PROCEDURE
REAGENT PREPARATION
1. If stored at ≤ 15°C, allow reagent to come to room temperature before equilibrating at room temperature for one hour. Alternatively, swirl the bottle for 30 seconds in a water bath at 60° ± 1°C and cool to room temperature and allow to sit for 10 minutes. Once equilibrated, LDRCheck-D reagent may be stored at room temperature.
2. Before pipetting sample, mix the LDRCheck-D reagent by swirling and gently inverting three times to mix. Do not vortex.

EQUIPMENT PREPARATION
1. Prepare the Hologic lumimeter for operation. Make sure there is sufficient volume of Detection Reagent(s) to run the instrument. Add reagent to the instrument.
2. Print out the instrument parameters as described in the Operator's Manual. On AccuDR read and record instrument parameters.
3. Prepare the lumimeter to run raw data using a second count time as described in the instrument's operator's manual.

CHEMICAL BACKGROUND CHECK
Perform a chemical background check on the Detection Reagents.
1. Run the instrument wash protocol as described in the Operator's Manual.
2. Load ten empty tubes for a one second count time in data mode. The exterior of each tube is to be wiped with a water dampened, lint free wipe before reading in the lumimeter.
3. If the average reading is < 150 RLU for Leader 1, Leader 50 or Leader 100, use a clean, lint free wipe for AccuDR, then repeat the Detection Reagents with a fresh set of detection reagents and repeat the test. If the result is not, notify your Hologic distributor/Customer Service.

D. MEASURING LDRCheck-D REAGENT
LDRCheck-D reagent should be read on the lumimeter a minimum of every two weeks.
1. Precision pipette ten replicates of LDRCheck-D reagent (100 µL per tube) into clean polypropylene tubes.
2. Read the LDRCheck-D samples plus 10 empty tubes using one second count time on raw data.
4. When the analysis is complete, remove the tubes from the lumimeter.

PROCEDURAL NOTES
A. Use calibrated micropipets to pipette LDRCheck-D reagent.
B. Pipette all samples to the bottom of the tube. Do not splash on the side of the tube in order to obtain accurate results.
C. The exterior of each tube is to be wiped with a water dampened, lint free wipe before being read in the lumimeter.
D. The situation of LDRCheck-D reagent should appear as a clear liquid with no clumps. If microbial contamination is apparent, discard bottle. Obtain a new bottle of LDRCheck-D reagent and proceed.

A. In any reagent system, excess powder on some golves may cause contamination of opened reagent or reaction tubes. Hologic recommends that customers experiencing difficulty with the test avoid using this type of laboratory glove. Using powdered gloves (no talcum powder) will avoid this difficulty.

RESULTS

A. COEFFICIENT OF VARIATION (NCV)
The Coefficient of Variation (NCV) for LDRCheck-D reagent should be ≤ 6.0%. If the NCV is greater than 6.0%, then prepare and read a new set of LDRCheck-D reagents. If the NCV is greater than 6.0%, then contact your Hologic distributor/Customer Service.
* For AccuDR calculate the NCV as follows:
NCV = Standard Deviation of LDRCheck-D in RLU / PLU × 100%

B. CALCULATION OF RESULTS
1. DAY 0 - INITIAL OBSERVED/EXPECTED RATIO
Calculate the Observed Value / Expected Value using the equation below:
Observed Value = (Avg. RLU) / LDRCheck-D × Avg. RLU / empty tubes
Expected Value = (Avg. RLU) / LDRCheck-D × Avg. RLU / empty tubes

2. SUBSEQUENT TEST DATES - OBSERVED / DAY 0 RATIO
Calculate the Observed Value / Day 0 Value ratio using the equation below:
Observed Value = (Avg. RLU) / LDRCheck-D × Avg. RLU / empty tubes
Day 0 Value = (Avg. RLU) / LDRCheck-D × Avg. RLU / empty tubes (observed on Day 0)
* For AccuDR to convert PLU to RLU: RLU = PLU × 0.031

QUALITY CONTROL RESULTS AND ACCEPTABILITY
A. DAY 0 - INITIAL OBSERVED/EXPECTED RATIO
The Observed/Expected (OE) ratio must be within the range indicated below. If the OE ratio is within this range, the instrument is functioning properly. If the OE range is still not met, perform a second quality control check. If the OE range is still not met, contact your Hologic distributor/Customer Service.

B. SUBSEQUENT TEST DATES - OBSERVED / DAY 0 RATIO
The Observed / Day 0 ratio must be between 0.90 and 1.10. If the OE ratio is within this range, the instrument is functioning properly. If the OE range is not met, review all instructions and calculations and perform a second quality control check. If the OE range is still not met, contact your Hologic distributor/Customer Service.

C. This procedure should be followed whenever a new bottle of LDRCheck-D reagent is to be opened and used for checking instrument calibration.

LIMITATIONS

A. LDRCheck-D reagent contains EDTA. The reagent is bacteriostatic, however the presence of EDTA may affect the RLU value of the reagent. If contamination is observed or suspected, the reagent should not be used and should be replaced by a previously unopened bottle of LDRCheck-D reagent.
B. Failure to follow the instructions given in this insert may result in erroneous results.

LDRCheck-D
(bioMérieux ref. 39318, Hologic Cat. No. 301042F/1042F)
Zur Qualitätskontrolle der Hologic Lumimeter AccuDR, Leader 50, Leader 100 und Leader 1

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