HOLOGIC®

AccuProbe®

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(bioMérieux ref. 39305 /
Hologic Cat. No. 102800)

CULTURE IDENTIFICATION REAGENT KIT

For Export Use Only

(bioMérieux ref. 39305 / Hologic Cat. No. 102800)

INTENDED USE

The reagents in the ACCUPROBE CULTURE IDENTIFICATION REAGENT KIT are to be used with the ACCUPROBE CULTURE IDENTIFICATION REAGENT TESTS according to the instructions outlined in the respective Package Inserts. Reagents are provided for 200 bacterial or 100 mycobacterial or fungal determinations per kit.

REAGENTS AND MATERIALS PROVIDED

Note: For information on any hazard and precautionary statements that may be associated with reagents, refer to the Safety Data Sheet Library at www.hologic.com/sds.

CULTURE IDENTIFICATION KIT
(bioMérieux ref. 39305 / Hologic Cat. No. 102800)

200 Bacterial or 100 Mycobacterial/Fungal Tests

Reagent 1 (Lysis Reagent) (1) Buffered solution containing 0.04% sodium azide	1 x 10 mL
Reagent 2 (Hybridization Buffer) (2) Buffered solution	1 x 10 mL
Reagent 3 (Selection Reagent) (3) Buffered solution	1 x 60 mL

WARNINGS AND PRECAUTIONS

For *in vitro* diagnostic use.

Reagents in this kit contain sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. Upon disposal of the reagents, always dilute the material with a large volume of water to prevent azide build-up in the plumbing.

STORAGE AND HANDLING REQUIREMENTS

The reagents contained in the ACCUPROBE CULTURE IDENTIFICATION REAGENT KIT are to be stored at 2° to 25°C and are stable until the expiration date indicated.

DO NOT FREEZE THESE REAGENTS.

Note: Reagent 2 (Hybridization Buffer) must be at room temperature before use.

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©1990 - 2018 Hologic, Inc. All rights reserved.
Hologic, Inc.