

ThinPrep® Stain Rinse Solution

Instructions for Use



INTENDED USE

ThinPrep® Rinse Solution is intended to be used in a Papanicolaou staining procedure in conjunction with ThinPrep Nuclear Stain, ThinPrep Bluing Solution, ThinPrep Orange G Solution, and ThinPrep EA Solution.

SUMMARY AND EXPLANATION

ThinPrep Rinse Solution is an aqueous, detergent solution that removes excess ThinPrep Nuclear Stain from cellular components on slide preparations for cytologic evaluation.

PRINCIPLES OF PROCEDURE

When used in conjunction with the ThinPrep staining solutions and cytologic slide preparations, ThinPrep Rinse Solution will remove excess ThinPrep Nuclear Stain from cellular components, enhancing nuclear staining and clarifying cytoplasm.

COMPOSITION

An aqueous solution of propriety surfactant.

WARNINGS

For In Vitro Diagnostic Use.

Not for external or internal use in humans or animals. May cause irritation on contact.

PRECAUTIONS

Avoid contact with eyes.

Wear gloves and eye protection when handling solution.

Keep container closed when not in direct use.

Results with non-ThinPrep cytology samples have not been evaluated.

Do not use if primary packaging is damaged.

PRETREATMENT

No reconstitution, mixing or dilution is required.

STORAGE

Store unused ThinPrep Rinse Solution at 15°C (59°F) to 30°C (86°F) in the container provided. Do not use ThinPrep Rinse Solution beyond the expiration date marked on the container.

APPEARANCE AND INTEGRITY

Clear, non-sterile solution.



ThinPrep® Stain Rinse Solution

SPECIMEN COLLECTION AND PREPARATION

For use with gynecologic cytology samples that have been processed on any ThinPrep Processor.

PROCEDURE

See the ThinPrep Stain User's Manual for specific staining protocols.

LIMITATIONS OF PROCEDURE

Must be used according to the instructions provided in the ThinPrep Stain User's Manual. ThinPrep Rinse Solution cannot be substituted with any other solution.

PERFORMANCE CHARACTERISTICS

When used as directed on cytologic slide preparations, ThinPrep Rinse Solution removes excess ThinPrep Nuclear Stain from cellular components.

DISPOSAL

Dispose in accordance with all applicable regulations.

FIRST AID MEASURES

No ill effects are anticipated. See www.hologic.com/sds for the entire Safety Data Sheet.



Manufacturer



Authorised representative in the European Community



Batch code



Consult instructions for use



Catalogue number



Temperature limitation



In vitro diagnostic medical device



Use by



Hologic, Inc. • 250 Campus Drive • Marlborough, MA 01752 USA
Tel: +1 508 263-2900 • www.hologic.com

 Hologic BVBA • Da Vinciaan 5 • 1930 Zaventem • Belgium

Australian Sponsor: Hologic (Australia) Pty Ltd • Suite 302, Level 3 • 2 Lyon Park Road
Macquarie Park NSW 2113 • Australia • Tel: 02 9888 8000