

Artwork and Signature File for:

MAN-04965 “AFFIRM PRONE BIOPSY ANTIVIRUS SOFTWARE INSTALLATION”

Artwork consists of:

- Six (6) 8 ½ inch x 11 inch sheet(s) attached.

REV AUTHORED BY	DATE	HOLOGIC®	SIGNATURES ON FILE	
J. SHARP	09/29/16			
REV DRAFTED BY	DATE	TITLE	DOCUMENT NUMBER	REV
J. SHARP	09/29/16		AW, AFFIRM PRONE BIOPSY ANTIVIRUS SOFTWARE INSTALLATION	AW-15412
PROPRIETARY INFORMATION: The content of this document is the exclusive property of Hologic and may not, without prior written permission of Hologic, be reproduced, copied or used for any purpose whatsoever.		ARTWORK	SIZE A	SHEET 1 OF 1

1.1 Introduction

1.1.1 Purpose

To install antivirus software on the Affirm™ Prone Biopsy System product.

1.1.2 Scope

This document applies to all Affirm Prone Biopsy System products with version 1.0.1 software and above.

1.1.3 Estimated Time

Installation of antivirus products takes approximately 30 minutes to complete including configuration.

1.1.4 Reference List

Hologic has verified only the products and versions that are in the following list. These products and versions do not interfere with the operation of the system. Hologic does not recommend installing any other product or version. Hologic cannot guarantee the effectiveness of these products in the prevention of malicious software.

This document provides instructions for the following products.

- Symantec™ Endpoint Protection 12.x
- McAfee® Enterprise VirusScan 8.8.x
- Sophos® Endpoint Security & Control 10
- Trend Micro® OfficeScan® 11



NOTE

The customer must provide these products. Load only the client program and only one antivirus program per system. Refer to the appropriate section of the installation guide.

1.1.5 Definitions

- **LiveUpdate** – A feature that allows servers and clients to retrieve updates from an internal server or Symantec’s official LiveUpdate server.
- **Managed** – The client system is configured to send virus alerts and retrieve virus updates from an internal parent Symantec server.
- **Real-Time** – Real-time scanning of each file that is loaded in RAM.
- **SmartScan** – A scanning technique that scans the header of each file to determine its true file extension and to identify possible malicious code.
- **Unmanaged** – The clients do not connect to the network nor do they have a parent server with which they communicate. These clients must download their own virus definition updates.

1.2 Customer Preparation Checklist

Before the installation, note the following:

- Hologic does not supply antivirus software. It is the responsibility of the customer to procure the software and associated licenses.

1.3 Preinstallation Checklist

Prior to the installation, review the following:

- Ensure that no existing antivirus software is loaded on the workstation prior to installation.
- Ensure that the installer has the proper serial keys and associated licenses for the product to install.

1.4 Install the Software

Follow the instructions for the software installation:

- [Installing Symantec Endpoint Protection v12.x as Unmanaged Client](#) on page 3
- [Installing Symantec Endpoint Protection v12.x in a Managed Setting](#) on page 4
- [Installing McAfee Enterprise VirusScan v8.8](#) on page 5
- [Installing Sophos Endpoint Security and Control 10](#) on page 6
- [Installing Trend Micro OfficeScan 11](#) on page 6

1.4.1 Installing Symantec Endpoint Protection v12.x as Unmanaged Client

1. On the Acquisition Workstation, press and hold down the **Ctrl** key while selecting the **Shutdown** button. Continue holding down the **Ctrl** key while selecting the **Exit** button to close and shut down the Acquisition Workstation application program properly.
2. Install the software:
 - a. Insert the media from Symantec.
 - b. Browse to the media and double-click **setup.exe** to start the installation process.
 - c. When the Symantec window appears, click **Install Symantec Endpoint Protection Client**.
 - d. A Welcome window appears. Click **Next**.
 - e. Click **I accept the terms in the license agreement**, then click **Next**.
 - f. Select **Unmanaged client** radio button, then click **Next**.
 - g. Select **Custom**, then click **Next**.
 - h. Click the drop-down box in front of **Application and Device Control** and select **x** to disable this component from installing.
 - i. Click the drop-down box in front of **Network Threat Protection** and select **x** to disable this component from installing.
 - j. Click **Next**.
 - k. Ensure **Enable Auto-Protect** and **Run LiveUpdate** boxes are checked, then select **Next**.
 - l. Click **Install**.
 - m. Click **Finish**, then **Exit**.
 - n. Launch Symantec Endpoint Protection from **Start > All Programs**.
 - o. When prompted to start client service for Symantec Endpoint Protection, click **No**.
 - p. You are now looking at the Symantec Endpoint Protection dashboard. Ensure that there is no warning about out-of-date definitions. If there is a warning, work with your local IT group to configure the system before clicking **Fix**. (**Fix** requires an active internet connection to download the definition files from Symantec.)
3. Configure protection (balancing security with performance).
 - a. Select the **Change settings** tab on the left side pane of the Symantec Endpoint Protection dashboard.
 - b. Click **Configure Settings** next to the Antivirus and Antispyware Protection
 - c. Click the **File System Auto-Protect** tab, then click the **Selected** check box under "File Types" section. Click **OK** to close.

- d. Select the **Change Settings** tab on the left side pane of the Symantec Endpoint Protection dashboard.
 - e. Click **Configure Settings** to the right of Centralized Exceptions.
 - f. Click **Add > Security Risk Exceptions > Folder**.
4. Exclude the following directories and its subdirectories from scanning:
- C:\Aries
 - C:\Images



Note

If the drive letters do not apply to your configuration, substitute the appropriate drive letters.

1.4.2 Installing Symantec Endpoint Protection v12.x in a Managed Setting

Symantec Endpoint Protection in a managed configuration is an IT-centric product that is geared toward an enterprise with IT support staff. Hologic assumes that the customer:

- has the infrastructure already running,
- has the personnel with expertise to deploy and manage the antivirus product, and
- knows which files or directories on Hologic systems to exclude from scanning.

It is recommended that only the antivirus client is installed, and features such as firewall or application control are left uninstalled. Other features such as firewall or application control are not recommended, as they increase the risk of reduced productivity due to the high maintenance required. Install them at your own risk if you have the local personnel with the right skill set to configure and maintain them.

For optimal system performance, ensure that the installed antivirus software excludes the following directories and their subdirectories from scanning:

C:\Aries

C:\Images



Note

If the drive letters do not apply to your configuration, substitute the appropriate drive letters.

1.4.3 Installing McAfee Enterprise VirusScan v8.8

1. On the Acquisition Workstation, press and hold down the **Ctrl** key while selecting the **Shutdown** button. Continue holding down the **Ctrl** key while selecting the **Exit** button to close and shut down the Acquisition Workstation application program properly.
2. Install the software:
 - a. Insert the McAfee media.
 - b. Browse to the media and double-click **setupvse.exe** to start the installation process. The McAfee console appears.
 - c. Click **Next** at the McAfee VirusScan Enterprise Setup screen.
 - d. Select the appropriate licensing information, click **I accept**, and then click **OK**.
 - e. Select **Typical**, then click **Next**.
 - f. Choose **Protection Level Standard**, then click **Next**.
 - g. Click **Install**.
 - h. When installation is complete, unselect **Run On-Demand Scan**.
 - i. Click **Finish**.
 - j. If you are prompted to reboot, reboot the system.
3. Once the Acquisition Workstation completes booting, press and hold down the **Ctrl** key while selecting the **Shutdown** button. Continue holding down the **Ctrl** key while selecting the **Exit** button to close and shut down the Dimensions application program properly.
4. Exclude the following directories from scanning:
 - a. Double-click the **McAfee** shield in the system tray and select **Properties**.
 - b. In the left pane, click **All Processes**, then select the **Exclusions** tab.
 - c. Add the following directories and its subdirectories to Exclusions:
 - C:\Aries
 - C:\Images



Note

Do not forget to check the **Also exclude subfolders** box when you add directory to the exclusion.



Note

If the drive letters do not apply to your configuration, substitute the appropriate drive letters.

1.4.4 Installing Sophos Endpoint Security and Control 10

Sophos is an IT-centric product that is geared toward an enterprise with IT support staff. Hologic assumes that the customer:

- has the infrastructure already running,
- has the personnel with expertise to deploy and manage the antivirus product, and
- knows which files or directories on Hologic systems to exclude from scanning.

It is recommended that only the antivirus client is installed, and features such as firewall or application control are left uninstalled. Other features such as firewall or application control are not recommended, as they increase the risk of reduced productivity due to the high maintenance. Install them at your own risk if you have the local personnel with the right skill set to configure and maintain them.

For optimal system performance, ensure that the installed antivirus software excludes the following directories and their subdirectories from scanning:

C:\Aries

C:\Images



Note

If the drive letters do not apply to your configuration, substitute the appropriate drive letters.

1.4.5 Installing Trend Micro OfficeScan 11

Trend Micro OfficeScan antivirus is an IT-centric product that is geared toward an organization with dedicated IT staff. Hologic assumes that the customer:

- has the infrastructure already running,
- has the personnel with expertise to deploy and manage the antivirus product, and
- knows which files or directories on Hologic systems to exclude from scanning.

For optimal performance, be sure Trend Micro OfficeScan excludes the following directories and subdirectories from scanning:

C:\Aries

C:\Images



Note

If the drive letters do not apply to your configuration, substitute the appropriate drive letters.
