



# ThinPrep® Pap Test: specimen collection

## Training bulletin


The detection of cervical cancer and its precursors as well as other gynecologic abnormalities is the primary purpose of obtaining a cervical cell sample. The following guidelines are referenced from CLSI Document GP15-A3<sup>1</sup> and are recommended in the collection process for obtaining a ThinPrep® Pap Test (TPPT) specimen. In general, the guidelines state that it is important to obtain a specimen that is not obscured by blood, mucus, inflammatory exudate or lubricant.




Unsatisfactory specimen  
obscured by blood




Unsatisfactory specimen  
obscured by mucus



Unsatisfactory specimen  
obscured by inflammation



Unsatisfactory specimen  
obscured by lubricant



Satisfactory ThinPrep Pap Test specimen

### Patient information

- The patient should be tested two weeks after the first day of her last menstrual period, and definitely not when she is menstruating.  
Even though the TPPT reduces obscuring blood, clinical studies have demonstrated that excessive amounts of blood may still compromise the test and possibly lead to an unsatisfactory result.<sup>2</sup>
- The patient should not use vaginal medication, vaginal contraceptives, vaginal creams, vaginal jellies, or douches during the 48 hours before the exam.
- The patient should refrain from intercourse 48 hours prior to the exam<sup>3</sup>.

# ThinPrep® Pap Test: specimen collection

## Training bulletin

### Specimen collection preparation

- **Prepare the speculum.**

For patients without physical or physiological need for lubricant, use lukewarm water to warm and lubricate the speculum.

Water lubrication has the fewest risks to the quality of the Pap sample collected.<sup>4</sup>

When necessary, sparingly apply **carbomer-free** lubricant on the exterior of the speculum blades.

If lubricant is necessary due to patient discomfort or the use of a plastic speculum, sparingly apply a thin film of **carbomer-free** lubricant on the speculum's surface, avoiding the tip.

**Do not use an excessive amount of lubricant jelly to lubricate the speculum.**

Hologic® evaluated a variety of popular lubricants and found those containing carbomer or carbopol polymers (thickening agents) may interfere with obtaining a representative cervical sample or cause artifact in the alcohol-based transport medium.<sup>4</sup> Hologic® recognizes the varying availability of different types of lubricants and recommends that, if used, any lubricant should be applied sparingly.

- **Remove excess mucus or other discharge present before taking the sample. This should be gently removed with ring forceps holding a folded gauze pad.**

The excess cervical mucus is essentially devoid of meaningful cellular material and when present in the sample vial may yield a slide with little or no diagnostic material present.

- **Remove inflammatory exudate from the cervical canal before taking the sample. Remove by placing a dry 2-by-2 inch piece of gauze over the cervix and peeling it away after it absorbs the exudate or by using a dry procto swab or Scopette® swab.**

The excess inflammatory exudate is essentially devoid of diagnostic cellular material and, when present in the sample vial, may yield a slide with little or no diagnostic material present.

- **The cervix should not be cleaned by washing with saline or it may result in a relatively acellular specimen.**
- **The sample should be obtained before the application of acetic acid.**

### Specimen collection

Refer to Hologic's ThinPrep Pap Test Quick Reference Guides, Part No.s DS-05867-001 and DS-05720-001, and ThinPrep® Pap Test Specimen Collection Protocol Video<sup>3</sup>.

### Collection device rinsing

Refer to Hologic's ThinPrep Pap Test Quick Reference Guides, Part No.s DS-05867-001 and DS-05720-001, and ThinPrep® Pap Test Specimen Collection Protocol Video<sup>3</sup>

**hologic.com | diagnostic.solutions@hologic.com | +1.781.999.7300**

1. Papanicolaou Technique Approved Guidelines (CLSI Document GP15-A3) 2. Lee et al. Comparison of Conventional Papanicolaou Smears and a Fluid-Based, Thin-Layer System for Cervical Cancer Screening. *Ob Gyn* 1997; 90: 278-284. 3. Saslow D, et al. American Cancer Society guideline for the early detection of cervical neoplasia and cancer. *CA Cancer J Clin* 2002;52:342-62. 4. Lubricant Use during Pap Sample Collection, Part No. MISC-00579 Rev. 006 5. Specimen Collection Protocol Video, Part No. MOV-00052-001 Rev. 002.

MISC-02720-001 Rev. 001 ©2015 Hologic, Inc. All rights reserved. Hologic, Science of Sure, ThinPrep and associated logos are trademarks and/or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries. All other trademarks are the property of their respective owners. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your Hologic representative or write to [diagnostic.solutions@hologic.com](mailto:diagnostic.solutions@hologic.com).

