


Artwork consists of one (1) 8 ½ inch x 11 inch page.

REV AUTHORED BY P.KUPOVICH	DATE 2/13/08	HOLOGIC ® A FAMILY OF COMPANIES	Osteoporosis Assessment LORAD® Breast Cancer Detection DirectRay® Digital Imaging FLUOROSCAN® C-arm Imaging	 SIGNATURES ON FILE
REV DRAFTED BY G CORRIDORI	DATE 2/27/08			
PROPRIETARY: This document contains proprietary data of Hologic, Inc. No disclosure, reproduction or use of any part thereof may be made except by written permission from Hologic.		TITLE	DOCUMENT NUMBER	REV
		IHE INTEGRATION STATEMENT, APEX V2.3	AW-02360	001
REV. RELEASE DATE:	03/28/08	ARTWORK	SIZE A	SHEET 1 OF 1

IHE Integration Statement		Date
		13 February 2008
Vendor	Product Name	Version
Hologic, Inc.	Apex	2.3
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Acquisition Modality	Modality Image Storage
		Storage Commitment
		Patient Based Worklist Query
		Broad Based Worklist Query
Modality Performed Procedure Step		
Internet address for the vendor's IHE information:		
http://www.hologic.com/cc/iheinteg.htm		
Links to Standards Conformance Statements for the Implementation		
HL7	http://www.hologic.com/cc/dicom.htm	
DICOM	http://www.hologic.com/cc/dicom.htm	
Links to general information on IHE		
In North America: www.rsna.org/ihe	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j