


Artwork consists of one (1) 8 ½ inch x 11 inch pages.

REV AUTHORED BY S. KRZYSKO	DATE 03/04/09	<b>HOLOGIC</b> ® A FAMILY OF COMPANIES	Osteoporosis Assessment LORAD® Breast Cancer Detection DirectRay® Digital Imaging FLUOROSCAN® C-arm Imaging	 <b>SIGNATURES ON FILE</b>
REV DRAFTED BY C. PETERSON	DATE 03/16/09			
PROPRIETARY: This document contains proprietary data of Hologic, Inc. No disclosure, reproduction or use of any part thereof may be made except by written permission from Hologic.		TITLE <b>APEX V2.4.2 IHE INTEGRATION STATEMENT</b>	DOCUMENT NUMBER <b>AW-02815</b>	REV <b>002</b>
REV. RELEASE DATE:	03/13/09	ARTWORK	SIZE A	SHEET 1 OF 1

<b>IHE Integration Statement</b>		<b>Date</b>
		2 March 2009
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>
Hologic, Inc.	Apex	2.4.2
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
Scheduled Workflow	Acquisition Modality	Modality Image Storage
		Storage Commitment
		Patient Based Worklist Query
		Broad Based Worklist Query
Modality Performed Procedure Step		
<b>Internet address for the vendor's IHE information:</b>		
<a href="http://www.hologic.com/cc/iheinteg.htm">http://www.hologic.com/cc/iheinteg.htm</a>		
<b>Links to Standards Conformance Statements for the Implementation</b>		
<b>HL7</b>	<a href="http://www.hologic.com/cc/dicom.htm">http://www.hologic.com/cc/dicom.htm</a>	
<b>DICOM</b>	<a href="http://www.hologic.com/cc/dicom.htm">http://www.hologic.com/cc/dicom.htm</a>	
<b>Links to general information on IHE</b>		
In North America: <a href="http://www.rsna.org/ihe">www.rsna.org/ihe</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>