


Artwork consists of one (1) 8 ½ inch x 11 inch page.

REV AUTHORED BY PETER LEWIS	DATE 4/1/11	HOLOGIC ® Osteoporosis Assessment LORAD® Breast Cancer Detection DirectRay® Digital Imaging A FAMILY OF COMPANIES FLUOROSCAN® C-arm Imaging	 SIGNATURES ON FILE	
REV DRAFTED BY SBROWMACIOCI	DATE 4/1/11			
PROPRIETARY: This document contains proprietary data of Hologic, Inc. No disclosure, reproduction or use of any part thereof may be made except by written permission from Hologic.		TITLE IHE INTEGRATION STATEMENT, APEX V3.3	DOCUMENT NUMBER AW-06400	REV 001
REV. RELEASE DATE: 05/03/11	ARTWORK		SIZE A	SHEET 1 OF 1

IHE Integration Statement		Date
		1 April 2011
Vendor	Product Name	Version
Hologic, Inc.	Apex	3.3
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Acquisition Modality	Modality Image Storage
		Storage Commitment
		Patient Based Worklist Query
		Broad Based Worklist Query
Modality Performed Procedure Step		
Internet address for the vendor's IHE information:		
http://www.hologic.com/cc/iheinteg.htm		
Links to Standards Conformance Statements for the Implementation		
HL7	http://www.hologic.com/cc/dicom.htm	
DICOM	http://www.hologic.com/cc/dicom.htm	
Links to general information on IHE		
In North America: www.rsna.org/ihe	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j