

**IHE Integration Statement for Selenia v3.4.x**

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product</b>	<b>Version</b>	<b>Date</b>
Hologic	Selenia Full Field Digital Mammography Acquisition Workstation (AWS)	3.4.X	March 2009
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow	Acquisition Modality	Assisted Acquisition Protocol Setting	
		PPS Exception Management Option	
		Broad Worklist Query	
		Patient Based Worklist Query	
Patient Information Reconciliation	Acquisition Modality	None	
Mammography Image Profile	Acquisition Modality	None	
Access to Radiology Information	Image Display	None	
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	Not applicable		
<b>DICOM</b>	<a href="http://hologic.com/cc/dicom.htm">http://hologic.com/cc/dicom.htm</a>		
<b>Links for general information on IHE</b>			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	