

## ATEC® Introducer Localization System with the SureSight™ Obturator

The SureSight™ obturator simplifies the targeting of lesions under MRI guidance. In the targeting scans, an area of increased signal intensity mimics the location of the sampling aperture on the ATEC® biopsy needle. This provides the physician a visual confirmation of where the needle aperture will be during the biopsy.



### The SureSight Obturator Provides Outstanding Visual Confirmation



SureSight in the sagittal view looks like a bullseye.



SureSight in the axial view shows where the needle aperture will be during biopsy.

# ATEC<sup>®</sup> with SureSight<sup>™</sup> Obturator


## Increased Confidence, Faster Results

By listening to our customers, we've developed another advancement in MRI biopsy making the ATEC MRI biopsy procedure even faster and easier.

- Higher confidence that the aperture is in the right location prior to biopsy
- Increased speed of the biopsy procedure

## SureSight Obturator Compared to the ATEC MRI Needle

### Standard



### Petite

