


Artwork consists of one (1) 8 ½ inch x 11 inch page.

REV AUTHORED BY S KRZYSKO	DATE 9/28/07	HOLOGIC ® Osteoporosis Assessment A FAMILY OF COMPANIES LORAD® Breast Cancer Detection DirectRay® Digital Imaging FLUOROSCAN® C-arm Imaging	 SIGNATURES ON FILE	
REV DRAFTED BY G CORRIDORI	DATE 10/15/07			
PROPRIETARY: This document contains proprietary data of Hologic, Inc. No disclosure, reproduction or use of any part thereof may be made except by written permission from Hologic.		TITLE	DOCUMENT NUMBER	REV
		IHE INTEGRATION STATEMENT, INSIGHT	AW-01079	002
REV. RELEASE DATE:	11/20/07	ARTWORK	SIZE A	SHEET 1 OF 1



1 IHE Integration Statement		1.1 Date
		28 September 2007
1.2 Vendor	1.3 Product Name	1.4 Version
Hologic, Inc.	Fluoroscan InSight	2.1
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	1.5 Options Implemented
Scheduled Workflow	Acquisition Modality	Modality Image Storage
		Storage Commitment
		Modality Worklist Query
		Modality Performed Procedure Step
Internet address for the vendor's IHE information: http://www.hologic.com/cc/dicom.htm		
1.6 Links to Standards Conformance Statements for the Implementation		
DICOM	http://www.hologic.com/cc/dicom.htm	
Links to general information on IHE		
In North America: www.rsna.org/ihe	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j