



eviva[™]
Stereotactic Breast Biopsy

The First Fully Integrated Solution in stereotactic breast biopsy

Eviva™ breast biopsy device incorporates firing, pain management and site marking into one integrated system delivering superior clinical results and patient benefits.



Unparalleled Patient Benefits

- **Quiet pneumatic firing** mechanism removed from adapter and incorporated into the device for patient comfort
- **4.5 second** automated tissue acquisition with continuous vacuum means less time on the table
- **Saline lavage** and constant aspiration helps ensure a core with every cycle ensuring a positive patient experience
- **Patent-pending Y-Valve** allows application of your choice of analgesic to ensure comfort without interruption during the biopsy procedure

Advancing Biopsy Innovation

- **End-deploy** biopsy site identification system with SecurMark™ for Eviva ensures accurate delivery
- **Core visualization** of each sample from handpiece to remote tissue filter
- **Modified profile** for superior targeting and lesion access
- **Newly designed tocar tip** allows improved penetration of dense tissue
- **Petite** device allows for compression to 16mm



Integrated Biopsy Site Identification Solution

SecurMark™ for Eviva Biopsy Site Identifier

Our innovative system ensures easy and accurate deployment during a stereotactic biopsy.

- No pull back for deployment
- Rigid end-deploy delivery system
- Marker is comprised of two components
 - Bio-absorbable suture-like netting
 - Bio-compatible titanium
- Two shapes available



Eviva Introducer with SecurMark for Eviva



First & Second Shape



Eviva with Open Aperture

SSL-00022 (02/09) © Hologic 2009. All rights reserved.
Printed in USA. Specifications are subject to change without prior notice. Hologic, ATEC, Eviva, MultiCare, SecurMark, the Hologic and Eviva logos are trademarks and/or registered trademarks of Hologic, Inc. in the United States and other countries.

HOLOGIC™

Corporate Headquarters

35 Crosby Drive,
Bedford, MA 01730-1401 USA
Tel: +1.781.999.7300
Sales: +1.781.999.7453
Fax: +1.781.280.0668
sales@hologic.com
www.hologic.com

Europe

Everest (Cross Point)
Leuvensesteenweg 250A
1800 Vilvoorde, Belgium
Tel: +32.2.711.4680
Fax: +32.2.725.2087

Asia Pacific

3/F, 21 Li Yuen Street West
Central, Hong Kong
Tel: +852.3102.9200
Fax: +1.781.280.0668