

## Gardnerella vaginalis Promoter T7 Oligo

### Analyte Specific Reagent (ASR).

Analytical and performance characteristics are not established.

### Intended Use

Analyte Specific Reagent. Analytical and performance characteristics are not established. This analyte specific reagent is intended to be used as a component of a Lab Developed Test (LDT) used exclusively by laboratories. Laboratories can develop their own LDTs using one or more Analyte Specific Reagents and/or general laboratory reagents and/or general laboratory instruments for *in vitro* diagnostic purposes.

### Materials Provided

Gardnerella vaginalis Promoter T7 Oligo (Cat No. PRD-03862)

Component	Volume	Concentration	Description
Tube	1 x 0.28 mL	19.4 nmol/mL solution	<i>DNA oligonucleotide primer for Gardnerella vaginalis 16S rRNA, includes promoter sequence for T7 RNA polymerase</i>

### Warnings and Precautions

- For Laboratory Use only. Analytical and performance characteristics are not established.
- Use routine laboratory precautions. Wear disposable, powderless gloves, protective eye wear, and laboratory coats when handling specimens and kit reagents. Wash hands thoroughly after handling reagents.
- Dispose of all materials that have come in contact with reagents according to local, state, and federal regulations.

### Product Specifications

#### Quality Control

This analyte specific reagent was manufactured and released in accordance with the Hologic quality control and quality assurance procedures. This reagent has been analytically tested to verify oligonucleotide mass and concentration. However, it is the responsibility of each laboratory to develop and validate their own tests for clinical diagnostic purposes.

#### Storage and Stability

Please refer to product labels for information on storage and stability.



Hologic, Inc.  
10210 Genetic Center Drive  
San Diego, CA 92121 USA

Customer Support: +1 800 442 9892  
customersupport@hologic.com

Technical Support: +1 888 484 4747  
molecularsupport@hologic.com

For more contact information visit [www.hologic.com](http://www.hologic.com).

Hologic and associated logos are trademarks and/or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries.

All other trademarks that may appear in this package insert are the property of their respective owners.

This product may be covered by one or more U.S. patents identified at [www.hologic.com/patents](http://www.hologic.com/patents).

©2018 Hologic, Inc. All rights reserved.

AW-16913 Rev. 001  
2018-01