



Date: June 25, 2019

Author: Software Engineering

System: SELENIA DIMENSIONS / 3DIMENSIONS **Subsystem:** Microsoft
Windows OS

Subject: Microsoft Windows OS Patch release June 2019

Purpose:

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

List of patches/KB's validated for June 2019:

KB4503292 KB4503269 KB4503259

List of all patches/KB's validated on the SELENIA DIMENSIONS / 3DIMENSIONS System:

KB976932	KB4480960	KB4103712
KB4503292	KB4474419	KB4096418
KB4503269	KB4471328	KB4093108
KB4503259	KB4468323	KB4074837
KB4499175	KB4467106	KB4074587
KB4498206	KB4462915	KB4056897
KB4495606	KB4457145	KB4054521
KB4495588	KB4457044	KB4048960
KB4493448	KB4457035	KB4041678
KB4490628	KB4343899	KB4038779
KB4489885	KB4339284	KB4034679
KB4486564	KB4338823	KB4025337
KB4483187	KB4284867	KB4022722

KB4019990	KB3074543	KB2943357
KB4019263	KB3072305	KB2937610
KB4015546	KB3071756	KB2931356
KB4012212	KB3067903	KB2912390
KB3205394	KB3061518	KB2911501
KB3197867	KB3060716	KB2900986
KB3192391	KB3059317	KB2894854
KB3177467	KB3055642	KB2894844
KB3172605	KB3046269	KB2893294
KB3161949	KB3046017	KB2888049
KB3159398	KB3045685	KB2884256
KB3156016	KB3037574	KB2882822
KB3150220	KB3035126	KB2871997
KB3139914	KB3030377	KB2868038
KB3139398	KB3023215	KB2864202
KB3138910	KB3022777	KB2862335
KB3138612	KB3021674	KB2862330
KB3127220	KB3019978	KB2862152
KB3126587	KB3018238	KB2861698
KB3124275	KB3011780	KB2840631
KB3110329	KB3010788	KB2834140
KB3109560	KB2992611	KB2813430
KB3109103	KB2991963	KB2813347
KB3108664	KB2984972	KB2807986
KB3108371	KB2978742	KB2789645
KB3101722	KB2977292	KB2786081
KB3097989	KB2973201	KB2770660
KB3092601	KB2973112	KB2758857
KB3086255	KB2972211	KB2742599
KB3080079	KB2972100	KB2736422
KB3075220	KB2968294	KB2731771

KB2729452

KB2479943

KB2729094

KB2727528

KB2719033

KB2716513

KB2706045

KB2705219

KB2698365

KB2690533

KB2685939

KB2670838

KB2667402

KB2656356

KB2654428

KB2653956

KB2639308

KB2631813

KB2621440

KB2620704

KB2604115

KB2585542

KB2579686

KB2564958

KB2560656

KB2552343

KB2533623

KB2533552

KB2532531

KB2506212

KB2506014

KB2491683