

IHE Integration Statement		Date
		25 June 2013
Vendor	Product Name	Version
Hologic, Inc.	APEX	5.5
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Acquisition Modality	Modality Image Storage
		Storage Commitment
		Patient Based Worklist Query
		Broad Based Worklist Query
		Modality Performed Procedure Step
Internet address for the vendor's IHE information:		
http://www.hologic.com/en/product-support/bone-densitometry/discovery/		
Links to Standards Conformance Statements for the Implementation		
HL7	http://www.hologic.com/en/product-support/bone-densitometry/discovery/	
DICOM	http://www.hologic.com/en/product-support/bone-densitometry/discovery/	
Links to general information on IHE		
In North America: www.rsna.org/ihe.aspx	In Europe: www.ihe-europe.net	In Japan: http://www.ihe-j.org/en/