

**VZV Probe****Intended Use**

Analyte Specific Reagent (ASR). Analytical and performance characteristics are not established. This ASR is intended to be used as a component of a laboratory developed test (LDT) used exclusively by laboratories. Laboratories may develop their own LDTs using one or more ASRs and/or general laboratory reagents for *in vitro* diagnostic purposes.

**Product Description**

The Varicella zoster virus (VZV) Probe ASR contains an oligonucleotide designed to detect the ORF28 gene of VZV. The oligonucleotide contains a 5' FAM label and a 3' BHQ-1 label.

**Materials Provided**

VZV Probe (Cat. No. PRD-04995)

Component	Quantity	Concentration
Tubes	1 x 0.53 mL	25 µmol/L

**Warnings and Precautions**

- A. Analytical and performance characteristics are not established.
- B. Use routine laboratory precautions. Wear disposable, powderless gloves, protective eye wear, and laboratory coats when handling specimens and kit reagents. Wash hands thoroughly after handling reagents.
- C. Dispose of all materials that have come in contact with reagents according to local, state, and federal regulations.

**Quality Control**

This analyte specific reagent was manufactured and released in accordance with the Hologic quality control and quality assurance procedures.

**Storage and Handling**

Upon receipt, store material at -15 to -85 °C. Multiple freeze/thaw cycles should be minimized by appropriate aliquoting and should not exceed 3 freeze/thaw cycles.



Hologic, Inc.  
10210 Genetic Center Drive  
San Diego, CA 92121 USA

Customer Support: +1 800 442 9892  
customersupport@hologic.com

Technical Support: +1 888 484 4747  
molecularsupport@hologic.com

For more contact information visit [www.hologic.com](http://www.hologic.com).

Hologic and associated logos are trademarks or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries.

“Black Hole Quencher” and “BHQ” are trademarks of and licensed from Biosearch Technologies, Inc., Novato, CA. The Black Hole Quencher (BHQ) dye technology is the subject of existing or pending patents and is licensed and sold under agreement with Biosearch Technologies, Inc.

All other trademarks that may appear in this package insert are the property of their respective owners.

This product may be covered by one or more U.S. patents identified at [www.hologic.com/patents](http://www.hologic.com/patents).

© 2018 Hologic, Inc. All rights reserved.

AW-18452 Rev. 001  
2018-12