

IHE Integration Statement for Selenia® Dimensions®/3Dimensions™ Systems			
Vendor	Product	Version	Date
Hologic®, Inc.	Selenia Dimensions/3Dimensions Acquisition Workstation (AWS)	1.9/2.0	September 2017
This product implements all transactions required in the IHE Radiology Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Access to Radiology Information	Image Display	None	
Digital Breast Tomosynthesis	Acquisition Modality	For Presentation Breast Projection X-Ray Images	
		For Processing Breast Projection X-Ray Images	
		Media Creation	
Digital Breast Tomosynthesis	Print Composer	None	
Mammography Image	Acquisition Modality	None	
Mammography Image	Print Composer	None	
Patient Information Reconciliation	Acquisition Modality	None	
Portable Data for Imaging	Portable Media Creator	None	
Radiation Exposure Monitoring	Acquisition Modality	None	
Scheduled Workflow	Acquisition Modality	Assisted Acquisition Protocol Setting	
		PPS Exception Management Option	
		Broad Worklist Query	
		Patient Based Worklist Query	
Scheduled Workflow	Image Display	None	
Links to Standards Conformance Statements			
HL7	Not applicable		
DICOM	www.hologic.com/product-support/digital-mammography/dimensions/		
Links for general information on IHE			
In North America: www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	