

CTB-01183

Date: April 10, 2023

Author: Service Engineering

Product: Horizon

Subsystem: APEX

Subject: APEX 5.6.1.3 rev 007 Software Release Notes

APEX™ Software Release Notes

APEX 5.X is the operating system for the Horizon® series of densitometers. For differences in results from previous versions of APEX, refer to the Statement of Equivalency (MAN-09016).

**Note**

Refer to the Horizon User Guide (MAN-08072-002) for operating instructions related to new features.

Version 5.6.1.3 rev 007**Notes**

APEX version 5.6.1.3 is for Windows® 10 operating system only.

Resolved Issues

- Analysis Shutdown Fixed:
 - The “Display Region of Interest Coordinates” checkbox can now be checked again.
 - APEX software will no longer intermittently crash when auto-analysis is enabled, and/or when auto-compare is enabled.