

BKV Probe**Intended Use**

Analyte Specific Reagent (ASR). Analytical and performance characteristics are not established. This ASR is intended to be used as a component of a laboratory developed test (LDT) used exclusively by laboratories.

Product Description

The BK virus (BKV) Probe ASR contains an oligonucleotide designed to detect the VP2 gene of BKV. The oligonucleotide contains a 5' CAL Fluor Red 610 label and a 3' BHQ-2 label.

Materials Provided

BKV Probe (Cat. No. PRD-04992)

Component	Quantity	Concentration
Tubes	1 x 0.53 mL	25 µmol/L

Warnings and Precautions

- A. Analytical and performance characteristics are not established.
- B. Use routine laboratory precautions. Wear disposable, powderless gloves, protective eye wear, and laboratory coats when handling specimens and kit reagents. Wash hands thoroughly after handling reagents.
- C. Dispose of all materials that have come in contact with reagents according to local, state, and federal regulations.

Quality Control

This analyte specific reagent was manufactured and released in accordance with the Hologic quality control and quality assurance procedures.

Storage and Handling

Upon receipt, store material at -15 to -85 °C. Multiple freeze/thaw cycles should be minimized by appropriate aliquoting and should not exceed 3 freeze/thaw cycles.



Hologic, Inc.
10210 Genetic Center Drive
San Diego, CA 92121 USA

For country-specific Technical Support and Customer Service email address and telephone number, visit www.hologic.com/support.

For more contact information visit www.hologic.com.

Hologic and associated logos are trademarks or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries.

“Black Hole Quencher,” “BHQ,” and “CAL Fluor” are trademarks of and licensed from Biosearch Technologies, Inc., Novato, CA. The Black Hole Quencher (BHQ) and CAL Fluor dye technologies are the subject of existing or pending patents and are licensed and sold under agreement with Biosearch Technologies, Inc.

All other trademarks that may appear in this package insert are the property of their respective owners.

This product may be covered by one or more U.S. patents identified at www.hologic.com/patents.

© 2024 Hologic, Inc. All rights reserved.

AW-18449 Rev. 002
2024-04