

Revised Nameplate Label

This addendum modifies your Selenia® Dimensions®/3Dimensions™ product manuals for the system nameplate label on the rear of the Gantry. A new RFID circuit board used in the compression device was implemented, requiring a new FCC registration number. This FCC registration number is stated on the system nameplate label, requiring a new system nameplate label.

- For systems containing RFID circuit board PCB-00116, LBL-05704 is used as the system nameplate label (item 4 in Table 1 of your manuals).




Figure 1: System nameplate LBL-05704 for FCC registration of RFID circuit board PCB-00116

- For systems containing RFID circuit board PCB-01673, LBL-07018 is used as the system nameplate label (item 4 in Table 1 of your manuals).




Figure 2: System nameplate LBL-07018 for FCC registration of RFID circuit board PCB-01673

Place this addendum with your product manuals for future reference.