

Marker Compatibility

for Ultrasound Breast Biopsy

Biopsy Device	Marker Name	Shape	Metal	Bio-absorbable Material (PLA/PGA)	Deployment Method
ATEC® Device	SecurMark® for ATEC® Biopsy Site Marker				
	SMark-ATEC-13-09	"mini cork" 	Titanium	Yes	Needle Hub (short shaft)
	SMark-ATEC-2s-13-09	"top hat" 	Titanium	Yes	Needle Hub (short shaft)
	SMark-ATEC-13-12	"mini cork" 	Titanium	Yes	Needle Hub (short shaft)
	SMark-ATEC-2s-13-12	"top hat" 	Titanium	Yes	Needle Hub (short shaft)
	SMark-A-13-12-ss2	"infinity" 	Stainless Steel	Yes	Needle Hub (short shaft)
	TriMark® Biopsy Site Marker				
	TriMarkTD-13-09	"cork" 	Titanium	No	Needle Hub (short shaft)
	TriMarkTD-2s-13-09	"hourglass" 	Titanium	No	Needle Hub (short shaft)
	TriMarkTD-13-12	"cork" 	Titanium	No	Needle Hub (short shaft)
	TriMarkTD-2s-13-12	"hourglass" 	Titanium	No	Needle Hub (short shaft)
	ATEC®, Celero® and Sertera® Devices	SecurMark® Biopsy Site Marker			
SMark-Celero		"mini cork" 	Titanium	Yes	Introducer or Free Hand
SMark-2s-Celero		"top hat" 	Titanium	Yes	Introducer or Free Hand
SMark-U-ss1		"buckle" 	Stainless Steel	Yes	Introducer or Free Hand
SMark-U-ss2		"infinity" 	Stainless Steel	Yes	Introducer or Free Hand
SMark-U-ss3		"stoplight" 	Stainless Steel	Yes	Introducer or Free Hand
CeleroMark™ Biopsy Site Marker					
CeleroMark		"cork" 	Titanium	No	Introducer or Free Hand
CeleroMark-2s		"hourglass" 	Titanium	No	Introducer or Free Hand
Tumark® Biopsy Site Marker					
351230		"vision" 	Nitinol	No	Introducer or Free Hand
351220		"Q" 	Nitinol	No	Introducer or Free Hand
351210	"X" 	Nitinol	No	Introducer or Free Hand	

Note: Not all markers are available worldwide.

MISC-01696 Rev. 006 (3/17) © 2017 Hologic Inc. Printed in USA. All rights reserved. Specifications are subject to change without prior notice. Hologic, the Science of Sure, ATEC, Celero, CeleroMark, SecurMark, Sertera and TriMark are trademarks or registered trademarks of Hologic and/or Hologic subsidiaries in the United States and/or other countries. Tumark and associated logos are trademarks and/or registered trademarks of Somatex Medical Technologies GmbH. Hologic is an authorized distributor of Tumark Professional X, Q and Vision biopsy site markers. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshow, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative or write to womenshealth@hologic.com.