

# The first biopsy system to unite **tissue acquisition, real-time imaging, verification and post-biopsy handling.**

Brevera® Breast Biopsy System with CorLumina® Imaging Technology



# REAL POWER REAL TIME



**Are you looking for a way to improve biopsy productivity?  
Increase your confidence in tissue acquisition? Provide more  
compassionate patient care?**

The Hologic Brevera<sup>®</sup> breast biopsy system with CorLumina<sup>®</sup> imaging technology delivers. It's the world's first and only solution to combine **vacuum-assisted tissue acquisition, real-time imaging, verification** and **advanced post-biopsy handling – all in one, integrated system.**

There's no moving to another room for imaging. No waiting for sample verification. No interrupting your screening schedule. Instead, you have the detailed, real-time information you need to **deliver an improved patient experience, enhance the workflow,** and **increase accuracy.**

**The Hologic Brevera breast biopsy system with CorLumina imaging technology will change the way you think about performing breast biopsies.**



**Improve the patient experience.**

Fast, accurate, and streamlined procedures mean less time under compression and can result in a more positive and compassionate biopsy experience for you and your patients.



**Enhance workflow.**

An intuitive user interface, real-time imaging, and automated specimen collection and separation work together to save facilities an average of 13 minutes per procedure.<sup>4</sup>




**Increase accuracy.**

Real-time imaging delivers a wealth of information at the point of care – so you can make informed clinical decisions with confidence.

# Give your Patients a more Compassionate Biopsy Experience

Because the patient experience is in Hologic's DNA, we wanted to create an efficient biopsy system to help reduce the anxiety and stress of having a biopsy. That's why we designed the Brevera® breast biopsy system with CorLumina® imaging technology.



## More Comfortable

**9 out of every 10** patients felt less uncomfortable with their Brevera breast biopsy procedure than they expected.<sup>1</sup>



## Less Invasive

Hologic's premier stereotactic biopsy solution offers patients a **positive cosmetic outcome** post-biopsy with maximized core retrieval and minimal incision size.<sup>2</sup>



# 92%

## Faster

In **over 92%** of biopsy procedures performed with the Brevera system, patients felt the biopsy experience was faster than they expected.<sup>1</sup>

## Fewer steps for shorter procedures.

The Brevera system eases the biopsy journey – for you and your patients.

### Conventional Stereotactic Breast Biopsy



### Brevera System Integrated Breast Biopsy





## Image accuracy, in hand.

### Biopsy Patients

- **Spend less time** under compression to enhance the patient experience.
- **Ensure accurate transfer of samples** from radiology to pathology for peace of mind.

### Physicians and Technologists

- **Confidently verify** excised tissue using CorLumina real-time imaging technology.
- **Quickly identify** calcifications with built-in image enhancement tools on a high-resolution, touch-screen monitor.
- **Potential to shorten procedures an average of up to 25% or more<sup>5</sup>** to free up your time and can result in a more positive patient experience.
- **Automatically separate, identify, and analyze** target samples, then transfer seamlessly to pathology with minimal handling to help maintain core integrity.

### Breast Imaging Facilities

- **Streamline workflow** to free up physicians and technicians to perform other important tasks.
- **Increase availability** of mammography systems to take on additional screenings.
- **Integrate** with your upright or prone biopsy systems to maximize your investment.
- **Lower the cost** of handling waste<sup>6</sup> and simplify inventory management.

6  
Manually  
Image  
Samples

7  
Place  
Marker

8  
Prep  
Samples  
for Pathology

Saving an average of  
**13 minutes per procedure<sup>4</sup>**



# Transforming 2D and 3D™ breast biopsy – in real-time.<sup>7</sup>

The Brevera system streamlines the entire breast biopsy process from start to finish – with real-time imaging for instant verification and automated post-biopsy specimen handling, and the potential to cut procedure time by up to 25% or more<sup>5</sup> to transform the patient experience and boost productivity.

## The future is in your hands.

The Brevera breast biopsy system features a reusable device driver and disposable needles to provide new functionality, simplify storage, and improve waste management.<sup>6</sup> The remote system operation puts control where it belongs – in the hands of the physician.



Brevera Biopsy Device

Integrated pain management and saline lavage — so you can administer local anesthetic to reduce pain and cleanse the biopsy site during each biopsy cycle.



**Control system operation and monitor performance instantly.**

**Quickly verify samples on the easy-to-use, 15-inch touchscreen monitor using a comprehensive suite of imaging tools.**

**Annotate, print, and export procedure information for a higher level of accountability.**

**Access patient worklist to help streamline workflow.**

**Tilt, swivel, rotate, and adjust the height of the Brevera system's dual monitors to suit the way you work.**


**Collect and image tissue in real time with our proprietary CorLumina imaging technology.**



Take advantage of wireless connectivity for advanced image sharing and transfer of patient records – with a single touch.



Redesigned disposable needle for 54% less waste<sup>2</sup> and better inventory management.




Variable aperture settings to help access challenging lesions.



Integrated site marking simplifies operation.



Remote operation frees technologists to focus on patients while giving physicians full control of the procedure.



Accurately separate and identify samples with a 12-chamber tissue filter to improve core handling and sample integrity.



In-line specimen imaging provides real-time feedback for more information at the point of care.

# Better Together.



Make fast, accurate decisions at the point of care, streamline workflow, and deliver an exceptional patient experience

 **94%**

Over 94% of patients reported a **faster biopsy procedure** than they expected when performed on an Affirm Prone® biopsy system with the Brevera breast biopsy System.<sup>1</sup>

 **90%**

On average, in almost 90% of procedures, patients reported a **less uncomfortable biopsy experience** when using the Brevera system together with the Affirm Prone biopsy system.<sup>3</sup>



Paired together, the **Brevera® biopsy system** and the **Affirm™ Prone biopsy guidance system** accelerate the speed and efficiency of your prone biopsy procedures to further reduce your patients' time under compression<sup>1</sup>.

To learn more, talk with your Hologic representative. Or visit [BreveraBiopsy.com](http://BreveraBiopsy.com)

[www.hologic.com](http://www.hologic.com) | [BSHInsideSales@hologic.com](mailto:BSHInsideSales@hologic.com) | **+1.800.442.9892**

**REFERENCES:** 1. Results from, "A Prospective, Block Stratified Clinical Trial to Evaluate the Performance and Operation of the Brevera® Breast Biopsy System" (NCT03300206) (US) 2. MISC-04810 3. Results from, "A Prospective Study to Evaluate the Performance and Operation of the Brevera® Breast Biopsy System" (NCT03456583) (EU) 4. Brevera Pulse Survey, Quantitative report. Inspired Insights. April 2019. 5. 2015 Kadence International survey of 200 healthcare professionals. 6. Compared to Eviva® device. Hologic Data on file 2017. 7. When used in conjunction with the Affirm® breast biopsy system.

**Hologic Headquarters  
United States / Latin America**  
250 Campus Drive  
Marlborough, MA 01752  
USA Tel: +1.508.263.2900  
Sales: +1.781.999.7453  
Fax: +1.781.280.0668  
Email: [info@hologic.com](mailto:info@hologic.com)

**Hologic Europe**  
The Corporate Village  
Da Vincilaan 5  
Building Caprese, 8th floor  
1935 Zaventem  
Belgium  
Tel: +32.2.711.4680  
Fax: +32.2.725.2087

**Hologic Asia Pacific**  
One Temasek Avenue  
21-04 Millenia Tower  
Singapore 03192  
Tel: +65 6206 5833

**Hologic Australia**  
Suite 302, Level 3  
2 Lyon Park Road  
Macquarie Park NSW 2113  
Australia  
Tel: +61.2.9888.8000  
Fax: +61.2.9870.7555

PB-00510 Rev.007 (08/20) ©2020 Hologic Inc. All rights reserved. Hologic, 3D, Affirm, Brevera, CorLumina, Eviva, The Science of Sure, and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative.