



Artwork and Signature File for:

MAN-00469, "IHE INTEGRATION STATEMENT, SV-DX 3-0-X"

Artwork consists of:

- One (1) 8 ½ inch x 11 inch page.

REV AUTHORED BY	DATE		DANBURY, CT 06810		
J. SHARP	2/27/2006				
REV DRAFTED BY	DATE				
B. MOSES	3/1/2006	TITLE		DOCUMENT NUMBER	REV
PROPRIETARY INFORMATION: The content of this document is the exclusive property of Lorad and may not, without prior written permission of Lorad, be reproduced, copied or used for any purpose whatsoever.		AW, IHE INTEGRATION STATEMENT, SV-DX 3-0-X		AW-01271	001
		ARTWORK		SIZE A	SHEET 1 OF 1

IHE Integration Statement for SecurView_{DX} v3-0-X

IHE Integration Statement			
Vendor	Product	Version	Date
Hologic	SecurView Diagnostic Workstation	3-0-X	Feb 01, 2006
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Image Display	None	
Access to Radiology Information	Image Display	None	
Links to Standards Conformance Statements for the Implementation			
HL7	Not applicable		
DICOM	http://hologic.com/cc/dicom.htm		
Links for general information on IHE			
In North America: www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	