



Artwork and Signature File for:

MAN-00896, "IHE INTEGRATION STATEMENT, SECURVIEW-DX 6-X"

Artwork consists of:

- (1) 8 ½ inch x 11 inch pages.

REV AUTHORED BY	DATE		DANBURY, CT 06810	
A.FAND	05/23/08			
REV DRAFTED BY	DATE			
A.FAND	05/23/08			
PROPRIETARY INFORMATION: The content of this document is the exclusive property of Lorad and may not, without prior written permission of Lorad, be reproduced, copied or used for any purpose whatsoever.	TITLE		DOCUMENT NUMBER	REV
	AW, IHE INTEG. STMT, SV-DX 6-X		AW-02329	002
	ARTWORK		SIZE A	SHEET 1 OF 1

IHE Integration Statement for SecurView_{DX} v6-x

IHE Integration Statement			
Vendor	Product	Version	Date
Hologic, Inc.	SecurView _{DX} Diagnostic Workstation	6-x	May 2008
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Image Display	None	
Access to Radiology Information	Image Display	None	
Mammography Image	Image Display	None	
Mammography Image	Print Composer	None	
Consistent Presentation of Images	Evidence Creator	None	
Portable Data for Imaging	Portable Media Creator	None	
Links to Standards Conformance Statements for the Implementation			
HL7	Not applicable		
DICOM	http://hologic.com/cc/dicom.htm		
Links for general information on IHE			
In North America: www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	