

<b>IHE Integration Statement for Selenia</b>			
<b>Vendor</b>	<b>Product</b>	<b>Version</b>	<b>Date</b>
Hologic, Inc.	Selenia FFDM Acquisition Workstation (AWS)	5.1	February 2013
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Access to Radiology Information	Image Display	None	
Mammography Image Profile	Acquisition Modality	None	
Patient Information Reconciliation	Acquisition Modality	None	
Scheduled Workflow	Acquisition Modality	Assisted Acquisition Protocol Setting	
		PPS Exception Management Option	
		Broad Worklist Query	
		Patient Based Worklist Query	
Scheduled Workflow	Image Display	None	
<b>Links to Standards Conformance Statements</b>			
<b>HL7</b>	Not applicable		
<b>DICOM</b>	<a href="http://www.hologic.com/product-support/digital-mammography/selenia">www.hologic.com/product-support/digital-mammography/selenia</a>		
<b>Links for general information on IHE</b>			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	