



**Date:** 5/17/2019

**Author:** Software Engineering

**System:** Apex

**Subsystem:** Microsoft

Windows OS 7

**Subject:** Microsoft Windows 7 OS Patch release April 2019

---

**Purpose:**

To release the list of validated Microsoft critical patches on the Apex systems.

Note: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the "Sentinel Run-Time 7.80.0.0 " which can be located in the Windows Control Panel "Program and Features". If needed you can install the patch by ordering "APEX Service Dongle Patch Release" (ASY-12610). Please contact Hologic Support if you have any questions.

<b>System Type</b>	<b>Software Versions</b>
Horizon	APEX v5.5 – v5.6.0.5
Discovery	APEX v4.0 – v4.6.0.3
Workstation	APEX v4.0 – v5.6.0.5

**List of patches/KB's validated for April 2019:**

KB4493472	KB4493448	KB4493435	KB4490628
KB4490128	KB4489885	KB4486564	KB4486459
KB4483458	KB4483451	KB4480960	KB4474419

**List of all patches/KB's validated on the Apex System:**

KB976932	KB4493472	KB4493448	KB4493435
KB4490628	KB4490128	KB4489885	KB4486564
KB4486459	KB4483458	KB4483451	KB4483187
KB4480960	KB4474419	KB4471328	KB4467106
KB4462915	KB4457145	KB4457044	KB4457035
KB4343899	KB4338823	KB4284867	KB4103712
KB4096418	KB4093108	KB4074587	KB4056897
KB4054521	KB4048960	KB4041678	KB4019990
KB3177467	KB3172605	KB3161949	KB3159398
KB3156016	KB3150220	KB3139914	KB3139398
KB3138910	KB3138612	KB3127220	KB3126587
KB3124275	KB3122648	KB3110329	KB3109560
KB3109103	KB3108664	KB3108371	KB3101722
KB3097989	KB3093513	KB3086255	KB3078601
KB3075226	KB3074543	KB3072305	KB3071756
KB3067903	KB3061518	KB3060716	KB3059317
KB3055642	KB3046269	KB3046017	KB3045685
KB3042058	KB3037574	KB3035126	KB3031432
KB3030377	KB3023215	KB3022777	KB3021674
KB3020388	KB3011780	KB3010788	KB3004375
KB2992611	KB2991963	KB2984972	KB2978742
KB2977292	KB2973201	KB2973112	KB2972100
KB2968294	KB2943357	KB2937610	KB2931356
KB2923545	KB2912390	KB2911501	KB2894854
KB2894844	KB2893294	KB2884256	KB2871997
KB2868038	KB2864202	KB2862335	KB2862330
KB2862152	KB2861698	KB2857650	KB2852386
KB2847927	KB2840631	KB2834140	KB2830477

KB2813430	KB2813347	KB2807986	KB2786081
KB2770660	KB2742599	KB2736422	KB2729094
KB2727528	KB2705219	KB2698365	KB2690533
KB2670838	KB2667402	KB2654428	KB2639308
KB2631813	KB2621440	KB2620704	KB2604115
KB2579686	KB2574819	KB2564958	KB2560656
KB2552343	KB2533552	KB2532531	KB2506212
KB2491683	KB2479943		