

Aptima® BV and CV/TV Assays

■ An objective, comprehensive and accurate method for diagnosing the cause of vaginitis.

Vaginitis Is the Leading Reason for OBGYN Visits¹



Most women experience an episode of vaginitis in their lifetime.²



10 million clinical visits each year are due to symptomatic vaginitis.¹



90% of vaginitis is caused by:³

- Bacterial vaginosis (BV)
- Candida vaginitis (CV)
- *Trichomonas vaginalis* (TV)

Mixed Infections and Overlapping Symptoms Make Diagnosis a Challenge⁴

There are distinct treatments for each etiology, so the correct diagnosis is key to proper care.

37%

■ **37%** of women with BV are also infected with TV and/or *Candida* species.⁵

More than
50%

■ **More than 50%** of women diagnosed with BV experience recurrent symptoms within 12 months.⁶



NAAT Results in More Accurate Diagnosis²

Detect Up to

5x

More Infections

- **Detects 3 times more mixed infection cases** than clinical diagnosis with wet mount, culture, and Amsel's criteria.⁷
- **Detects mixed infections more frequently** than either clinical evaluation or probe testing.⁷
- **Detects 3 to 5 times more TV infections** than wet-mount microscopy.⁸

Test with exceptional performance. Aptima® BV and Aptima® CV/TV assays are highly sensitive NAATs that deliver accurate diagnosis of vaginitis.

Assay/Target	Aptima Multitest Vaginal Swab	Sensitivity (%)	Specificity (%)
Aptima® BV Assay	Bacterial vaginosis ⁹	Clinician-Collected	95.0
		Patient-Collected	97.3
Aptima® CV/TV Assay	<i>Candida spp.</i> ¹⁰	Clinician-Collected	91.7
		Patient-Collected	92.9
	<i>C. glabrata</i> ¹⁰	Clinician-Collected	84.7
		Patient-Collected	86.2
<i>Trichomonas vaginalis</i> ¹⁰	Clinician-Collected	96.5	
	Patient-Collected	97.1	

Streamline Testing for Sexual and Vaginal Health

The Aptima® portfolio powers efficient, streamlined testing with assay and collection consolidation on scalable automation.

Aptima® Multitest Swab Specimen Collection Kit

- Detect up to **7 infections** and disease states from just one vaginal swab.⁹⁻¹³
- **One-step**, direct-load workflow on the automated Panther® system for streamlined testing.
- Can pull up to **4 aliquots** from specimen transfer tube.
- **Better experience** for providers and patients with simplified collection and decreased callbacks.



References: 1. Kent HL. Epidemiology of Vaginitis. Am J Obstet Gynecol. 1991 Oct;165(4 Pt 2):1168-76. 2. Goje OL. Advancing the Diagnosis of Vaginitis. CLPMag. 2020. 3. Paladine HL. Vaginitis: Diagnosis and Treatment. Am Fam Physician. 2018 Mar 1;97(5):321-329. Sobel J. Approach to females with symptoms of vaginitis. UpToDate. Updated May 25, 2021. Accessed August 2021. www.uptodate.com 4. Anderson MR, Klink K, Cohrsen A. Evaluation of vaginal complaints. JAMA. 2004;291(11):1368–79. 5. Belley-Montfort L, Lebed J, Smith B, et al. Sensitivity of the Amsel's Criteria compared to the Nugent Score in absence and in presence of Trichomonas vaginalis (TV) and/or Candida SPP among women with symptomatic vaginitis/vaginosis. Sex Transm Infect. first published as 10.1136/sextrans-2015-052126.290 on 18 May 2015 6. Bradshaw CS, et al. High Recurrence Rates of Bacterial Vaginosis over the Course of 12 Months after Oral Metronidazole Therapy and Factors Associated with Recurrence. J Infect Dis. 2006;193(11):1478–1486. doi.org/10.1086/503780 7. Schwabke JR, et al. Diagnostic Performance of a Molecular Test Versus Clinician Assessment of Vaginitis. J Clin Microbiol. 2018; 56. doi:10.1128/JCM.00252-18 8. CDC. 2015 Sexually Transmitted Diseases Treatment Guidelines. https://www.cdc.gov/std/tg2015/trichomoniasis.htm. Updated June 4, 2015. 9. Aptima BV assay. US package insert AW-23712. Hologic, Inc., 2022. 10. Aptima CV/TV assay. US package insert AW-23713. Hologic, Inc., 2023. 11. Aptima Combo 2 Assay. US package insert AW-25929. Hologic, Inc., 2023. 12. Aptima Mycoplasma genitalium assay. US package insert AW-17946. Hologic, Inc., 2022. 13. Aptima Trichomonas vaginalis Assay, Panther system. US package insert AW-25942. Hologic, Inc., 2023.