

UNLOCK POTENTIAL: WHY UPGRADE?

MultiCare® Platinum System to Affirm® Prone Biopsy System



KEEP TECHNOLOGY UP TO DATE

- Exceptional 3D biopsy imaging^{1,2}
- Updated Software Version 1.1.3.0

ADVANCED IMAGING

- 12cm x14cm detector is 5.5x larger to visualize more tissue
- Biopsy lesions that are only visible under 3D imaging

ENHANCED WORKFLOW

- Translucent paddles help to see breast tissue for positioning
- Improved biopsy workflow with programmed needle parameters to automate calculations and reduce opportunity for error

To learn more, contact your Account Representative today!

References:

1. Compared to MultiCare Platinum system. Based on the fact that we have used the same detector technology as the Dimensions and introduction of tomo capability.
2. Schrading S, Martine D, Dirrachs T, et al. "Digital Breast Tomosynthesis-guided Vacuum-assisted Breast Biopsy: Initial Experiences and Comparison with Prone Stereotactic Vacuum-assisted Biopsy." Radiology. 2014 Nov 12. [Epub ahead of print].

MISC-09674 Rev 001 (04/24) ©2024 Hologic, Inc. All rights reserved. Printed in the USA. Specifications subject to change without notice. Hologic, MultiCare Platinum, Affirm and associated logos or trademarks are registered trademarks of Hologic, Inc. and/or subsidiaries in the United States and/or other countries. All other trademarks, registered trademarks, and product names are the property of their respective owners. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative or write to womenshealth@hologic.com.