



Dear Customer:

Thank you for your inquiry regarding continued software support for those Breast and Skeletal Health Hologic products which utilize Microsoft Windows 10 as their operating system. Windows IoT Enterprise was specially designed by Microsoft for fixed purpose devices, such as Medical Devices, that require tight change control and long-life support. Hologic leverages this special purpose operating system to achieve stability and reliability of the device, ensuring that operating system changes are limited to bug fixes and security fixes only.

While Microsoft has announced an end of support date of October 2025 for some Windows 10 operating systems, please note, Hologic does not utilize any of these versions of Windows 10 within any of our Breast and Skeletal Health products. Instead, Hologic products utilize various IoT Enterprise LTSC versions of Windows 10 which are not included in Microsoft's October 2025 end of support announcement. The specific versions of Windows 10 used in Breast and Skeletal Health Hologic products, as well as their currently scheduled end of support dates, are provided below:

<b><u>Product</u></b>	<b><u>Operating System</u></b>	<b><u>Currently Scheduled End-of-Support Date<sup>1</sup></u></b>
Trident HD	Windows 10 IoT Enterprise LTSC 2016	October 13, 2026
Affirm Prone Biopsy System Brevera Cenova Faxitron CoreVision/PathVision Dimensions/3Dimensions Horizon DXA Insight/Insight FD Securview Workstation Faxitron Vision DXA	Windows 10 IoT Enterprise LTSC 2019	January 9, 2029
SecurView Manager	Windows Server Standard 2016	January 12, 2027

---

<sup>1</sup> All end-of-support dates for the operating systems listed in this letter are determined by Microsoft. For the most up-to-date information regarding scheduled end-of-support announcements for Windows products, please contact Microsoft.

# HOLOGIC®

Advanced Workflow Manager (AWM)	Windows Server 2019	January 9, 2029
------------------------------------	---------------------	-----------------

As technology changes, we will continue working to ensure Hologic-manufactured equipment remains up to date with the latest software updates to meet our customers' needs. Hologic will begin evaluating Windows 11 IoT Enterprise LTSC for our Breast and Skeletal Health products when it becomes available later this year. If Hologic decides to switch to Windows 11 IoT Enterprise LTSC then we will provide our customers with ample time to prepare ahead of any scheduled end-of-support announcements.

Should you have any questions, please do not hesitate to reach out to your service or sales representative. Thank you for your continued partnership.

Sincerely,



Sr. Director, Global Service Marketing