



Marker Compatibility for Stereotactic Breast Biopsy



Biopsy Device	Catalog Number (SKU)	Gauge	Length	Shape	Metal	Bio-absorbable Carrier	Deployment Method	# per box	
Brevera	TUMARK-BREV-P-V-02	18G	12.35cm	VISION	Nitinol	No	Introducer	10	
	TUMARK-BREV-P-X-02			X	Nitinol	No	Introducer	10	
	TUMARK-BREV-P-Q-02		12.55cm	Q	Nitinol	No	Introducer	10	
	TUMARK-BREV-S-V-02		12.75cm	VISION	Nitinol	No	Introducer	10	
	TUMARK-BREV-S-X-02			X	Nitinol	No	Introducer	10	
	TUMARK-BREV-S-Q-02		12.95cm	Q	Nitinol	No	Introducer	10	
Eviva and Brevera	Smark-Eviva-10	12G	10cm	Mini Cork	Ti	Yes	Introducer	10	
	Smark-Eviva-2s-10			Top Hat	Ti	Yes	Introducer	10	
	Smark-E13-ss1		13cm	Buckle	SS	Yes	Introducer	10	
	Smark-E13-ss2			Infinity	SS	Yes	Introducer	10	
	Smark-E13-ss3			Stoplight	SS	Yes	Introducer	10	
	Smark-Eviva-13			Mini Cork	Ti	Yes	Introducer	10	
	Smark-Eviva-2s-13		Top Hat	Ti	Yes	Introducer	10		
	TriMark-Eviva-10		15G	10cm	Cork	Ti	No	Introducer	10
	TriMark-Eviva-2S-10	Hourglass			Ti	No	Introducer	10	
	TriMark-Eviva-13	13cm		Cork	Ti	No	Introducer	10	
	TriMark-Eviva-2S-13			Hourglass	Ti	No	Introducer	10	
	Eviva	TUMARK-E13-S-V-02	18G	12.95cm	VISION	Nitinol	No	Introducer	10
		TUMARK-E13-S-X-02			X	Nitinol	No	Introducer	10
		TUMARK-E13-S-Q-02		13.2cm	Q	Nitinol	No	Introducer	10
TUMARK-E13-P-V-02		13.3cm		VISION	Nitinol	No	Introducer	10	
TUMARK-E13-P-X-02				X	Nitinol	No	Introducer	10	
TUMARK-E13-P-Q-02		13.55cm		Q	Nitinol	No	Introducer	10	

Biopsy Device	Catalog Number (SKU)	Gauge	Length	Shape	Metal	Bio-absorbable Carrier	Deployment Method	# per box
ATEC	TriMarkTD-2S-13-09	9G	13cm	Hourglass	Ti	No	Needle Hub	10
	Smark-ATEC-13-09			Mini Cork	Ti	Yes	Needle Hub	10
	Smark-ATEC-2s-13-09			Top Hat	Ti	Yes	Needle Hub	10
	TriMark TD 1309			Cork	Ti	No	Needle Hub	10
	Smark-A3609-ss1	9G	36cm	Buckle	SS	Yes	Handpiece	10
	Smark-Atec-2s-36-09			Top Hat	Ti	Yes	Handpiece	10
	Smark-ATEC-36-09			Mini Cork	Ti	Yes	Handpiece	10
	TriMark TD 3609			Cork	Ti	No	Handpiece	10
	TriMarkTD-2S-36-09	12G	13cm	Hourglass	Ti	No	Handpiece	10
	TriMark TD 1312			Infinity	SS	Yes	Needle Hub	10
	TriMarkTD-2S-13-12			Cork	Ti	No	Needle Hub	10
	Smark-A1312-ss2			Hourglass	Ti	No	Needle Hub	10
	Smark-ATEC-13-12		Mini Cork	Ti	Yes	Needle Hub	10	
	Smark-ATEC-2s-13-12		Top Hat	Ti	Yes	Needle Hub	10	
	Smark-ATEC-36-12		36cm	Mini Cork	Ti	Yes	Handpiece	10
	TriMark TD 3612			Cork	Ti	No	Handpiece	10
	TriMarkTD-2S-36-12	Hourglass		Ti	No	Handpiece	10	

		
<p>Tumark® Marker</p> <p>Intelligently designed to provide long term ultrasound visibility, non-bio-absorbable, bio-compatible permanent marker with excellent visibility under mammography at deployment.¹</p>	<p>SecurMark® Marker</p> <p>The bio-absorbable suture-like net that carries a bio-compatible permanent marker, is designed to minimize movement within the biopsy cavity and enhance ultrasound visualization for three to four weeks post-biopsy.²</p>	<p>TriMark® Marker</p> <p>Reliable and simple bio-compatible Titanium metal markers, optimal for thin-breasted patients and superficial lesions.</p>

1. In 85% of marker placements, physicians stated that the ultrasound visibility was good to excellent upon deployment. In 56% of marker placements, physicians stated that ultrasound visibility was excellent upon deployment. Based on Tumark® Marker Data Collection Study DHM-06169, 3 clinicians at 3 hospitals for 90 marker placements. Hologic Data On File 2017 **2.** Based on Clinical Evaluation Report DHM-08585. Visibility is dependent on surrounding tissue, experience may vary.

www.hologic.com | 1.800.442.9892

Hologic Headquarters
United States / Latin America
 250 Campus Drive
 Marlborough, MA 01752
USA Tel: +1.508.263.2900
Sales: +1.781.999.7453
Fax: +1.781.280.0668
Email: info@hologic.com

Hologic Europe
 Everest (Cross Point)
 Leuvensesteenweg
 250A 1800
 Vilvoorde, Belgium
Tel: +32.2.711.4680
Fax: +32.2.725.2087

Hologic Asia Pacific Limited
 Unit 01-03A, 13/F
 No. 909 Cheung Sha Wan
 Road Cheung Sha Wan
 Kowloon, Hong Kong
Tel: +852 3748 7700
Fax: +852 3526 0723

Hologic Australia
 Suite 402, Level 4
 2 Lyon Park Road
 Macquarie Park NSW 2113
 Australia
Tel: +61.2.9888.8000
Fax: +61.2.9870.7555

MISC-01695 Rev.012 (6/24) ©2024 Hologic Inc. All rights reserved. Hologic, ATEC, Eviva, Brevera, SecurMark, TriMark, Tumark and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative.