

## Tuloc™ localization wire Solutions for Wire Localization

Introducing Tuloc™ localization wire, wire localization kits for interoperative breast lesion localization prior to surgery. Available for stereotactic, and ultrasound guided procedures, Tuloc's double hook anchor expands upon deployment providing a firm tissue anchor. In addition, it is designed to be repositioned, by simply pulling the wire back into the cannula and redeploying.

### Tuloc wire localization deployment systems

#### Tuloc



- With 20-Gauge puncture cannula
- Puncture cannula suitable for stereotactic surgery
- Sharp and echogenic cannula tip

#### Tuloc Premium



- With 20-Gauge puncture cannula and ergonomic handle
- For ultrasound-guided operation using one hand
- Sharp and echogenic cannula tip with ultrasonically visible marking

### Tuloc Localization Wire at a Glance

Marking for safe repositioning

Monofilament marking wire for intraoperative localization of findings

Double hook design with sharp tips



Marking wire made of implant material Nitinol



#### Secure anchoring

Double hook expands upon deployment and anchors firmly in the tissue



#### Repositioning

Wire can be pulled back into the cannula and released again



#### Single-handed operation

Tuloc Premium for ultrasound-guided operation using one hand

## The Tuloc localization wire kit consists of:

- A cannula with cm markings
- A marking wire, that is placed within the cannula shaft
- A fixation clip



Components	Material description
Cannula	Stainless steel 1.4301 / AISI 304
Cannula holder	PC
Marking wire	Nitinol
Tube	Stainless steel 1.4301 / AISI 304
Cover	LDPE + PP
Fixation clip	ABS

## The Tuloc Premium localization kit consists of:

- A cannula with cm markings and ultrasonic marking
- A handle with lever and slide button
- A marking wire, that is placed within the cannula shaft
- A fixation clip



Components	Material description
Cannula	Stainless steel 1.4301 / AISI 304
Marking wire	Nitinol
Tube	Stainless steel 1.4301 / AISI 304
Cover	LDPE
Cap	ABS
Nut screw	PP/TPE
Fixation clip	ABS

	Product code	Gauge	Cannula length
<b>Tuloc</b>	271660	20 G / 0.95 mm	1.97 inches
	271661	20 G / 0.95 mm	3.55 inches
<b>Tuloc Premium</b>	271880	20 G / 0.95 mm	1.97 inches
	271881	20 G / 0.95 mm	3.55 inches

DS-10143 Rev. 002 (10/23) Hologic Inc. ©2023 All rights reserved. Hologic, Tuloc, and the Science of Sure are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative.