

Marker Compatibility for Ultrasound Breast Biopsy

SecurMark® Biopsy Site Marker					Tumark® Biopsy Site Marker					CeleroMark, TriMark® Biopsy Site Markers	
-------------------------------	--	--	--	--	----------------------------	--	--	--	--	--	--



Biopsy Device	Catalog Number (SKU)	Gauge	Length	Shape	Metal	Bio-absorbable Carrier	Deployment Method	# per box
Celero, Sertera, ATEC with introducer & other core needle devices	CeleroMark	15G	10cm	Cork	Ti	No	End Deploy (Introducer or Free Hand)	10
	CeleroMark-2S			Hourglass	Ti	No	End Deploy (Introducer or Free Hand)	10
	Smark-Celero			Mini Cork	Ti	Yes	End Deploy (Introducer or Free Hand)	10
	Smark-Celero-2S			Top Hat	Ti	Yes	End Deploy (Introducer or Free Hand)	10
	Smark-U-ss1			Buckle	SS	Yes	End Deploy (Introducer or Free Hand)	10
	Smark-U-ss2			Infinity	SS	Yes	End Deploy (Introducer or Free Hand)	10
	Smark-U-ss3			Stoptlight	SS	Yes	End Deploy (Introducer or Free Hand)	10
	351260	18G	10cm	Conic	Nitinol	No	End Deploy (Introducer or Free Hand)	10
	351210			X	Nitinol	No	End Deploy (Introducer or Free Hand)	10
	351220			Q	Nitinol	No	End Deploy (Introducer or Free Hand)	10
	351230			Vision	Nitinol	No	End Deploy (Introducer or Free Hand)	10
	271510			Eye	Nitinol	No	End Deploy (Introducer or Free Hand)	10
	351280			U	Nitinol	No	End Deploy (Introducer or Free Hand)	10
	351262		12cm	Conic	Nitinol	No	End Deploy (Introducer or Free Hand)	10
	351212			X	Nitinol	No	End Deploy (Introducer or Free Hand)	10
	351222			Q	Nitinol	No	End Deploy (Introducer or Free Hand)	10
	351232			Vision	Nitinol	No	End Deploy (Introducer or Free Hand)	10
	271512			Eye	Nitinol	No	End Deploy (Introducer or Free Hand)	10
	351282			U	Nitinol	No	End Deploy (Introducer or Free Hand)	10
	ATEC side deploy			Smark-A1312-ss2	12G	13cm	Infinity	SS
Smark-ATEC-13-12		Mini Cork	Ti	Yes			Side Deploy (Needle hub)	10
Smark-ATEC-2s-13-12		Top Hat	Ti	Yes			Side Deploy (Needle hub)	10
TriMark TD 1312		Cork	Ti	No			Side Deploy (Needle hub)	10
TriMarkTD-2S-13-12		Hourglass	Ti	No			Side Deploy (Needle hub)	10
Smark-ATEC-13-09		9G	13cm	Mini Cork	Ti	Yes	Side Deploy (Needle hub)	10
Smark-ATEC-2s-13-09				Top Hat	Ti	Yes	Side Deploy (Needle hub)	10
TriMarkTD-2S-13-09				Hourglass	Ti	No	Side Deploy (Needle hub)	10
TriMark TD 1309				Cork	Ti	No	Side Deploy (Needle hub)	10

Marker Compatibility for Ultrasound Breast Biopsy

<p>Tumark® Marker</p> <p>Intelligently designed to provide long term ultrasound visibility, non-bio-absorbable, bio-compatible permanent marker with excellent visibility under mammography at deployment.¹</p>	<p>SecurMark® Marker</p> <p>The bio-absorbable suture-like net that carries a bio-compatible permanent marker, is designed to minimize movement within the biopsy cavity and enhance ultrasound visualization for three to four weeks post-biopsy.²</p>	<p>CeleroMark, TriMark® Markers</p> <p>Reliable and simple bio-compatible Titanium metal markers.</p>

1. In 85% of marker placements, physicians stated that the ultrasound visibility was good to excellent upon deployment. In 56% of marker placements, physicians stated that ultrasound visibility was excellent upon deployment. Based on Tumark® Marker Data Collection Study DHM-06169, 3 clinicians at 3 hospitals for 90 marker placements. **2.** Based on Clinical Evaluation Report DHM-08585. Visibility is dependent on surrounding tissue, experience may vary.

www.hologic.com | BSHSalesSupportUS@Hologic.com | + 1.800.442.9892

Hologic Headquarters
United States / Latin America
 250 Campus Drive
 Marlborough, MA 01752
USA Tel: +1.508.263.2900
Sales: +1.781.999.7453
Fax: +1.781.280.0668
Email: info@hologic.com

Hologic Europe
 Everest (Cross Point)
 Leuvensesteenweg
 250A 1800
 Vilvoorde, Belgium
Tel: +32.2.711.4680
Fax: +32.2.725.2087

Hologic Asia Pacific Limited
 Unit 01-03A, 13/F
 No. 909 Cheung Sha Wan
 Road Cheung Sha Wan
 Kowloon, Hong Kong
Tel: +852 3748 7700
Fax: +852 3526 0723

Hologic Australia
 Suite 402, Level 4
 2 Lyon Park Road
 Macquarie Park NSW 2113
 Australia
Tel: +61.2.9888.8000
Fax: +61.2.9870.7555