



**Date:** 04/23/2021

**Author:** Software Engineering

**System:** Apex

**Subsystem:** Microsoft Windows OS 7

**Subject:** Microsoft Windows 7 OS Patch release April 2021

---

---

**Purpose:**

To release the list of validated Microsoft critical patches on the Apex systems.

Note 1: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the “Sentinel Run-Time 7.80.0.0” which can be located in the Windows Control Panel “Program and Features”. If needed you can install the patch by ordering “APEX Service Dongle Patch Release” (ASY-12610). Please contact Hologic Support if you have any questions.

Note 2:

- Please check if January 2020 Windows 7 Monthly Patch Disk is installed on system, if not then first Install January 2020 Windows 7 Monthly Patch Disk (SFW-03075).
- Apply SFW-03338 Rev. 001 (Windows 7 ESU License Installation) to the system prior to installation of January 2021 Windows 7 Monthly Patch Disk (SFW-03375) to allow for updates to be installed, due to end of support for Microsoft Windows 7.
- Apply SFW-03338 Rev. 003 (Windows 7 ESU License Installation) to the system prior to installation of April 2021 Windows 7 Monthly Patch Disk (SFW-03430) to allow for updates to be installed, due to end of support for Microsoft Windows 7.

<b>System Type</b>	<b>Software Versions</b>
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

**List of patches/KB's validated for April 2021:**

KB5001392

KB5000851

KB4601363

KB5001335

KB5000800

KB4601275

**List of all patches/KB's validated on the Apex System:**

KB976932	KB4556843	KB4486459
KB5001392	KB4550965	KB4483187
KB5001335	KB4541500	KB4480960
KB5000851	KB4537829	KB4474419
KB5000800	KB4537813	KB4471328
KB4601363	KB4536952	KB4467106
KB4601275	KB4534314	KB4462915
KB4598289	KB4534310	KB4457145
KB4598279	KB4534251	KB4457044
KB4597239	KB4532945	KB4457035
KB4592510	KB4532932	KB4343899
KB4592503	KB4530692	KB4338823
KB4592471	KB4525233	KB4284867
KB4586827	KB4520003	KB4103712
KB4586805	KB4519108	KB4096418
KB4586768	KB4516033	KB4093108
KB4585204	KB4512486	KB4074587
KB4580970	KB4507704	KB4057116
KB4580387	KB4507456	KB4056897
KB4578952	KB4507004	KB4054521
KB4578623	KB4503269	KB4048960
KB4577053	KB4501226	KB4041678
KB4575903	KB4499175	KB4038779
KB4571719	KB4493448	KB4034679
KB4566371	KB4490628	KB4025337
KB4565539	KB4490128	KB4022722
KB4561669	KB4489885	KB4022619
KB4557900	KB4486564	KB4019990

KB4019263	KB3092601	KB2984972
KB4018073	KB3086255	KB2978742
KB4015546	KB3080079	KB2977292
KB4012212	KB3075226	KB2973201
KB3205394	KB3075220	KB2973112
KB3197867	KB3074543	KB2972211
KB3192391	KB3072305	KB2972100
KB3172605	KB3071756	KB2968294
KB3161949	KB3067903	KB2943357
KB3159398	KB3061518	KB2937610
KB3156016	KB3060716	KB2931356
KB3155178	KB3059317	KB2923545
KB3150220	KB3046269	KB2912390
KB3139914	KB3046017	KB2911501
KB3139398	KB3045685	KB2900986
KB3138910	KB3042058	KB2894854
KB3138612	KB3037574	KB2894844
KB3133977	KB3035126	KB2893294
KB3127220	KB3031432	KB2888049
KB3126587	KB3030377	KB2884256
KB3124275	KB3023215	KB2882822
KB3122648	KB3022777	KB2871997
KB3115858	KB3021674	KB2868038
KB3110329	KB3020388	KB2864202
KB3109560	KB3019978	KB2862330
KB3109103	KB3011780	KB2862152
KB3108664	KB3010788	KB2861698
KB3108371	KB3004375	KB2857650
KB3101722	KB2992611	KB2852386
KB3097989	KB2991963	KB2847927
KB3093513	KB2984976	KB2840631

KB2834140	KB2564958
KB2830477	KB2560656
KB2813430	KB2552343
KB2813347	KB2533623
KB2807986	KB2533552
KB2786081	KB2532531
KB2770660	KB2506212
KB2758857	KB2491683
KB2742599	KB2479943
KB2736422	KB2467173
KB2731771	
KB2729094	
KB2727528	
KB2705219	
KB2698365	
KB2690533	
KB2685939	
KB2670838	
KB2667402	
KB2656356	
KB2654428	
KB2653956	
KB2639308	
KB2631813	
KB2621440	
KB2620704	
KB2604115	
KB2585542	
KB2579686	
KB2574819	
KB2565063	