



Date: May 2, 2023

Author: Software Engineering

System: InSight

Subsystem: Microsoft Windows 10 Enterprise 2019 LTSC

Subject: Microsoft Windows 10 Enterprise 2019 LTSC Patch release April 2023

Purpose:

To release the list of validated Microsoft critical patches on the InSight systems.

| System Type | Software Versions |
|--------------------|--------------------------|
| Workstation | InSight v6.1.4 |

List of patches/KB's validated for April 2023:

KB5025229

List of all patches/KB's validated on the Insight System:

| | | |
|-----------|-----------|-----------|
| KB890830 | KB5022511 | KB4589208 |
| KB5025229 | KB5012170 | KB4577586 |
| KB5023702 | KB5005112 | |
| KB5022840 | KB5004424 | |