



Date: May 2, 2023

Author: Software Engineering

System: SecurView

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release April 2023

Purpose:

To release the list of validated Microsoft critical patches on the SecurView System.

List of patches/KB's validated for April 2023:

KB5025229

List of all patches/KB's validated on the SecurView System:

KB890830	KB5018419	KB5010351
KB5025229	KB5017315	KB5009557
KB5023702	KB5016623	KB5008873
KB5022840	KB5015811	KB5008218
KB5022511	KB5014692	KB5007206
KB5022286	KB5013941	KB5005112
KB5021237	KB5013641	KB4589208
KB5021126	KB5012647	KB4577586
KB5020866	KB5012170	KB4535680
KB5020627	KB5012128	
KB5019966	KB5011503	