



Date: May 11, 2026

Author: Software Engineering

System: AWM

Subsystem: Windows Server 2019 & Software Security Patches

Subject: Windows Server 2019 & Software Security Patch release April 2026

Purpose:

To release the list of validated Microsoft critical patches on the AWM System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for April 2026:

KB5082413	KB5082123	KB5068402
-----------	-----------	-----------

List of all patches/KB's validated on the AWM System:

KB5068402	KB5073723	KB5066586
KB5082123	KB5071544	KB5049608
KB5082413	KB5065225	KB5063877
KB5078752	KB5068791	KB5062557
KB5075904	KB5066143	KB5060531

KB5058392	KB5026362	KB5006669
KB5055519	KB5025229	KB5005698
KB5050533	KB5023702	KB5005573
KB890830	KB5022840	KB5005112
KB5053596	KB5022511	KB5005043
KB5052000	KB5022286	KB5004238
KB5050008	KB5021237	KB5003638
KB5049608	KB5021126	KB5003197
KB5048661	KB5020866	KB5001402
KB5046615	KB5020627	KB5001347
KB5044277	KB5019966	KB5001078
KB5044016	KB5019964	KB5000803
KB5043050	KB5018419	KB5000685
KB5042215	KB5018411	KB4601392
KB5041578	KB5017396	KB4601318
KB5040430	KB5017305	KB4598243
KB5039879	KB5017095	KB4593226
KB5039217	KB5016884	KB4589210
KB5037765	KB5016622	KB4589208
KB5036896	KB5016058	KB4586830
KB5036604	KB5015808	KB4580346
KB5035849	KB5014702	KB4577586
KB5034768	KB5014553	KB4577015
KB5034127	KB5014026	KB4576750
KB5033904	KB5013952	KB4571694
KB5033371	KB5013641	KB4565912
KB5032196	KB5012596	KB4565511
KB5031984	KB5012170	KB4562561
KB5031361	KB5011570	KB4561616
KB5030214	KB5011495	KB4557397
KB5029931	KB5010359	KB4551853
KB5029376	KB5009557	KB4550929
KB5029247	KB5009546	KB4549947
KB5028168	KB5008873	KB4540723
KB5027222	KB5008207	KB4539571
KB5027131	KB5007192	KB4537759

KB4535680

KB4514366

KB4505225

KB4494175

KB4494174

KB4132216

KB2565063

KB2467173