



Date: May 12, 2026

Author: Software Engineering

System: InSight

Subsystem: Microsoft Windows 10 Enterprise 2016 LTSC

Subject: Microsoft Windows 10 Enterprise 2016 LTSC Patch release April 2026

Purpose:

To release the list of validated Microsoft critical patches on the InSight systems.

System Type	Software Versions
Workstation	InSight v6.1.4

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for April 2026:

KB5082089	KB5082411	KB5082198
-----------	-----------	-----------

List of all patches/KB's validated on the Insight System:

KB5082089	KB5068864	KB5060954
KB5082198	KB5070247	KB5061010
KB5082411	KB5002789	KB5058524
KB5002848	KB5066136	KB5058383
KB5078938	KB5066584	KB5055661
KB5002839	KB5066836	KB5055521
KB5075902	KB5065427	KB5049614
KB5075999	KB5065687	KB5002702
KB5002829	KB5002763	KB5041773
KB5073447	KB5063871	KB5041576
KB5073722	KB5062560	KB5040562
KB5071543	KB5062799	KB5040434

KB5039885	KB5022838	KB5004424
KB5039334	KB5022503	KB5004238
KB5039214	KB5022289	KB5003638
KB5037763	KB5021235	KB5003197
KB5037016	KB5020873	KB5002662
KB5036899	KB5020614	KB5002619
KB5036609	KB5019964	KB5002542
KB5035962	KB5018411	KB5002520
KB5035855	KB5017396	KB5002497
KB5035849	KB5017305	KB5002464
KB5034862	KB5017095	KB5002406
KB5034767	KB5016622	KB5002369
KB5034119	KB5016058	KB5002323
KB5033910	KB5015808	KB5002223
KB5033904	KB5014702	KB5001402
KB5033373	KB5014026	KB5001347
KB5032391	KB5013952	KB5001078
KB5032197	KB5013625	KB5000803
KB5031989	KB5012596	KB4601392
KB5031362	KB5012170	KB4601318
KB5030504	KB5011570	KB4598243
KB5030213	KB5011495	KB4593226
KB5029924	KB5010359	KB4589210
KB5029242	KB5009546	KB4589208
KB5028169	KB5008207	KB4586830
KB5027219	KB5007192	KB4580346
KB5027123	KB5006669	KB4580325
KB5026363	KB5005698	KB4577586
KB5025228	KB5005573	KB4577015
KB5023788	KB5005112	KB4576750
KB5023697	KB5005043	KB4571694

KB4565912

KB4565511

KB4562561

KB4561616

KB4561600

KB4556813

KB4550994

KB4550929

KB4540723

KB4540670

KB4537764

KB4537759

KB4535680

KB4534271

KB4530689

KB4525236

KB4520724

KB4516115

KB4494175

KB4132216

KB2467173