



**Date:** May 12, 2026

**Author:** Software Engineering

**System:** SELENIA DIMENSIONS / 3DIMENSIONS

**Subsystem:** Windows 10 LTSC 2021 & Software Security Patches

**Subject:** Windows 10 LTSC 2021 & Software Security Patch release April 2026

---

**Purpose:**

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

**Latest SQLSERVER Versions Validated:**

**Note:** Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for April 2026:**

KB5082426

KB5082200

**List of all patches/KB's validated on the SELANIA DIMENSIONS /  
3DIMENSIONS System:**

KB5082200

KB5078885

KB5071546

KB5082426

KB5075912

KB5068402

KB5073724

KB5068781

KB5065225

KB5066135

KB5066791

KB5065429

KB5063709

KB5049613

KB5012170

KB5050533

KB5062554