



Date: May 12, 2026

Author: Software Engineering

System: SecurView

Subsystem: Windows 10 LTSC 2016 (v1607) & Software Security Patches

Subject: Windows 10 LTSC 2016 (v1607) & Software Security Patch release April 2026

Purpose:

To release the list of validated Microsoft critical patches on the SecurView System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for April 2026:

KB5082089	KB5082198	KB5082411
-----------	-----------	-----------

List of all patches/KB's validated on the SecurView System:

	KB5075999	KB5070247
KB5082089	KB5073447	KB5066584
KB5082198	KB5073722	KB5066836
KB5082411	KB5071543	KB5065427
KB5078938	KB5066136	KB5065687
KB5075902	KB5068864	KB5063871

KB5062560	KB5036899	KB5017095
KB5062799	KB5035962	KB5016622
KB5061010	KB5035855	KB5016058
KB5060954	KB5034862	KB5015808
KB5058524	KB5034767	KB5014702
KB5058383	KB5034119	KB5014026
KB5055661	KB5033373	KB5013952
KB5055521	KB5032391	KB5012596
KB890830	KB5032197	KB5012170
KB5054006	KB5031362	KB5011570
KB5053594	KB5030504	KB5011495
KB5052006	KB5030213	KB5010359
KB5050109	KB5029242	KB5009546
KB5049993	KB5028169	KB5008207
KB5048671	KB5027219	KB5007192
KB5046612	KB5026363	KB5006669
KB5044293	KB5025228	KB5005698
KB5043124	KB5023788	KB5005573
KB5043051	KB5023697	KB5005112
KB5041773	KB5022838	KB5005043
KB5041576	KB5022286	KB5004238
KB5040562	KB5021235	KB5003638
KB5040434	KB5020866	KB5003197
KB5039334	KB5019964	KB5001402
KB5039214	KB5018411	KB5001347
KB5037763	KB5017396	KB5001078
KB5037016	KB5017305	KB5000803

KB4601393	KB4540670
KB4601392	KB4539571
KB4601345	KB4537764
KB4601318	KB4537759
KB4601060	KB4535680
KB4598243	KB4534271
KB4593226	KB4530689
KB4589210	KB4525236
KB4589208	KB4521858
KB4586830	KB4520724
KB4580346	KB4519998
KB4580325	KB4516115
KB4577586	KB4516044
KB4577015	KB4512574
KB4576750	KB4512517
KB4571694	KB4509091
KB4565912	KB4503308
KB4565511	KB4494175
KB4562561	KB4494174
KB4561616	KB4346087
KB4561600	KB4132216
KB4557397	KB2565063
KB4556813	KB2467173
KB4550994	
KB4550949	
KB4550929	
KB4540723	