



**Date:** 08/30/2020

**Author:** Software Engineering

**System:** Apex

**Subsystem:** Microsoft Windows OS 10

**Subject:** Microsoft Windows 10 OS Patch release August 2020

**Purpose:**

To release the list of validated Microsoft critical patches on the Apex systems.

Note 1: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the "Sentinel Run-Time 7.80.0.0" which can be located in the Windows Control Panel "Program and Features". If needed you can install the patch by ordering "APEX Service Dongle Patch Release" (ASY-12610). Please contact Hologic Support if you have any questions.

| System Type | Software Versions |
|-------------|-------------------|
| Horizon     | APEX v5.6.0.7     |
| Workstation | APEX v5.6.0.7     |

**List of patches/KB's validated for August 2020:**

KB4571694

**List of all patches/KB's validated on the Apex System:**

- |           |           |           |
|-----------|-----------|-----------|
| KB4571694 | KB4540723 | KB4516044 |
| KB4565912 | KB4540670 | KB4512574 |
| KB4565511 | KB4537764 | KB4505422 |
| KB4562561 | KB4537759 | KB4503537 |
| KB4561616 | KB4535288 | KB4503308 |
| KB4561600 | KB4534271 | KB4503267 |
| KB4556813 | KB4530689 | KB4498947 |
| KB4550994 | KB4520724 | KB4494175 |
| KB4550929 | KB4516115 | KB4346087 |

KB4132216

KB4018073

KB2565063

KB4022619

KB2565063

KB2467173